MULTITORQUE SURGICAL WRENCH

1. PRODUCT DESCRIPTION

REFERENCE	DESCRIPTION
LLMQ	SURGICAL MULTITORQUE WRENCH

Kit for multitorque surgical wrench consisting of:

- Multitorque surgical wrench
- Adjustment key
- "LUBRICANT FOR MULTITORQUE WRENCH" (Ref. LLMQ-OIL) USDA H1 approved
- Instructions for use

The SURGICAL MULTITORQUE WRENCH with torque adjustment, is a dental device for manually tightening and loosening screws, prosthetic elements and implants. It is a precision instrument that can be disassembled, supplied unsterilized. In order to guarantee its correct operation, the wrench must be disassembled, disinfected, cleaned, greased and sterilized after each use and following the instructions for use. Therefore, those who use this instrument must always read the instructions for use before handling it. The handling and use of the product are the user's responsibility. Before each use always carry out a complete revision of the tool to guarantee the precision of the tightening torque. This instrument is not a measuring device. The precision of the device is ± 3.5 Ncm with a confidence interval of 95%. Its useful lifetime is approximately 3 years from the manufacturer's delivery date, always depending on the use. Discard the product if there is any physical or functional sign of deterioration. Below are the structures that form the SURGICAL MULTITORQUE WRENCH, all of them made in stainless steel except for the washer (9) (polyphenylene sulfide PPS) and the brake (10) (teflon PTFE).

2. INTENDED USE

The SURGICAL MULTITORQUE WRENCH with torque adjustment, is a dental device used for tightening and loosening screws, prosthetic elements and dental implants.

3. INDICATIONS FOR USE

The SURGICAL MULTITORQUE WRENCH is indicated for tightening the implant in the maxillary bone to an exact and known torque as part of the implantological surgery. The torques that the multitorque wrench can be adjusted to are 20, 30, 40, 50, 60 and 70 Ncm.

4. INSTRUCTIONS FOR USE

By turning the tightening torque adjustment screw, the wrench can be adjusted to the desired tightening torque. To adjust the wrench to the desired tightening torque the tightening torque adjustment screw is turned. To correctly regulate the tightening torque, turn the regular clockwise to reach the desired tightening torque and locate it at the exact reference point. Always try to make the reference located on the handle form a line with those of the regulator. If a torque lower than that used must be adjusted, loosen two turns below the new desired torque and re-tighten to the desired value. Once the torque has been set, fit the adapter of the multitorque surgical wrench and proceed to the insertion of the implant. When the preset torque has been reached, the wrench will turn without applying torque (with a characteristic sound). The "IN" inscription on the cap (3) indicates the position of the wrench that enables tightening; by reversing the device, the "OUT" inscription is to apply the unscrewing function.

5. CONTRAINDICATIONS

There are no particular contraindications for this product except for those derived from the surgical technique in itself, (anaesthetic, antibiotic therapy, etc.) and dental implantology: A lack of adequate quality and quantity of bone, infection, a lack of hygiene or cooperation by the patient, and systemic diseases (diabetes, etc.) are potential causes that can hinder the subsequent recovery from the surgical intervention. See MA-003: BTI implant system.

6. WARNINGS

-  The use of this product does not entail additional precautions in children, pregnant women and lactating women to those already contemplated in general implantology.
-  Do not use the wrench with a torque value higher than 70 Ncm as the tool can suffer damage.
-  The use of the multitorque surgical wrench is restricted to odontologists, stomatologists and maxillofacial surgeons. BTI regularly offers ongoing training courses for an optimum performance of the product.

7. PRECAUTIONS

This instrument must not be used for applications other than those mentioned in the sections "Product description" and "Intended use". Do not use the wrench with a torque value higher than 70 Ncm as the tool can suffer damage. No further precautions are considered for special risk groups such as children, pregnant women and breast-fed infants, in addition to those already contemplated in general in oral implantology techniques. Nevertheless, these situations must be taken into account by the clinician. Before use the sterilization of the multitorque surgical wrenches is recommended. This instrument must not be used for applications other than those mentioned in "Indications for use" or with equipment that affects the anticipated efficacy of the device. Any deterioration of the tightening, ratchet and torque mechanisms must be monitored by the staff responsible for the use and maintenance of the device. During the assembly, it is essential to not mix up the components that belong to different wrenches as the parts are not interchangeable.

-  Do not store the wrench with the spring compressed, it must be stored with the torque at the minimum value. This device must not be sterilised with the packaging supplied by the manufacturer.

8. ADVERSE EFFECTS

Use of the LLMQ surgical multitorque wrench does not entail adverse effects in itself, but its use is indirectly related to the risks involved in an implantological surgical intervention. Therefore the same adverse effects as those described for the use of implants must be observed. See MA-003: BTI implant system.

9. STORAGE AND HANDLING; CLEANING, ASSEMBLY/DISASSEMBLY AND STERILISATION

- STORAGE AND HANDLING -
The BTI multitorque surgical wrenches do not require special conservation or storage conditions.

-  It is recommended that they are sterilised before their first use.

- CLEANING

If during use residue from the surgical procedure is produced (blood secretions, tissue remains), the wrench must always be disassembled (see figs. 2 to 4) and immersed in an appropriate disinfection bath, scrupulously observing the manufacturer's recommendations (neutral detergent with disinfectant action).

-  During the same cleaning operation instruments of different materials should never be mixed. Avoid all contact between the different parts of the wrench. The materials of the multitorque wrench are compatible with most disinfectants in the market, except for those derived from the following chemical substances: iodine, chromic acid, permanganates and dichromates, hydrogen peroxide, peroxyacetic and peroxybenzoic acids and strong acids.

The disinfection bath facilitates cleaning, as dry residue causes corrosion.

After cleaning, rinse the parts well with water and with a nylon brush scrub the internal and external surfaces of the different parts of the wrench.

- DISASSEMBLY

Figure 2: Remove the spring from the handle. Completely unscrew the regulator (7) and remove the spring/end cap assembly (6). If necessary, use the end of the regulator to extract the spring; this operation must be done with precaution to avoid damaging the hexagonal tip (9).

-  Do not separate the end cap (6) from the spring and the washer (10) of the regulator (7). These two elements must never be disassembled.

Figure 3: Remove the screw with the hexagonal tip. With the hexagonal tip (9), remove the screw (5) applying slight pressure to remove the cap (3) of the head (4).

Figure 4: Remove the ratchet and the sleeve from the head. The head (4) contains two elements that must be removed: the ratchet (2) and the sleeve (1). The head and the handle form an assembly that cannot be separated. The screw (5) is attached to the body.

- ASSEMBLY

Figure 5: Insert the sleeve and the ratchet in the head. To assemble the wrench it is necessary to insert the following two elements in the order indicated: the sleeve (1) and the ratchet (2).

Figure 6: Greasing point

Moderately lubricate the area of contact between the teeth of the sleeve (1) and the ratchet (2) (1) as indicated in figure 6. Remove the remains of lubricant from the external surface of the wrench as excess lubricant leaves remnants on the surface of the instrument when sterilising it. Only use the lubricant "LUBRICANT FOR MULTITORQUE WRENCH" (Ref. LLMQ-OIL) that is supplied with the wrench.

Figure 7: Screw the cap on the head

Once elements (1) and (2) have been placed, fit the cap (3) adjusting it with the head (4). Turn the wrench always having the cap in its place. Firmly tighten the screw (5) with the tool attached in the regulator (7).

Figure 8: Insert the spring/end cap assembly in the handle and screw on the regulator. The assembly of the wrench is finished when the spring/end cap assembly (6) is inserted in the handle (4) and the regulator (7) is screwed to the latter. Once the assembly has been completed, a functional assessment must be made. To do so activate the sleeve. The instrument works perfectly if the wrench produces a normal ratchet noise. Before using the wrench always revise the assembly and the functionality of the wrench.

Figure 9: With the adjustment key the high tightening values can be reached more easily.

- STERILISATION

The instruments must be subjected to sterilisation with water vapour at 134°C/273°F during 18 minutes. Apply the regulations in force of the country the instrument is used in. During the sterilisation the wrench must be totally assembled and regulated at its minimum tightening value. (fig. 5 to 8)

-  Sterilise the wrench according to the sterilisation cycles recommended by the autoclave manufacturer. We recommend the use of devices equipped with vacuum pumps (type B) to reduce the risk of air pockets forming. This recommendation is particularly important for hollow instruments and to guarantee perfect drying. The dry air steriliser is not recommended as it can accelerate the ageing of the spring and therefore cause a modification of the tightening torque..
-  Make sure that no oxidized instrument or material is sterilised, as it increases the risk of the rest of the instruments becoming oxidised when sterilised together.

10. MAINTENANCE

After 1 year we recommend having the wrench checked by the BTI technical service.

11. LABELLING OF THE CONTAINER / SYMBOLS USED

LOT	Symbol for "BATCH CODE". "Batch number" is synonymous. This symbol is accompanied by the manufacturer's batch code, which appears beside the symbol.
REF	Symbol for "CATALOGUE NUMBER". "Reference number", "new order number" are synonymous.
CE 0123	CE marking according to European Directive 93/42/CEE corrected by Directive 2007/47/EC in relation to healthcare products
Rx only	The sale of this device is restricted solely to licensed health professionals
	Symbol of "MANUFACTURER".
	Symbol of "PRECAUTION".
	Symbol of "See instructions of use".

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