
BTI APNiA®

---

|           |                                   |    |
|-----------|-----------------------------------|----|
| <b>en</b> | Information for the patient ..... | 4  |
| <b>es</b> | Información para el paciente..... | 12 |
| <b>de</b> | Patienteninformation.....         | 20 |
| <b>fr</b> | Information pour le patient.....  | 28 |
| <b>hr</b> | Informacije za pacijente .....    | 36 |
| <b>it</b> | Informazioni per il paziente..... | 44 |
| <b>pt</b> | Informação para o paciente .....  | 52 |
| <b>sl</b> | Informacije za bolnika .....      | 60 |

### PATIENT INFORMATION SHEET:

*Note: Information related to the medical device BTI APNiA® and MA-153.*

### PRESENTATION AND COMPOSITION

Together with this instruction for use MA152: Patient Guide the patient must receive the following components:

- 1 BTI APNiA module 1 battery LR06 AA
- 1 Pulse oximetry sensor

|                                                |                                                                                                |            |
|------------------------------------------------|------------------------------------------------------------------------------------------------|------------|
| DEU610ANDB7<br>APNiA SENSOR<br>Disposable Sp02 | Medical Cables<br>C/Ferrocarril del Puerto nº18, planta 1, oficinas 3-4,<br>29002 Málaga Spain | CE<br>1639 |
|------------------------------------------------|------------------------------------------------------------------------------------------------|------------|

- 1 Disposable nasal cannula

|                                                          |                                                                         |    |
|----------------------------------------------------------|-------------------------------------------------------------------------|----|
| NC-002/30c<br>Nasal pressure monitoring<br>cannula adult | CNSAC MedShop GmbH<br>Am Sonnenstuhl 63 D-97236 Randersacker<br>Germany | CE |
|----------------------------------------------------------|-------------------------------------------------------------------------|----|

- 1 Elastic strap
- 1 Case
- 1 Strap hook

### WHAT IS IT AND WHAT IT IS USED FOR?

The BTI APNiA unit measures and saves the oxygen saturation, air flow in respiration, heart rate, position of the body, respiratory rate, thoracic movement and snoring by means of different sensors. Thus, with the information stored and with the help of the diagnosis software provided with the unit, it is possible to detect events such as: obstructive, central and mixed apnoea, hypopnoea, desaturations, artifacts or snoring, and draw up a detailed report. The aim is to obtain a prior diagnosis of the patient's clinical situation, which must then be confirmed by a specialist.

It is indicated for use by medical staff skilled in the diagnosis of respiratory disorders in sleep, in hospitals or clinics. At the same time, the electronic BTI APNiA unit is also indicated for domestic use by the patient.

### HOW IT IS USED?

It is a portable device that can be used by the patient himself, on himself, in the home during the hours of sleep by following these steps:

It is advisable to do this standing up, just before doing the test.

### **1. Putting on the Unit:**

- 1.1 Adjust the unit to your thorax, leaving it at the front. If you wish to use the optional hook, insert the strap and place it on the right side of the unit. If you do not wish to use it, hook the strap directly to the right side of the device.

The purpose of this hook is to set the tension of the strap on the hook, and only use the other end to adjust the unit.

#### **MEN:**


Place the device at the centre of your thorax.

#### **WOMEN:**

Place the device under your breasts.

- 1.2 Take the other end of the strap, pass it behind your waist around your back, and insert it in the other groove (left) of the module, like you have done with the other end. The side with Velcro must always be outward.
- 1.3 Adjust the unit at the site indicated for men and women, with the connections area up. Fasten the Velcro without tightening it excessively so it is not uncomfortable. Bear in mind that the straps can loosen with movement during the night, so it could be necessary to readjust it.

Avoid placing it too low (for example, on the abdomen) or on one side, as the device has accelerometers that measure the body position and movements, and, consequently, is designed to be placed at the site indicated provided this is possible.

The device must be placed at the site indicated for men and women and with the side with the connectors upwards.


We recommend placing the module over your clothes and adjusting the belt so it is not too tight.

We recommend using the most suitable strap for each patient. In any event, for subjects with a large abdominal perimeter, the large strap size (No. 3) should be used.

Once switched on, you should try to spend the first minute of study seated quietly on the edge of the bed, to gather the baseline parameters.

**The module must be well secured to your body, but not too tight.**


## **2. Placing the sensors:**

Insert the cannula into the nose, pass the tubes behind the ears and adjust the strap under the chin, so that it is neither too loose nor too tight. Secure the position of the cannula with adhesive plaster, to prevent it moving during the night.

### **2.1 Nasal cannula**

Remove the protective cap of the nasal flow connector and put it in the case.

Then screw the tube of the cannula into the round, transparent output. It is important to turn the connector until it is firmly secured.

**IMPORTANT:** The cannula must be inserted with the two small tubes toward the nose, not upwards. If they cause **too much** discomfort, you can use scissors to cut them **very carefully** maintaining the shape.

## 2.2 Pulse oximeter

Connect the pulse oximeter sensor cable to the module in the metal connection. Insert it gently, without damaging the connection pins, and push it all the way in.

Place the other end of the pulse oximeter around your index, middle or ring finger according to the drawings on the sensor. Do not stretch the adhesive tape when you are affixing the digital sensor. This could cause incorrect readings or blistering of the skin.

Detach the adhesive film from the sensor and dispose of it. Place your fingernail against the sensor. Wrap the sensor around your finger.

Ensure that the light emitter and detector are aligned with an imaginary vertical line passing through your finger from top to bottom. For best results, use medical tape or Elastoplast or secure the cable independently of the sensor, preferably around the base of your finger.

The sensor must be firm, but without excessive compression. We recommend the concept of "in position but not tight", if it is too tight, it will cause discomfort during the night.

Ensure there are no circumstances that alter or restrict your cardiac pulse, as the signal from the pulse oximeter could be lost (e.g. Blood pressure

monitor sleeve, nail polish or artificial fingernails, venous pulse or catheter, etc.).

Inspect the skin where the sensor is applied at least every 6-8 hours (adults) or 2-4 hours (neonates) to ensure the sensor is correctly placed and the skin beneath it is intact. You may experience adverse reactions of skin sensitivity to the self-adhesive pulse oximetry sensor. In this case, remove the sensor immediately.

Ensure that the connectors are connected to the device to prevent the pulse oximeter and the nasal cannula becoming tangled, because in certain circumstances, this can cause constriction.

### SWITCHING ON:

Press the on/off button for 1 second until you hear a short beep, the lights come on and the unit starts the recording. The light on the finger sensor will come on.

Note: Use alkaline batteries. It is not possible to use rechargeable batteries under any circumstances. The recording of 8 hours of study is guaranteed with the use of new alkaline batteries.

Note: It is essential not to sleep next to a person who snores, so that the microphone only captures the sound of the person being studied. Otherwise, it is preferable to perform the study in a different bedroom.

Once switched on, you should try to spend the first minute of study seated quietly on the edge of the bed, to gather the baseline parameters.

### RECORDING:

From this moment the BTI APNIA unit is recording and will not exceed 8 hours. The device will automatically be turned off after this time or when the battery has run flat.

The following icons will help on the patient to understand the meaning of each one:

In green: correct operation.

In yellow: incorrect operation. Sensor not in place or incorrectly positioned.

## NASAL CANNULA:

## OXIMETRY:

## OVERALL CONDITION:

## DATA RECORDING:

1. The testing sequence of the sensors lasts 30 seconds. During this time the lights on the front screen (nasal cannula and oximetry) blink in green until it is detected that the sensors are properly positioned. At that moment they will turn green and remain fixed.
2. If after 30 seconds they do not detect a signal, the icon corresponding to the incorrectly positioned or missing sensor will turn yellow.
3. The overall condition light on the upper screen will start in yellow and will remain that way while there are sensors that are not correctly positioned. It will change to green when it detects a correct signal in all the sensors.
4. Once the Testing Sequence is complete (30 seconds), the unit will begin to save data, regardless of whether any sensor is incorrect.
5. The lights will stay on for 3 more minutes, and then they will turn off to avoid excessive battery use, but the unit will continue to record normally.
6. If at any time during the night the patient wishes to verify that the BTI APNIA unit is still recording, simply check that the red light on the finger sensor is on or press the on/off button once and all the lights of module will come on for 5 seconds, in the status (green/yellow) they are in at that precise moment.

**If, during the data recording, the nasal cannula or the pulse oximeter sensor is dislodged, put them back following the guidelines in the section “Positioning BTI APNIA”. The data recording continues even if the cannula or sensor is dislodged.**

**If you accidentally turn the unit off during recording, it can be switched on as explained in the section “switching on”.**

## SWITCHING OFF:

1. Press the on/off button for 5 seconds. While holding the button down, all the lights will blink in green, then you will hear a long beep, the lights will go out and the unit will stop recording.

## AFTER THE RECORDING:

1. Remove the nasal cannula and the pulse oximeter sensor from your finger.
2. Loosen the elastic strap that secures the unit to your chest.
3. Disconnect the nasal cannula and the sensor from the BTi APNiA module and replace the protective cap on the nasal cannula connection. Disconnect the nasal cannula and the sensor from the BTi APNiA module and put the protective cap back on the nasal flow connector.
4. Separate the nasal cannula and dispose of it with the conventional organic waste.
5. Keep the rest of the unit in the case, as it was when you received it. The case has a defined shape so you can position the unit correctly. Please observe this to avoid transport damage.
6. Finally, return the case and the pulse oximeter sensor to the corresponding doctor/health care staff (so they can dispose of the pulse oximeter at a recycling centre: Waste Collection Centre).

**Both the nasal cannula and the pulse oximeter cable are single-use accessories, so you must discard them after use.**

### **What problems could be caused by using it? Possible side effects**

In prior studies in which BTi APNiA has been used, adverse effects due to its use were not described. Nevertheless, you may experience adverse reactions of skin sensitivity to the self-adhesive pulse oximetry sensor. In these cases, remove the sensor immediately.

### **Contraindications**

No specific contraindications have been found regarding the use of this product.

### **Warnings**

- The device and the pulse oximeter must never be exposed to direct contact with water or high levels of moisture. Ensure that no liquids enter these products.
- The pulse oximeter and nasal flow connectors are very sensitive; they must be protected while the device is not in use. Never blow on the connectors or handle them.

Ensure that the cannula and the pulse oximeter sensor are in good condition. If you observe any defects, do not use it and contact your doctor. If you detect any type of adverse reaction, contact your doctor immediately.

### **Patient responsibilities and precautions**

- Use this medical device according to the instructions for use and precautions described.
- This unit is prepared for one patient, so it is for use exclusively by the patient for whom it has been prepared. Never use the unit on other people before returning it to the corresponding healthcare staff.

- The BTI APNiA® unit should be kept out of reach and sight of children.
- If you detect any type of adverse reaction, contact a doctor immediately.

### Storage and conservation

The unit must be stored in a dry and dust-free environment, away from animals and insects. Otherwise, the unit may malfunction.

The patient shall not perform any maintenance on the BTI APNiA unit. The health professional will be responsible for replacing the battery.

### Meaning of the symbols used in this manual

|                                                                                   |                              |                                                                                   |                                                                                     |
|-----------------------------------------------------------------------------------|------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  | Precaution                   |  | See instructions for use                                                            |
|  | Nasal flow                   |  | Pulse oximetry                                                                      |
|  | ON/OFF                       |  | Do not touch the terminals, electrostatic discharge sensitive device. (ESD warning) |
|  | CE marking                   |  | Type BF applied part                                                                |
|  | See the instructions for use |  | Manufacturer of the unit                                                            |
|  | Date of manufacture          |  | Medical device                                                                      |
|  | Unique device identifier     |                                                                                   |                                                                                     |

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