
STOP & GUIDE® KIT

en	3
es	10
de	17
fr	24
it	31
pt	38
cs	45
da	52
el	59
bg	67
fi	75
hu	82
lt	89
nb	96
sv	103

STOP & GUIDE® KIT

1. PRODUCT DESCRIPTION

Stop & Guide® is a surgical system with components to guide the placement of dental implants. It is a versatile system that can be used as a "stopper kit" in conventional (analog) surgery, or as a "guided surgery kit" to perform piloted, semi-guided, and fully guided surgery.

2. INTENDED USE


Surgical instrument system intended for the placement of dental implants.

3. INDICATIONS FOR USE

The system is designed to allow fully guided placement of implants with a diameter of less than 5.5 mm whose drilling protocol has a drill with a diameter of less than 4.80 mm as the largest diameter drill and whose implant holder does not exceed 5.5 mm in diameter. Using semi-guided surgery protocols, it is possible to place larger diameter implants using this system, as long as the largest diameter drill is less than 4.80 mm.

4. STOP & GUIDE® SYSTEM COMPONENTS

- The Stop & Guide® kit includes the following components:
- A set of 12 drivers of different lengths.
- A connector for the drivers to the surgical contra-angle.
- A punch drill for *flapless* surgery.
- A driver tester, which is compatible with all 12 drivers.
- Specific implant holder for guided surgery for the narrow internal implant range.
- A *Stop & guide®* box

PUNCH DRILL Ø5	STOP & GUIDE® DRIVER TESTER
FP5	CMCG
	

STOP & GUIDE® DRIVER CONNECTOR	STOP & GUIDE® DRIVER	STOP & GUIDE® NARROW INTERNAL IMPLANT HOLDER
CDCG	DCG (followed by driver number)	PCGE
		

STOP & GUIDE® BOX
CCG

5. GUIDED SURGERY PROCEDURE WITH THE STOP & GUIDE® SYSTEM

To perform any of the 3 types of static guided surgery (piloted, semi-guided, or fully guided) it is necessary to have:

- A radiographic study of the patient (CBCT) and a scan of the surface of the patient's mouth.
- A surgical guide for guided surgery prepared prior to surgery, designed using surgical planning software.
- Sleeves to prepare the surgical guide for Guided Surgery, which can be:
 - Ref: SCG1. Manufactured by steco-system-technik GmbH & Co. Distributed by BTI. The sales unit contains 5 units.
 - Ref: SCG. Manufactured by BTI. The sales unit contains 5 units.
- A Stop & guide® guided surgery kit.
- Punch drill for flapless surgery 5mm diameter ref. FP5. Manufactured by BTI.
- Kit with the necessary drill for the dimensions of the planned implants, as well as the drilling sequence which can be consulted in the CAT218 catalogue.
- A WI-75 E/KM surgical contra-angle compatible with the connector to which the drivers are coupled. Manufactured by W&H. Distributed by BTI.

5.1 SURGICAL GUIDE FEATURES

Surgical guides for guided surgery are designed and made using different types of software which make it possible to superimpose or merge (match) the information from the scanning of the surface of the patient's mouth (in an STL file) and the radiographic study using CBCT (in a DICOM file).

After planning the case, the design of the guide must allow the inclusion of metal mastersleeves of the following dimensions: 5mm in height and 5.55 mm in internal diameter (Ref. SCG1); and 5mm in height and 5.57mm in internal diameter (Ref. SCG). The sleeves must be positioned with an offset (height from the implant platform to the upper margin of the sleeve) of 9mm in height.

5.2 SLEEVE GUIDED SURGERY

CP4 Titanium cylindrical sleeves with the following dimensions:

Ref. SCG

Outside diameter (D): 6.35mm

Inside diameter (d): 5.57mm

Length (L): 5.0mm


Ref: SCG1

Outside diameter (D): 6.35mm

Inside diameter (d): 5.55mm

Length (L): 5.0mm

For more information on the sleeves Ref. SCG1 see the STECO "STECOGUIDE" user instructions on their website <https://steco.de/download/#StecoGuide>

5.3 CONTRA-ANGLE W&H WI-75 E/KM

Contra-angle compatible with any micromotor in accordance with ISO 3964 "Dentistry. Coupling dimensions for handpiece connectors". The connection is standard and non-lit, so it does not require any particular specification.

For more information on the WI-75 E/KM surgical contra-angle, see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

5.4 PUNCH DRILL

The punch drills are drills bit is cylindrical with a cutting edge of 3.5mm-diameter (Ref. FP35), 4mm-diameter (Ref. FP4) or 5mm-diameter (Ref. FP5) cutting edge to pierce the gum. The outer diameter allows it them to be inserted and guided through the sleeves of the guide.

Ref: FP35

Ref: FP4

Ref: FP5

Indications for use:

They should be connected to the surgical motor, do not require irrigation, and the recommended drilling speed is 50-75 rpm.

They should not be used when the driver connector is mounted on the contra-angle, as it will not be able to reach the gums.

These products rotate in a clockwise manner (in normal surgical motor mode).

5.5 DRIVER TESTER

The driver tester is an element consisting of a handle and a head into which a driver can be screwed. This element attached to the driver facilitates its insertion through the sleeve to check that the initial part (smaller diameter) of the driver passes completely through the sleeve.

If the crest is not flat or the implant has been planned in a subcrestal position, it will not stop and this distance must be measured with the help of laser marking and then set to the depth at which you want to drill.

Indications for use:

This element attached to the driver facilitates its insertion through the sleeve to check the inclination of the bone crest margin and the depth at which the implant is planned to be positioned (in case the implant is planned to be positioned subcrestal).

5.6 DRIVER CONNECTOR

An element that fits on the head of a compatible contra-angle, WI-75 E/KM, with a part on which it is possible to attach the drivers.

5.7 DRIVERS

The drivers are cylindrical elements, which once coupled (attached) to the driver connector at the top, concentrically surround the drill used to prepare the implant bed. The outer diameter of the drivers allows them to be inserted and guided through the sleeves of the guide, so that the drill rotate freely inside, thus preventing the drill bit from rubbing against the sleeve of the guide.

The surgical box has 12 numbered drivers (from 2 to 13) to facilitate their identification and with depth marks (laser marking) that help in the placement of implants when they are planned in a subcrestal position. All drivers have an identical 9 mm long part that is inserted through the sleeve of the surgical guide and that allows the drill to be guided. The outer diameter of this initial part of the driver is adjusted to the inner diameter of the sleeve in the surgical guide. From these 9 mm, the area of the driver closest to the contra-angle has a larger diameter, so it creates a cap at the upper part of the sleeve on the guide. This driver zone varies in length from mm to mm between consecutive drivers.

The internal diameter of the drivers is 4.80 mm, which allows drilling with any drill with a smaller diameter.

Indications for use:

When placing a driver and a drill in the contra-angle, the part of the drill that is outside the driver determines the effective drilling depth, since the initial end of the driver stops the advance in depth when it reaches the bone crest, so that only the part of the drill that is outside the driver comes into contact with the bone. Changing a driver for a shorter one increa-

ses the drilling length and vice versa. The laser marking on the drivers indicates this effective drilling depth.

DCG2	DCG3	DCG4	DCG5	DCG6	DCG7
2.5mm	3.5mm	4.5mm	5.5mm	6.5mm	7.5mm

DCG8	DCG9	DCG10	DCG11	DCG12	DCG13
8.5mm	9.5mm	10.5mm	11.5mm	12.5mm	13.5mm
	2mm	3mm	4mm	5mm	6mm

5.8 GUIDED SURGERY IMPLANT HOLDER

The guided surgery implant holder is an implant holder whose outer diameter is adjusted to the inner diameter of the different sleeves in the guide, which allows the implant to be guided during its insertion. Its height is 9 mm, the same as the offset of the sleeve, which also provides control of the depth of the final positioning of the implant.

Indications for use:

The guided surgery implant holder is attached to the implant, when the aim is to perform a fully-guided surgery protocol. This implant holder is compatible with the prosthetic connection of BTI narrow internal implants (3.5 mm prosthetic connection).

6. INTENDED USER AND INTENDED PATIENT GROUP

The *Stop & Guide*® surgical system should be used by clinicians with specific training in dental implantology, guided surgery for implants, and in BTI dental products.

The group of patients envisaged are partially or totally toothless patients who are going to undergo oral implantology techniques and therefore oral or maxillofacial surgery.

7. CONTRAINDICATIONS

The *Stop & Guide*® surgical system is not suitable for use in patients who are allergic or hypersensitive to titanium.

The *Stop & Guide*® surgical system is not indicated when there is a limitation of oral opening in the patient and/or mesio/distal space of less than 6mm.

The *Stop & Guide*® surgical system should not be used in children who have not completed their growth and development stage.

For more information on the sleeves see the STECO "STECOGUIDE" user instructions on their website

<https://steco.de/download/#StecoGuide>

For more information on the WI-75 E/KM surgical contra-angle, see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

8. WARNINGS

The components of the *Stop & Guide*® surgical system as well as the SCG sleeve and the FP35 and FP4 punch drills must be sterilised prior to use in patients.

Lack of adequate bone quality and quantity, infections, lack of hygiene or cooperation on the part of the patient, teeth grinding and generalised diseases (diabetes, etc.) are potential causes that can aggravate the subsequent recovery from the surgical intervention.

In case of ingestion of the piece by the patient, refer them to the emergency department of the hospital for appropriate treatment.

All products must only be used with the corresponding original BTI components. Otherwise, the parts could be damaged.

The system is optimised for BTI drill, as in other systems compatibility with the marks on the drill cannot be ensured.

For more information on the sleeve Ref. SCG1 see the STECO "STECOGUIDE" user instructions on their website <https://steco.de/download/#StecoGuide>

For more information on the WI-75 E/KM surgical contra-angle, see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

9. CAUTIONS

NUMBER OF USES

BTI recommends a maximum of 13 uses for all drills, including the punch drill. This value serves as a reference, as the true service life may be different depending on the type of application and/or characteristics of the material being worked with (hardness of the bone or mould). For the rest of the kit components manufactured by BTI, a maximum of 50 uses is recommended.

In any case, instruments should be discarded as soon as any physical deterioration that may affect their proper functionality is noticed, e.g., instruments with damaged, deformed, and nicked blades, as they cause vibrations and imperfections in the preparation margins, as well as uneven surfaces.

The *sleeves* (SCG1 and SCG) are single-use products and therefore should not be reused. For more information on the sleeve

SCG1 see the STECO "STECOGUIDE" user instructions on their website <https://steco.de/download/#StecoGuide>

For more information on the WI-75 E/KM surgical contra-angle, see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

10. ADVERSE EFFECTS

The surgical procedure carries certain risks, such as inflammation of the implant site, loose sutures, short-term tenderness of the area, bruising, bleeding, loss of sensation in the lower lip and chin, nosebleeds, trismus, sinusitis, peri-implantitis, periodontitis, gingivitis, fistula, mucositis, anesthesia, paresthesia, dysesthesia, local pain, and bone resorption of the maxillary or mandibular crest. Although these effects are most likely only temporary, in some isolated cases the loss of sensation has been permanent. An infection can also occur around the implant, but can usually be resolved with local treatment.

For more information on the sleeve Ref. SCG1 see the STECO "STECOGUIDE" user instructions on their website <https://steco.de/download/#StecoGuide>

For more information on the WI-75 E/KM surgical contra-angle see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

11. INSTRUCTIONS FOR USE

1. Ensure that the surgical guide you have prepared after the case planning sits correctly in the patient's mouth.
2. If you want to perform a flapless surgery, place the appropriate punch drill in the counter-angle, without having mounted the driver connector, and start drilling in a clockwise manner at a low speed, 50-75 rpm until you reach the bone crest. If you perform the procedure by opening a flap, this step is not necessary.
3. Place a driver (preferably long) in the driver tester and pass it through the sleeve.
When the bone crest is flat and the planned implant positioning is juxtacrestal, the initial part of the driver passes completely through the appropriate sleeve until it reaches its upper part. When the crest is not flat or the desired positioning of the implant is subcrestal, the driver will not stop at the top of the sleeve and its laser marks (marked every 0.5 mm) will allow this distance to be estimated. The distance measured with the laser marking in this situation is the measurement to be added to the effective milling length, by selecting a shorter Driver (with a greater number), as explained in greater detail below.
4. Place the driver connector on the counter-angle (reference: WI-75 E/KM) ensuring the correct alignment between the head and the washer. Ensure that it is correctly mounted and that it does not show any lateral or vertical movement. The friction that can be felt when inserting it into

the drive head is normal to ensure correct guiding.

5. Place the longest driver (No. 2) and the largest diameter long drill that can pass through it (Ø4.75 BTI drill); rotate it at low speed and check the correct concentricity of all components.
6. Remove the Ø4.75 drill and driver no. 2.
7. Place the long starting drill and then the driver corresponding to the drilling depth you want to reach.
8. Remove the driver and then the long starter drill.
9. Continue with the drilling protocol corresponding to the implant to be placed and adapted to the patient's bone quality. Refer to CAT218 for BTI INTERNA implants for details. To do this, in addition to the drills, you must exchange the drivers depending on the depth to which you want to drill.

	DCG2	DCG3	DCG4	DCG5	DCG6	DCG7
LONG DRILLS	2.5mm	3.5mm	4.5mm	5.5mm	6.5mm	7.5mm
SHORT DRILLS						

DCG8	DCG9	DCG10	DCG11	DCG12	DCG13
8.5mm	9.5mm	10.5mm	11.5mm	12.5mm	13.5mm
	2mm	3mm	4mm	5mm	6mm

For the **juxta-crestal positioning** of an implant, the driver matching the length of the implant to be placed is selected. E.g., if a 6.5 mm implant is to be placed, driver no. 6 is selected. The driver system is optimised for use with *BTI Long drills*, so with other types of drills this numbering correspondence may not exist (depending on the length of the drills used). After placing the driver, check that the depth markings of the drill match perfectly with the lower end of the driver*. The drill are exchanged from the starting drill to the final diameter drill, following the usual drilling protocols and adapted to the quality of the patient's bone. Additionally, the 5 shortest drivers can also be used with the *BTI Short drills*, allowing a drilling depth range between 2.5 and 6.5 mm. For this reason, the drivers between number 9 and 13 have an additional number in their marking that corresponds to the drilling length when a short drill is used. E.g., the 11/4 driver indicates that it allows a drilling depth of 11.5 mm with a long drill and 4.5 mm when a short drill is used.

*Since there is some manufacturing tolerance in the length of the drills, it is important to confirm that the lower end of the driver matches the desired depth mark on the drill every time the drill is changed. The system for connection between the driver and contra-angle allows to overcome discrepancies of up to 1 mm, simply by not connecting the driver all the way.

For **sub-crestal positioning**, as for juxta-crestal positioning, the appropriate driver is initially selected for the length of the implant to be placed. Again, if a 6.5 mm implant is to be placed, driver no. 6 is selected and it is verified that the part of the drill that protrudes from the driver (and that determines the drilling depth) correctly matches the depth marks on the drill.

Since the default *offset* is 9 mm (distance between the implant platform and the top of the sleeve in the guide), in cases where the implant has been positioned below the ridge margin, when drilling, the initial end of the driver will encounter the crest before it can be fully inserted into the guide. The drivers are marked with laser marks every 0.5 mm to be able to measure this missing distance to be able to fully insert the driver in the sleeve (this measurement has been verified prior to starting the drilling protocol with the Driver tester). This measurement must be added to the drilling length (the part of the drill that is outside the driver) in order to correctly drill and position the subcrestal implant. If a 6.5 mm long implant has been positioned subcrestal by 2 mm during planning, the laser mark that we will see when the lower end of driver no. 6 hits the crest will be 2 mm. To reach the planned positioning, this distance of 2 mm should be added to the length of the implant when calculating the drilling depth. For this purpose, a driver that is 2 mm shorter can be used (no. 8 in this specific situation).

10. Once the drilling protocol is completed, remove the driver.
11. Proceed with the replacement of the conventional implant holder (pre-assembled in the implant) by the guided surgery implant holder, which can be done in different ways:

11.1 If you wish to perform a half-guided procedure, also remove the surgical guide from the patient's mouth, and proceed with the placement of the implant as if it were a conventional surgery.

11.2 If you want to follow a fully guided protocol, you have several options depending on the length of the implant and bone density (quality):

11.2.1 When standard length implants are inserted, and bone quality is favourable, insertion of the first 2-3 turns of the implant through the sleeve (without removing the guide) can be started as in a conventional surgery, with the pre-assembled implant holder and the implant holder connector of the implant. With the first coils inserted, the conventional implant holder is removed and the guided surgery implant holder is placed (**ensuring it sits perfectly in the connection**). The implant continues to be placed until the upper end

(with a larger diameter) of the guided surgery implant holder makes contact with the sleeve. The implant will have been positioned at this time in the planned direction and depth (9 mm *offset*). Remove the implant holder for guided surgery and the surgical guide.

11.2.2 When extra-short implants need to be placed or the bone quality is not favourable, another possible way to make this change is to start by placing the implant in a conventional way through the sleeve of the guide (as in the situation indicated above) until the first 2-3 turns are inserted and a degree of stability is reached in the implant which allows the change to be made. At this time, the guide is carefully removed, the conventional implant holder is removed, and then the guide is replaced. Through the guide, the guided surgery implant holder is connected (**ensuring it sits perfectly in the connection**) and the insertion is completed. Remove the implant holder for guided surgery and the surgical guide.

11.2.3 Another simple way to exchange the implant holders, in situations of low bone quality or when it is complex to perform the above procedures, is to start in the same way as indicated above (through the sleeve) and remove the guide after the insertion of the first 2-3 turns. Without the surgical guide in the mouth, the conventional implant holder is removed and replaced by the guided surgery implant holder. Next, the implant (of which only the first turns have been inserted) is disinserted next to the guided surgery implant holder. Then, the surgical guide is placed back in the mouth and the implant is completely inserted with the guided surgery implant holder through the sleeve in the guide. Remove the implant holder for guided surgery and the surgical guide.

12. Following these indications avoids the need for additional elements for the exchange of the implant holders and prevents contact of the implant surface with elements other than the patient's own bone.

11.1 REPROCESSING

The components of the *Stop & guide*® surgical system, the SCG sleeve and the punch drills are NON-STERILE when sold and must be sterilised prior to use. For products manufactured by BTI, the reprocessing protocol described in the CAT246 cleaning, disinfection and sterilisation guide is recommended.

The driver connector must be disassembled to be sterilised (head and washer separately).

For more information on the sleeve Ref. SCG1 see the STECO "STECOGUIDE" user instructions on their website <https://steco.de/download/#StecoGuide>

For more information on the WI-75 E/KM surgical contra-angle, see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

12. MAINTENANCE, STORAGE AND REMOVAL

The components of the *Stop & guide*® surgical system, the SCG sleeve and the punch drills do not require special storage conditions.

For more information on the sleeves see the STECO "STECOGUIDE" user instructions on their website <https://steco.de/download/#StecoGuide>

For more information on the WI-75 E/KM surgical contra-angle see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

WASTE DISPOSAL:

It is recommended to discard safely potentially contaminated or no longer usable medical devices as dental health care waste in accordance with local regulation and government legislation.

13. SERIOUS INCIDENT WARNING

If a serious incident occurs while or after using this product, please notify the manufacturer and your national authority. The manufacturer's contact details are: <https://www.qualifiedperson@bti-implant.es>.

14. OTHER INFORMATION

You must only accept devices with intact packaging and factory labelling. Contact your dealer if the container is opened or modified.

It is the responsibility of the user to examine the product and check that its components are in a suitable condition for the intended use. If the instrument does not show wear, it can be used.

For more information on the sleeve Ref. SCG1 see the STECO "STECOGUIDE" user instructions on their website <https://steco.de/download/#StecoGuide>

For more information on the WI-75 E/KM surgical contra-angle, see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

14.1 EUDAMED

When the European database on medical devices is available, the summary of safety and clinical performance of this medical device will be made available through Eudamed <https://ec.europa.eu/tools/eudamed>. To do this, it will be necessary to enter the UDI-DI code, available on the implant label.

The UDI code corresponds to a series of numerical or alphanumeric characters that allow the traceability of the medical device and is shown on the label in ICAD format and in a humanly readable manner.

For more information on the sleeve Ref. SCG1 see the STECO "STECOGUIDE" user instructions on their website <https://steco.de/download/#StecoGuide>

For more information on the WI-75 E/KM surgical contra-angle see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

15. SYMBOLS

For a description of the symbols that appear on product labels and in these instructions, please refer to the eMA087 guide.

For more information on the sleeves see the STECO "STECOGUIDE" user instructions on their website <https://steco.de/download/#StecoGuide>

For more information on the WI-75 E/KM surgical contra-angle see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

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