
IMPLANT EXTRACTION KIT

bg	Български	КОМПЛЕКТ ЗА ИЗВАЖДАНЕ НА ИМПЛАНТИ _____	(Страница 3)
cs	Čeština	SADA PRO EXTRAKCI IMPLANTÁTŮ _____	(Stránka 7)
da	Dansk	SÆTTIL UDTRÆKNING AF IMPLANTATER _____	(Side 11)
de	Deutsch	IMPLANTAT-ENTFERNUNGSSET _____	(Seite 15)
el	Ελληνικά	KIT ΑΦΑΙΡΕΣΗΣ ΕΜΦΥΤΕΥΜΑΤΩΝ _____	(Σελίδα 19)
en	English	SQUARE BTI PILLARS AND INTERFACES _____	(Page 23)
es	Español	KIT DE EXTRACCIÓN DE IMPLANTES _____	(Página 27)
fi	Suomen	IMPLANTTIEN IRROTUSSARJA _____	(Sivu 31)
fr	Français	KIT D'EXTRACTION D'IMPLANTS _____	(Page 35)
hu	Magyar	IMPLANTÁTUMOK KIVONÁSÁRA SZOLGÁLÓ KÉSZLET _____	(Oldal 39)
it	Italiano	KIT DI ESTRAZIONE IMPIANTI _____	(Pagina 43)
lt	Lietuvių	IMPLANTŲ IŠTRAUKIMO RINKINYS _____	(Puslapis 47)
nb	Norsk Bokmål	IMPLANTATFJERNINGSSETT _____	(Side 51)
pt	Português	KIT DE EXTRAÇÃO DE IMPLANTES _____	(Página 55)
sv	Svenska	IMPLANTATUTTAGNINGSSATS _____	(Sida 59)

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IMPLANT EXTRACTION KIT

1. PRODUCT DESCRIPTION

The implant extraction kit (KEXIM) is a presentation that combines several elements necessary for use in the extraction of implants. The system offers a non-invasive method for the clean and simple extraction of a wide range of implants of different diameters and lengths surgically osseointegrated into the maxillary or mandibular bone. The kit incorporates an autoclavable box designed to facilitate sterilisation of the material and tidy storage until the next use.

COMPONENTS

- **WRENCH (LLT200):** Counter-torque wrench for use in the extraction of dental implants.
 - **EXTRACTORS (EXIMxx/INEXIMxx):** Indicated for the extraction of osseointegrated implants in the maxillary or mandibular bone when they are no longer functional or when they are unnecessary in a new prosthetic situation.
- * EXTRACTORS (EXIM8 and EXIM9) are not included in the implant extraction kit (KEXIM), but their mode of use is the same.
- **EXTENSORS (ELLCaxx):** Extenders to be used in combination with hand instruments during surgical and prosthetic procedures for the insertion or extraction of BTI dental implants.
 - **HANDLE (MDPT):** Handle for hand instruments used for the placement of screw-retained prosthetic components.
 - **DRILLS (DRL-E / DRL-25):** Drills for removal of fractured or rounded BTI External (DRL-E) connection screws and explantation of osseointegrated implants. DRL-25 Used in conjunction with the INEXIML extractor.
 - **BOX (CEXIM):** Box for holding surgical instruments in their specific place on the base. Can be autoclaved. Not designed to maintain sterility. Not designed for instrument reprocessing.

2. INTENDED USE

Implant extraction kit intended for the removal of dental implants.

3. INDICATIONS FOR USE

The system is indicated for the removal of osseointegrated implants when they are no longer functional or are unnecessary in a new prosthetic situation. In most cases, these products avoid the use of the conventional trephine drill extraction protocol.

4. INTENDED USER AND INTENDED PATIENT GROUP

Intended user: Professional healthcare personnel in hospital or clinical settings with specific training in oral implantology and/or prosthodontics and training in BTI dental products.
Intended patient group: Adult males and females who require

the removal of dental implants and who do not have one or more of the conditions or pathologies listed in the contraindications.

5. CONTRAINDICATIONS

Situations in which the product should not be used under any circumstances because the risk is too high.

- Children who have not completed their growth and development stage.

Situations where use of the product may be possible under certain conditions, but with caution and medical supervision. It should always be assessed whether the benefit outweighs the risk.

- Pregnant or lactating women.

6. GENERAL WARNINGS

- **EXTRACTORS** and **DRILLS** are intended for one-off clinical use. The product should not be reused or re-sterilised after use in a clinical procedure, as there is a risk of infection, cross-contamination or alteration of the mechanical properties.
- If the product has been unpacked, but not used during the surgical procedure, it may be reprocessed for up to 6 cycles, following the cleaning, disinfection and sterilisation protocol described in the CAT246 guide. For reprocessing, it is essential to have previously removed the original protective packaging.
- The **EXTENSORS**, **HANDLE** and **KEY** are multiple clinical use products that must be reprocessed between each use exclusively according to the protocol described in the CAT246 guide. BTI recommends a maximum of 50 uses for **EXTENSORS** and **HANDLE**. For the **KEY** the number of uses is not limited by its use or by the cleaning, disinfection and sterilisation cycles. However, incorrect maintenance may alter the functionality of the device.
- Use only original BTI components. The use of non-original components may cause damage to the part or functional failure.
- Factors such as insufficient bone quality, infections, bruxism, poor hygiene, poor patient cooperation or systemic diseases (e.g. diabetes) may compromise recovery.
- In case of accidental ingestion of a component, go immediately to a hospital emergency department.

7. PRECAUTIONS

- Visually check the part before and after each use and cycle to ensure that it is still suitable and shows no signs of rust, degradation or breakage that may affect the safety or functionality of the product.
- Some references have specific compatibilities or additional considerations. See **INSTRUCTIONS FOR USE** for more information.

- Handle all biological specimens and sharps in accordance with the biosafety policies of the facility.

8. ADVERSE EFFECTS

All materials used in the manufacture of the components of the IMPLANTS EXTRACTION KIT are biocompatible.

If you observe any type of adverse effect, please notify the manufacturer at the following email address: **technicalmanager@bti-implant.es**. By reporting adverse effects, you can help provide more information about the safety of this medical device.


9. INSTRUCTIONS FOR USE

PRELIMINARY WARNINGS AND PRECAUTIONS


- Process ALL components prior to use only according to the protocol described in the CAT246 cleaning, disinfection and sterilisation guide.
- Sterilisation of the LLT200 must always be carried out with the key in the 'On' position.
- It is essential to carefully follow the instructions for use to avoid the application of lateral loads or leverage during tightening, as well as not to exceed the maximum torque specified, as this could lead to breakage of the parts or injury to the bone.
- When following the trephine protocol, the implant extractor must be removed from the implant before inserting the trephine drill.
- When using the DRL-E and DRL-25 drills, it is important to:
 - Vacuum the waste generated to avoid the risk of accidental ingestion.
 - Avoid overheating by controlling the milling time.
 - Carry out the milling in successive stages, gradually increasing the depth.
 - Refrain from milling in a single borehole.
 - Apply moderate pressure and maintain copious irrigation during the procedure.
 - Drills should be carefully inserted into the implant in "REVERSE" mode at 600 rpm prior to extractor placement. This procedure improves the efficiency and safety of the explantation.
- Extractor No. 0 (EXIME) must not be used with the LLT200 wrench included in the kit. It must be used exclusively with the LLMQ wrench (not included), without exceeding a torque of 70 Ncm.
- If, when using an extractor No. 1, 1A or 1B, it does not contact the upper conical part of the implant connection and only contacts the threaded area, the LLMQ wrench with a maximum torque of 70 Ncm must also be used, instead of the LLT200 wrench.
- Under no circumstances should you force the LLT200 counter-torque wrench or use additional tools to attempt to exceed the torque limit of 200 Ncm. Once this torque has been reached, the wrench mechanism is activated to prevent damage to tissue or instruments. Use the plastic tube provided to reverse the position of the mechanism to allow for reuse.

Wrench that has been disengaged upon reaching 200 Ncm

To reverse this disengagement and be able to use it again, insert it into the plastic reset tube.


Apply pressure to return the wrench to its original position (engaged).

Extract the wrench from the plastic tube.


9.1. GENERAL USE

The first step is to insert the extractor along its axial axis into the implant connection in an anticlockwise direction. To do this, use the transport handle which, when attached to the extender, becomes a tool similar to a manual screwdriver. This will provide you with the necessary precision to generate a good fit of the extractor in the desired position and axis.

Once the puller is in position, remove the handle and insert the 200 Ncm counter-torque wrench. Once connected, start turning it anticlockwise continuously and always maintain the axial angle of the key-extractor-implant to avoid bending movements that could fracture the implant or the extractor. Before the wrench disengages, maintain tension until the implant begins to pull out, as osseointegration is broken. Once this has occurred, continue to remove the implant gradually with the wrench until complete removal.

If the LLT200 counter-torque wrench jumps (turning movement over its joint), this indicates that the torque limit has been reached.

Once the implant has been removed, complete preservation of the implant site allows a new implant to be inserted in the same operation. To do this, it is recommended to drill the extraction site and then insert an implant with a larger diameter than the one removed, provided that there is no large crater. If after holding the tension for 20 seconds in an anticlockwise direction the implant has not started to come out, repeat up to three times. If it still does not come out, apply torque until the key disengages and, if this still does not work, have the key disengage up to three times before using the extraction kit with the trephine drill.

9.2. USE WITH TREPHINE


With implants that have been integrated for a long time and are firmly anchored, it is sometimes possible to disengage the key without having removed the implant.

In this case, drill around the implant with the trephine drills (not included in the kit) to a depth of 2-5 mm and then repeat the implant extraction. Drilling the coronal area of the implant drastically reduces the extraction torque of the implant, a fact based on the corticalisation that the bone undergoes in the first few turns of the implant, due to the stress it has had to withstand.


9.3. USE WITH COMPLEX CONNECTIONS REQUIRING DRILLING

The protocol for the use of the extraction drills is as follows:

1. Drill the internal housing of the implant connection with the DRL-E or DRL-25 drill (depending on the connection) in "REVERSE" (reverse rotation) at a speed of 600 rpm.
2. Drill as deep as possible so that the extractor can be correctly inserted and fixed in the most robust area of the thread (start of the connection).


3. Insert the appropriate puller depending on the connection or situation. Pre-drilling facilitates the insertion of the puller and increases its breaking torque. In the case of extractors 1-A and 1-B, make sure that the attachment to the implant is made with the most robust threaded area. Therefore, once the extractor 1A has been inserted to the end, continue drilling if you notice that the conical area with the largest diameter does not secure the extractor. If this is not possible, insert extractor 1B.

4. Continue the conventional extraction procedure with the counter-torque wrench until the implant is removed.

9.4. USE WITH DAMAGED CONNECTION

Sometimes, due to the improper insertion of the extractor or other circumstances, the connection may be damaged. When this occurs with a BTI INTERNA implant, it can be removed using the DRL-25 drill and the INEXIML extractor (no. 2). In these cases, drill the internal thread of the implant with the DRL-25 at 600 rpm until the 5 mm depth mark coincides with the platform.

Then insert the INEXIML implant extractor (no. 2) and continue the conventional extraction process with the counter-torque wrench until the implant is removed.

10. MAINTENANCE, TRANSPORT AND STORAGE

- The LLT200 counter-torque wrench should be checked by BTI technical service every 30 extractions. Contact your distributor for service management.
- Do not expose to direct sunlight.
- Store in a dry environment. Exposure to moisture may affect the integrity of the container.
- The product should not be removed from its original packaging or handled until ready for use.

- Do not use after the expiry date.
- Inspect the product before use to check that the container is not damaged. Do not use if the container is damaged.
- Contact your distributor if the container is open or modified.

11. DISPOSAL

- Dispose of all used materials as biological waste.
- Wear gloves during processing.
- Follow current regulations for infectious waste.

12. NOTIFICATION OF SERIOUS INCIDENTS

In accordance with Regulation (EU) 2017/745 (MDR), if a serious incident occurs during the use of this product or as a result of its use, the user or patient should immediately report it to the following address: technicalmanager@bti-implant.es

13. ADDITIONAL INFORMATION

BTI conducts regular training courses to optimise product performance.

The UDI code corresponds to a series of numeric or alphanumeric characters that allow traceability of the medical device and is displayed on the label in ICAD (Automatic Identification and Data Capture) format and in HRI (Human Readable Interpretation) format.

14. SYMBOLS

For a description of the symbols that appear on product labels and in these instructions, see MA087.

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