
IMPLANT EXTRACTION KIT

en	3
es	9
de	15
fr	21
it	27
pt	33
cs	39
da	45
et	51
ro	57
el	63
fi	69
hu	75
lt	81
lv	87
nb	93
sv	99

IMPLANT EXTRACTION KIT

1. PRODUCT DESCRIPTION

Implant extraction kit (KEXIM) is a presentation that combines several elements needed for use in implant removal. The system offers a non-invasive method for the clean, simple extraction of a wide range of implants of different diameters and lengths, surgically osseointegrated in the maxillary or mandibular bone. The kit includes extractors, extenders, drills, a wrench, a handle and a sterilisation box.

2. INTENDED USE

System intended for the implant removal.

3. INDICATIONS FOR USE

The indication for use of the system is the removal of osseointegrated implants when they are no longer functional or unnecessary in a new prosthetic situation. BTI has developed extractors for both internal and external connection implants. These products will, in most cases, avoid the use of the traditional trephine drill extraction protocol.

4. KEXIM SYSTEM COMPONENTS

The Implant extraction kit includes the following components:

- **EXTENDERS (ELLCAX).** A set of extenders to be used in combination with manual instruments during surgical and prosthetic procedures. The extenders are the intermediate parts that connect the torque wrench with the implant extractor. They have the same morphology as the conventional ratchet wrench extenders widely used in the placement of implants. They are sold in three lengths to adapt to all clinical situations.
- **HANDLE (MDPT).** A handle to allow/facilitate the manual handling of surgical instruments or prosthetic components. The transport handle (MDPT) enables us, in combination with the extenders, to fit the extractor correctly in the connection, making it easier to adjust.
- **DRILLS (DRL-E / DRL-25).** Drills for removal fractured or rounded screws and explantation of osseointegrated implants with or without rounded conexions. Drilling the interior of the implant increasing the depth of the connection, facilitating a deeper insertion of the extractor that secures its grip throughout its whole geometry.
- **An implant extraction kit box (CEXIM).** Contains the necessary instruments organized according to the technique to be performed. They consist of a lid and a base with holes for the placement of instruments and is used as sterilization container and to keep surgical instruments organized and stored.

200 NCM COUNTER-TORQUE WRENCH (LLT200)

EXTRACTORS (EXIMxx / INEXIMxx)

EXTENDERS (ELLCAX) AND HANDLE (MDPT)

DRILLS (DRL-E / DRL-25)

KEXIM BOX (CEXIM)

5. INTENDED USER AND INTENDED PATIENT GROUP

Implant extraction kit (KEXIM) must be used by dental health-care professionals with specific training in dental implantology and in BTI dental products.

The group of patients envisaged are partially or totally edentulous patients who are going to undergo oral implantology techniques and therefore oral or maxillofacial surgery.

6. CONTRAINDICATIONS

Implant extraction kit (KEXIM) should not be used in children who have not completed their growth and development stage.

7. WARNINGS

 Extractors (EXIMxx / INEXIMxx) and drills (DRL-E / DRL-25) are single-use products and are supplied NON-STERILE. They must be reprocessed before use according to section 10.5.

 Single-use extractors and drills that have not been used clinically may be reprocessed up to 6 times if the protective packaging has been removed and they have been organised and stored in the box provided.

 Extenders (ELLCAX), handle (MDPT) and torque wrench (LLT200) are multiple-use products and are supplied NON-STERILE. They must be reprocessed before use according to section 10.5.

 The kit box (CEXIM) is not designed to maintain sterility and is not intended for instrument reprocessing; it is intended only for sterilization.

 The extractor number 0 (EXIME) must not be used with the wrench included in the kit (LLT200). It must be used with the wrench LLMQ (not included in the kit) and must never exceed the torque of 70 Ncm.

 If the extractor (no. 1, 1A or 1B) does not make contact with the connection of the implant in the upper tapered part, in other words it only makes contact in the area of the thread, you must use the wrench LLMQ (not included in the kit) at a torque of 70 Ncm instead of the wrench LLT200.

 It is important to follow the instructions described in order not to exert leverage (lateral loads) in the process of applying torque as the parts may break. The torque limit must not be exceeded.

 If the LLT200 counter-torque wrench jumps (turning movement on its joint), this indicates that the torque limit has been reached. Never force the wrench or use any other element to exceed it.

 Drills must be inserted carefully in "REVERSE" mode at 600 rpm prior to inserting the extractor.

 The drilling must never be done in one single perforation; it must be done in different stages, gradually increasing the depth each time.

 Drilling must be done carefully, without applying excessive pressure and with abundant irrigation to avoid overheating. In addition, it is recommended to vacuum the remains of the drilled screw to avoid ingesting them.

 Lack of adequate bone quality and quantity, infections, lack of hygiene or cooperation on the part of the patient, teeth grinding and generalised diseases (diabetes, etc.) are potential causes that can aggravate the subsequent recovery from the surgical intervention.

 In case of ingestion of the piece by the patient, refer them to the emergency department of the hospital for appropriate treatment.

 All products must only be used with the corresponding original BTI components. Otherwise, the parts could be damaged.

 BTI does not guarantee the suitability of BTI implant extractors with all implants available on the market.

8. PRECAUTIONS

NUMBER OF USES

BTI recommends a maximum of 50 uses for extenders and handle. In any case, instruments should be discarded as soon as any physical deterioration that may affect their proper functionality is noticed, e.g., instruments with damaged, deformed, and nicked blades, as they cause vibrations and imperfections in the preparation margins, as well as uneven surfaces.

For LT200, BTI recommends having the wrench checked by the BTI technical service – SAT after 30 extractions. See section 11.


9. ADVERSE EFFECTS

The surgical procedure carries certain risks, such as inflammation of the implant site, loose sutures, short-term tenderness of the area, bruising, bleeding, loss of sensation in the lower lip and chin, nosebleeds, trismus, sinusitis, peri-implantitis, periodontitis, gingivitis, fistula, mucositis, anesthesia, paresthesia, dysesthesia, local pain, and bone resorption of the maxillary or mandibular crest. Although these effects are most likely only temporary, in some isolated cases the loss of sensation has been permanent. An infection can also occur around the implant, but can usually be resolved with local treatment.

10. INSTRUCTIONS FOR USE


10.1 GENERAL USE

The first step is to insert the extractor following its axial shaft into the connection of the implant, anti-clockwise. For this purpose, use the transport handle that when joined to the extender becomes a tool similar to a manual screwdriver that will provide the accuracy necessary to generate a good fit of the extractor in the position and the axis desired.

Once the extractor is in position, remove the handle and insert the 200 Ncm counter-torque wrench. Once connected, start turning it anticlockwise continuously and always maintaining the axial angle of the wrench-extractor-implant, to avoid flexion movements that could fracture the implant or the extractor.

Before the wrench disengages, maintain the tension until the implant begins to be extracted, as the osseointegration is broken. Once this has occurred, continue removing the implant gradually with the wrench until it is fully extracted.

Once the implant has been extracted, the complete preservation of the bed where it was located makes it possible to

insert a new implant in the same surgery. For this purpose, we recommend drilling the extraction area and then inserting an implant with a larger diameter than the one extracted, provided there is no large crater.

If after maintaining the tension for 20 seconds anticlockwise the implant has not begun to be extracted, repeat up to three times. If it still does not come out, apply torque until the wrench disengages and, if this still does not work, make the wrench disengage up to three times before using the extraction kit with the trephine.

10.2 USE WITH TREPHINING

On occasions, when faced with implants that have been integrated for a long period of time and are firmly anchored, the wrench can be disengaged without having extracted the implant.

In this case, trephine around the implant with the trephine drills (not included in the kit) to a depth of 2-5 mm and then repeat the removal of the implant. Trephining the coronal area of the implant drastically reduces the removal torque of the implant, a fact based on the corticalization that the bone suffers in the first turns of the implant, due to the stress it has had to withstand.

The extractor must be removed from the implant before inserting the trephine.

10.3 USE WITH COMPLEX CONNECTIONS THAT REQUIRE DRILLING


The protocol for the use of the extraction drills is as follows:

1. Drill the internal housing of the implant connection with the DRL-E or DRL-25 drill (depending on the connection) in reverse function at a speed of 600 rpm.
2. Drill as deep as possible to ensure the extractor can be correctly inserted and secured in the most robust area of the thread (start of the connection).

3. Insert the corresponding extractor based on the connection or situation. The prior drilling facilitates the insertion

of the extractor, increasing its breakage torque. In the case of 1-A and 1-B extractors, make sure that the attachment to the implant is done with the most robust thread area. Therefore, once the extractor 1A has been inserted to the end, if you see that the conical area with the greater diameter does not secure the extractor, continue drilling. If this is not possible, insert extractor 1B.


4. Continue with the conventional extraction process with the counter-torque wrench until the implant is extracted.

10.4 USE WITH DAMAGED BTI CONNECTION

Sometimes, due to the poor insertion of the extractor or other circumstances, the connection may become damaged. When this happens with a BTI INTERNA implant, it may be extracted using the drill DRL-25 and the extractor INEXIML (No. 2).

In these cases, drill the internal thread of the implant with the DRL-25 at 600 rpm until the 5 mm depth mark coincides with the platform.

After that, insert the implant extractor INEXIML (no. 2) and continue with the conventional extraction process with the

counter-torque wrench until the implant is extracted.

10.5 REPROCESSING FOR NON-STERILE AND/OR MULTIPLE-USE DEVICES

BTI only recommends the reprocessing protocol described in the CAT246 cleaning, disinfection, and sterilisation guide.


 LLT200 must be sterilized in engaged position. When the wrench reaches the torque, use the plastic tube which is included to reverse the disengaged.

Wrench that has been disengaged upon reaching 200 Ncm.

To reverse this disengagement and be able to use it again, insert it into the plastic reset tube.

Apply pressure to return the wrench to its original position (engaged).

Extract the wrench from the plastic tube.

11. MAINTENANCE, STORAGE AND REMOVAL

The components of the implant extraction kit (KEXIM) do not require special storage conditions.

It is recommended to discard safely potentially contaminated or no longer usable medical devices as dental health care waste in accordance with local regulation and government legislation.

For verification and maintenance of the LLT200 by BTI technical service – SAT, please contact your distributor.

12. SERIOUS INCIDENT REPORTING

If a serious incident occurs while or after using this product, please notify the manufacturer and your national authority. The manufacturer's contact details are:
<https://www.qualifiedperson@bti-implant.es>.

13. OTHER INFORMATION

It is the responsibility of the user to accept devices with the packaging and labelling intact. Contact your distributor if the packaging is open or you suspect that it has been tampered with.

It is the responsibility of the user to examine the product and check that its components are in a suitable condition for the intended use. If the instrument does not show wear, it can be used.

BTI conducts regular training courses for optimal product performance.

EUDAMED

When the European database on medical devices is available, the summary of safety and clinical performance of this medi-

cal device will be made available through Eudamed <https://ec.europa.eu/tools/eudamed>. To do this, it will be necessary to enter the UDI-DI code, available on the device label.

The UDI code corresponds to a series of numerical or alpha-numeric characters that allow the traceability of the medical device and is shown on the label in ICAD format and in a humanly readable manner.

14. SYMBOLS

For a description of the symbols that appear on product labels and in these instructions, please refer to the eMA087 guide.

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