
# HEALING ABUTMENTS, HEALING SCREWS AND HEALING CAPS

**WARNING:**  
Please, read this manual carefully before connecting the centrifuge.

<b>en</b>	HEALING ABUTMENTS, HEALING SCREWS AND HEALING CAPS.....	(Page 4)
<b>es</b>	PILARES DE CICATRIZACIÓN, TORNILLOS DE CIERRE Y TAPAS DE PROTECCIÓN.....	(Página 7)
<b>de</b>	GINGIVAFORMER, DECKSCHRAUBEN UND SCHUTZKAPPEN.....	(Seite 10)
<b>fr</b>	PILIERS DE CICATRISATION, VIS DE FERMETURE ET CAPUCHONS DE PROTECTION .....	(Page 13)
<b>it</b>	PILASTRI DI CICATRIZZAZIONE, VITI DI SERRAGGIO E CAPPUCCI DI CICATRIZZAZIONE.....	(Pagina 16)
<b>pt</b>	PILARES DE CICATRIZAÇÃO, PARAFUSOS DE FECHO E TAMPAS DE PROTEÇÃO .....	(Pagina 19)
<b>ro</b>	BONTURI DE CICATRIZARE, ȘURUBURI DE STRĂNGERE ȘI CAPE DE VINDECARE.....	(Pagină 22)
<b>cs</b>	HOJÍČÍ ABUTMENTY, HOJÍČÍ ŠROUBY A HOJÍČÍ ČEPIČKY.....	(Stránka 25)
<b>sv</b>	LÄKNINGSFÖRANKRINGAR, LÄKNINGSSKRUVAR OCH LÄKNINGSKÅPOR.....	(Sida 28)
<b>no</b>	HELBREDENDE BROSTØTTER, SKRUER OG HETTER .....	(Side 31)
<b>lv</b>	DZIEDINOŠIE BALSTI, DZIEDINOŠĀS SKRŪVES UN DZIEDINOŠIE UZGAĻI.....	(34. lpp)
<b>lt</b>	GIJIMO GALVUTĖS, GIJIMO SRAIGTAI IR GIJIMO DANGTELIAI.....	(Puslapis 37)
<b>hu</b>	GYÓGYULÓFELÉPÍTMÉNYEK, GYÓGYULÓCSAVAROK ÉS ÍNYFORMÁZÓK .....	(Oldal 40)
<b>da</b>	TANDKØDSFORMER, HELINGSSKRUER OG HELINGSHÆTTER.....	(Side 43)

# HEALING ABUTMENTS, HEALING SCREWS AND HEALING CAPS

## 1. PRODUCT DESCRIPTION

BTI healing abutments, healing screws and healing caps are provisional attachments of commercially pure titanium with a polish that allows the soft tissues to heal around them. Some healing caps contain biocompatible peek (polyetheretherketone).

## 2. INTENDED USE

Healing abutments, healing screws and healing caps are temporary single use structures used to model the gingival tissue during the healing process and to protect the implant hexagon or internal thread, and to protect the tongue from rubbing.

Based on this indication usage of these is as follows:

### HEALING ABUTMENTS

Indication: They are screw-mounted into the implant, in the first surgery for implants fitted with one surgery, or when the second surgery is performed, in the case of two-stage implants.

*\*BTI recommends using the narrowest abutment, provided this is possible, to gain gingival width, so the papilla will have greater irrigation and consequently, less recession.*

### HEALING SCREW

Indication: They are screwed into the implant for the protection of the platform and connection in the case of a two-stage surgery.

### HEALING CAPS

Indication: They are screwed into the transepithelial abutment, and are used from when this is fitted until the fitting of the definitive or provisional prosthesis, to prevent the gum from growing over the abutment and protection of the tongue.

## 3. INSTRUCTIONS FOR USE

The healing abutments are screwed directly into the implant, with the large hexagonal connection, at a torque of 20 Ncm.

The healing screws are screwed directly into the implant.

The connection of the healing screws in this manual uses a small hexagon, and they are screwed at 10 Ncm.

The healing caps are screwed directly into the transepithelial abutments, using the large hexagonal connection, at a torque of 10 Ncm.

Interna healing abutments are supplied sterile. Some healing abutments have a small groove which can be used if the hexagon of the abutment is stripped, so it can be removed with a flat screwdriver.

The Externa and Tiny® healing abutments, healing screws and healing caps are supplied non-sterile.

 Sterilization before usage is recommended (at 134 °C during 6 min, with a pre-vacuum autoclave).

## 4. CONTRAINDICATIONS

There are no particular contraindications for this product except for those derived from the surgical technique in itself, (anaesthetic, antibiotherapy, etc.) and dental implantology: A lack of adequate quality and quantity of bone, infections, a lack of hygiene or cooperation by the patient, and systemic diseases (diabetes, etc.) are potential causes that can hinder the subsequent recovery from the surgical intervention. See IFU: BTI implant system.

BTI healing abutments, healing screws and healing caps are not suitable for use in patients who are allergic or hypersensitive to the material used (commercially pure titanium).

## 5. WARNING

 The use of this product does not entail additional precautions in children, pregnant women and lactating women to those already contemplated in general implantology.

 The reuse and/or sterilization of healing abutments, healing screws and healing caps may cause risk of infection and/or cross contamination.

 The use of healing abutments, healing screws and healing caps is restricted to odontologists, stomatologists and maxillofacial surgeons.

 BTI regularly offers ongoing training courses for an optimum performance of the product.

## 6. PRECAUTIONS

No further precautions are considered for special risk groups such as children, pregnant women and breast-fed infants, in addition to those already contemplated in general in oral implantology techniques. Nevertheless, these situations must be taken into account by the clinician.

## 7. ADVERSE EFFECTS

Not applicable. In any event those related to a surgical implant operation.

## 8. MRI SAFETY INFORMATION

 **MR Conditional**

 **Warning: The RF safety of the device has not been tested. The patient may only be imaged by landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the RF coil.**

Non-clinical worst case configuration review has determined that materials used in the BTI Interna Dental Implant UnicCa® Systems are safe from MR-induced displacement force for a given spatial gradient and also from MR-induced torque perspective. Therefore, a patient with this device can be scanned safely in an MR system under the following conditions:

<b>Device Name</b>	BTI Dental Implant System UnicCa®
<b>Static Magnetic Field Strength (B0)</b>	≤ 3.0T
<b>Maximum Spatial Field Gradient</b>	30 T/m (3,000 gauss/cm)
<b>RF Excitation</b>	Circularly Polarized (CP)
<b>RF Transmit Coil Type</b>	For body transmit coil, landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the coil.  Extremity T/R coils permitted.  Excludes Head T/R coil.
<b>Operating Mode</b>	Normal Operating Mode in the allowed imaging zone
<b>Maximum Whole-Body SAR</b>	2 W/kg (Normal Operating Mode)

<b>Maximum Head SAR</b>	Not evaluated for head landmark
<b>Scan Duration</b>	No specific constraints due to implant heating

## 9. LABELING OF THE CONTAINER / SYMBOLS USED

Symbol	Symbol Title and description	Symbol source
	BATCH CODE. Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1 (5.1.5)
	CATALOGUE NUMBER. Indicates the manufacturer's catalogue number so that the medical device can be identified	ISO 15223-1 (5.1.6)
	DO NOT REUSE. Indicates a medical device that is intended for one single use only.	ISO 15223-1 (5.4.2)
	USE-BY DATE. Indicates the date after which the medical device is not to be used.	ISO 15223-1 (5.1.4)
	STERILIZED USING IRRADIATION. Indicates a medical device that has been sterilized using irradiation	ISO 15223-1 (5.2.4)
	SINGLE STERILE BARRIER SYSTEM. Indicates a single sterile barrier system	ISO 15223-1 (5.2.11)
	CE MARKING according to European Directive 93/42/CEE corrected by Directive 2007/47/EC regarding medical devices (European requirement)	Council directive 93/42/CEE amended by Directive 2007/47/EC
<b>Rx only</b>	Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.	21 CFR 801.109(b) (1)
	MANUFACTURER. Indicates the medical device manufacturer	ISO 15223-1 (5.1.1)
	DO NOT USE IF THE CONTAINER IS DAMAGED. AND CONSULT INSTRUCTIONS FOR USE. Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information	ISO 15223-1 (5.2.8)
	DO NOT RE-STERILISE. Indicates a medical device that is not to be resterilized	ISO 15223-1 (5.2.6)

	<p><b>CAUTION.</b> Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.</p>	<p>ISO 15223-1 (5.4.4)</p>
	<p><b>CONSULT INSTRUCTIONS OF USE OR CONSULT ELECTRONIC INSTRUCTIONS FOR USE.</b> Indicates the need for the user to consult the instructions for use.</p>	<p>ISO 15223-1 (5.4.3)</p>
	<p><b>DATE DE MANUFACTURE.</b> Indicates the date when the medical device was manufactured</p>	<p>ISO 15223-1 (5.1.3)</p>
	<p><b>NON-STERILE.</b> Indicates a medical device that has not been subjected to a sterilization process</p>	<p>ISO 15223-1 (5.2.7)</p>
	<p><b>MEDICAL DEVICE.</b> Indicates the item is a medical device</p>	<p>ISO 15223-1 (5.7.7)</p>
	<p><b>MR CONDITIONAL.</b> An item with demonstrated safety in the MR environment within defined conditions including conditions for the static magnetic field, the time-varying gradient magnetic fields and the radiofrequency fields.</p>	<p>ASTM F2503 (FIG. 6)</p>

**NOTE:** The devices should be accepted only if the factory packaging and labeling arrive intact. Contact your distributor if the package is opened or altered.

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