
HEALING ABUTMENTS, HEALING SCREWS AND HEALING CAPS

en	HEALING ABUTMENTS, HEALING SCREWS AND HEALING CAPS.....	2)
es	PILARES DE CICATRIZACIÓN, TORNILLOS DE CIERRE Y TAPAS DE PROTECCIÓN.....	5)
de	GINGIVAFORMER, VERSCHLUSSSCHRAUBEN UND SCHUTZKAPPEN.....	8)
fr	PILERS DE CICATRISATION, VIS DE FERMETURE ET COUVERCLES DE PROTECTION.....	11)
cs	HOJIVÉ ABUTMENTY, UZAVÍRACÍ ŠROUBY A OCHRANNÉ UZÁVĚRY.....	14)
it	PILASTRI DI CICATRIZZAZIONE, VITI DI CHIUSURA E COPERTURE DI PROTEZIONE.....	17)
pt	PILARES DE CICATRIZAÇÃO, PARAFUSOS DE FECHO E TAMPAS DE PROTEÇÃO.....	20)
ro	BONTURI DE CICATRIZARE, ȘURUBURI DE ÎNCHIDERE ȘI CAPACE DE PROTECȚIE.....	23)
da	HEALING ABUTMENTS, LUKKESKRUER OG HELINGSHÆTTER.....	26)
hu	FEDŐCSAVAR, ZÁRÓCSAVAROK ÉS VÉDŐSAPKA.....	29)
lv	DZĪŠANAS BALSTI, NOSLĒDZOŠĀS SKRŪVES UN AIZSARGUZGAĻI.....	32)
lt	GIJIMO ATRAMOS, FIKSAVIMO VARŽTAI IR APSAUGINIAI DANGTELIAI.....	35)
no	HELENDE BROSTØTTER, LUKKESKRUER OG HELINGSHETTER.....	38)
sv	LÄKANDE DIREKT, LÅSSKRUVAR OCH SKYDDSKAPSLAR.....	41)
fi	PARANEMISJATKEET, KIINNITYSRUUVIT JA SUOJAKUORET.....	44)
et	TUGIPOSTID HAAVA PARANDAMISEKS, LUKUSTUSKRUVID JA KAITSEKORGID.....	47)
el	ΘΕΡΑΠΕΥΤΙΚΑ ΣΤΗΡΙΓΜΑΤΑ ΕΠΟΥΛΩΣΗΣ, ΒΙΔΕΣ ΑΣΦΑΛΙΣΗΣ ΚΑΙ ΠΡΟΣΤΑΤΕΥΤΙΚΑ ΚΑΛΥΜΜΑΤΑ.....	50)

HEALING ABUTMENTS, HEALING SCREWS AND HEALING CAPS

1. PRODUCT DESCRIPTION

HEALING ABUTMENTS are commercially pure titanium abutments that are screwed to implants. Healing abutments have an extra smooth surface finish to promote healing of the surrounding soft tissue. They have an EXTERNA or INTERNA connection. EXTERNA healing abutments are not supplied sterile and are single-use only. INTERNA healing abutments are supplied sterile and are single-use only.

HEALING SCREWS are single-use temporary screws made of commercially pure titanium that close the implant connection to protect the implant during two-stage procedures. Healing screws are not supplied sterile.

HEALING CAPS are single-use provisional structures to protect the abutment and the patient's mouth prior to the fitting of a prosthesis. The head of the crown has a hexagonal insert that allows it to be easily transferred into position. They are made of commercially pure titanium although some caps contain bio-compatible PEEK (polyetheretherketone). Healing caps are not supplied sterile.

2. INTENDED USE

Healing abutments, healing screws, and healing caps are single-use temporary structures to shape the gingival tissue during its healing process, protect the implant hexagon or internal thread, as well as to protect the tongue from chafing.

3. INDICATIONS FOR USE

HEALING ABUTMENTS

Indication: To be screwed to the implant during the first procedure in the case of one-stage implants, or during the second procedure in the case of two-stage implants. In the case of reduced healing abutments, which have a narrow supragingival geometry, these are indicated for thin gingival biotypes, when it is necessary to obtain more tissue thickness during the healing time.

**BTI recommends using the narrowest abutment, whenever possible, in order to gain gingival width, so that the papilla will have*

more irrigation and consequently less retraction.

HEALING SCREW

Indication: To be screwed to the implant to protect the plate and implant connection for two-stage procedures.

HEALING CAPS

Indication: To be screwed to the transepithelial; they are used from the placement of the transepithelial to the placement of the definitive or provisional prosthesis, to prevent the gum from growing over the transepithelial and to protect the tongue from chafing.

4. INTENDED USER AND INTENDED PATIENT GROUP

Healing abutments, healing screws and healing caps must be used by healthcare professionals with specific training in oral implantology and training in BTI dental products.

The target group of patients are partially or totally edentulous patients who are to undergo oral implantology techniques and thus oral or maxillofacial surgery.

5. CONTRAINDICATIONS

BTI healing abutments, healing screws and healing caps are not intended for use in children who have not completed their growth and developmental stages.

BTI healing abutments, healing screws and healing caps are not suitable for use in patients who are allergic or hypersensitive to titanium.

6. WARNINGS

Prosthetic components must only be used by suitably trained professionals. The prosthetic restoration techniques that are required to use these products correctly are very specialised and complex. The use of an improper technique can lead to implant damage and subsequent failure.

Sterile products must not be used if their original packaging is damaged as they may have lost their sterility.

 Do not use the products after the expiry date stated on the packaging, as the sterility and/or integrity of the product cannot be ensured after this date.

 EXTERNA healing abutments, healing screws and healing caps must be sterilised prior to use in patients. See section 8.2 STERILISATION.

 Reuse and/or re-sterilisation of healing abutments, healing screws and healing caps may result in risk of infection and/or cross-contamination. Furthermore, it may cause a change in mechanical characteristics and/or dimensions leading to implant damage and subsequent failure.

 Apply the specified torque to the product. Applying a higher torque than specified may result in thread breakage of the piece. Conversely, if a lower torque is applied, loosening of the product may occur.

 In case of ingestion of the piece by the patient, refer the patient to the emergency department of the hospital for appropriate treatment.

7. PRECAUTIONS

 This product must only be used with the corresponding original BTI components. Otherwise, the piece could be damaged.

8. ADVERSE EFFECTS

All materials used in the manufacture of BTI's prosthetic and laboratory components are biocompatible, however, some patients may experience allergic reactions or irritation due to undiagnosed sensitivity to titanium. This reaction may result in eczema, gingival hyperplasia, and pain / inflammation / erythema with granulomatous reaction around the contact surface.

Pregnancy and lactation: No additional precautionary measures to those observed in general oral implant surgery are considered, however, these special conditions must be taken into account by the clinician prior to placing the product.

9. INSTRUCTIONS FOR USE

9.1. HOW TO USE

Healing abutments are screwed directly to the implant by means of the large hexagonal connection at a tightening torque of 20 Ncm. Some abutments incorporate a small slot to be used in cases where the hex of the abutment is rounded so that it can be unscrewed with a flat screwdriver. Healing screws are screwed directly to the implant at a tighten-

ing torque of 10 Ncm. The Healing screw connection is a small hex.

Healing caps are screwed onto the transepithelials, using the large hex connection to a torque of 10 Ncm.

9.2. STERILISATION

INTERNA healing abutments are supplied sterile.

EXTERNA healing abutments, healing screws and healing caps are supplied unsterilised and must be sterilised prior to use in patients.

BTI only recommends the sterilisation protocol described in the CAT-246 cleaning, disinfection, and sterilisation guide.

10. MRI SAFETY INFORMATION

NOTE: This section does not apply to the Canadian market.

 This symbol does not apply to the Canadian market.

 RF safety of the BTI Implant System: implants + prosthetic components has not been proven. The patient can only be imaged by marking a reference point at least 30 cm away from the implant or by ensuring that the implant is outside the RF coil.

The non-clinical review of the worst-case scenario has determined that the materials used in the BTI Dental Implant Systems are safe from the point of view of MR-induced displacement force for a given spatial gradient and also from the perspective of MR-induced torsion. Therefore, a patient with this device can be safely scanned on an MRI system under the following conditions:

Static magnetic field intensity (B0)	≤ 3.0 T
Maximum spatial gradient of the field	30 T/m (3,000 gauss/cm)
RF Excitation	Circular polarisation (CP)
RF transmission coil type	In the case of the body coil, a reference point must be marked at least 30 cm from the implant or ensure that the implant is outside the coil. Allowed end T/R coils Excluding the head T/R coil.
Operating mode	Normal mode of operation in the permitted image area
Maximum SAR throughout the body	2 W/kg (Normal operating mode)

Maximum SAR in the head	Head reference point not assessed
Scan duration	There are no specific limitations due to the heating of the implants

11. MAINTENANCE, STORAGE, AND REMOVAL

Prosthetic components do not require special storage conditions as long as they are kept in their original packaging.

WASTE DISPOSAL:

1. Wear gloves to reduce the risk of infections (latex or nitrile).
2. Medical devices containing contaminated blood must be disposed of in containers suitable for potentially infectious waste.
3. General hygiene directives and legal regulations governing the proper disposal of infectious material must be complied with.

12. NOTICE REGARDING SERIOUS INCIDENTS

If a serious incident occurs during the use of this product or as a result of its use, notify the manufacturer and your national authority. The manufacturer's contact details are: <https://www.qualifiedperson@bti-implant.es>.

13. FURTHER INFORMATION

It is the responsibility of the user to examine the product and check that its components are in a suitable condition for the intended use.

BTI conducts regular training courses for optimal product performance.

EUDAMED

When the European database on medical devices is available, the summary of safety and clinical performance of this medical device will be made available through EUDAMED <https://ec.europa.eu/tools/eudamed>. This will require entering the UDI-DI code, which is available on the product label.

The UDI code corresponds to a series of numeric or alphanumeric characters that allow the traceability of the medical device and is displayed on the label in ICAD format and in a human readable form.

14. LABELLING OF THE PACKAGING / SYMBOLS USED

For a description of the symbols that appear on product labels and in these instructions, please refer to the MA-087 guide.

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