

PROSTHETIC AND LABORATORY COMPONENTS

en		4
es		8
de		13
fr		18
it		23
pt		27
cs		32
ro		36

PROSTHETIC AND LABORATORY COMPONENTS

1. PRODUCT DESCRIPTION

BTI prosthetic and laboratory components are attachments that enable the preparation of single or multiple dental prostheses, whether provisional or definitive.

2. INDICATIONS FOR USE

BTI prosthetic and laboratory components are used in the partial or complete restoration of the teeth in partially or fully edentulate patients, respectively.

The prosthetic components are used as accessories to dental implants, either as a component of support or fixation for the different fixed or removable prosthetic solutions that make it possible to improve or restore the patient's chewing function.

The laboratory components, however, are used as auxiliary elements in the manufacture of the different prosthetic solutions.

3. INSTRUCTIONS FOR USE

All the prosthetic and laboratory components by BTI are supplied NON-STERILE. It is recommended that the products that are used directly in the patient are sterilised before use at 134 °C for 6 min (with pre-vacuum autoclave).

Products used in the prosthetic laboratory do not require sterilisation prior to use.

All the prosthetic and laboratory components of BTI are single-use products, so they must not be reused or re-sterilised.

Section 9 indicates the specific instructions for use for each of the product families that make up the BTI pros-

thetic and laboratory components.

4. CONTRAINDICATIONS

There are no particular contraindications for this product except for those derived from the surgical technique in itself, (anaesthetic, antibiotherapy, etc.) and dental implantology: A lack of adequate quality and quantity of bone, infections, a lack of hygiene or cooperation by the patient, and systemic diseases (diabetes, etc.) are potential causes that can hinder the subsequent recovery from the surgical intervention. See the Instructions for Use of BTI implants.

5. WARNINGS

 The reuse and/or sterilisation of prosthetic components may cause risk of infection and/or cross contamination.

 The reuse and/or re-sterilisation of the prosthetic components can lead to a loss of the mechanical characteristics and/or dimensions that may cause the prosthetic restoration to fail.

6. PRECAUTIONS

The prosthetic and laboratory products must be used only by professionals with sufficient training. The prosthetic restoration techniques required to use these products correctly are very specialised and complex. Inadequate technique can cause the prosthesis to fail.

7. ADVERSE EFFECTS

All the materials used in the manufacture of the BTI

prosthetic and laboratory components are biocompatible, nevertheless some patients may present allergic reactions or hypersensitivity to some of the materials and their components.

The materials in which the different prosthetic components are manufactured are indicated in section 9.

8. STORAGE, CONSERVATION AND HANDLING

BTI prosthetic and laboratory components do not require special storage conditions.

9. SPECIFIC INSTRUCTIONS FOR USE

a. Aesthetic abutment posts /SQUARE

Aesthetic abutment posts for Zirconium, with different gingival heights and bio profiles (except height 0) for Interna/Externa/Multi-Im/Unit connections and all the BTI platforms. They incorporate a new anti-rotational system that facilitates the prosthetic work as it is not necessary to machine the male of the connection. They are made in commercially pure titanium.

- With casting cylinder: They are used for conventional techniques of prosthesis preparation by casting in conjunction with the plastic casting cylinder for single prostheses.
- Without casting cylinder: They are used as a machined base in titanium on which the zirconium cylinders can be modelled with systems for manual drilling or CAD-CAM systems for the preparation of prostheses or pillars with aesthetic materials such as Zirconium.

In both cases, the titanium component and the part made of zirconium will be cemented in the laboratory.

b. Aesthetic temporary abutments

Plastic provisional cylinders (PEEK) in white, indicated for preparing provisional prostheses in aesthetic areas. They can be loaded with acrylic materials without the need to hide the grey of the titanium.

Provisional abutments are not for prolonged (more than 180 days) or definitive use.

They are screwed into the implant with the titanium hexagonal screw at 20 Ncm and at 20 Ncm on the MULTI-IM® transepithelial or UNIT® with the titanium screw.

In the laboratory titanium hexagonal screws should be used so as not to damage the definitive gold screw.

The drilling in the laboratory must be performed at low

revolutions. Excessive overheating that could melt the plastic must be avoided. In cases in which it is necessary to use cement on these parts, it is advisable to avoid polyurethane-based cements, as they may present problems of setting.

They are not indicated in patients with bruxism or in situations with a high demands for functional load.

c. Ceramic-On® ceramisable abutment posts

Abutment posts with profile suitable for direct ceramisation for screw-mounting on implants. These enable the restoration of multiple and individual teeth. They are screwed directly into the implant and used with cemented prostheses. The impression is taken directly from the implant (indirect technique).

They are screwed into the implant with the hexagonal titanium screw or with the gold square at 35 Ncm.

They are indicated for multiple and single restorations. When the height of the tissues is 1 mm or greater, it is possible to use angled ceramic abutments. Crowns with a ceramic shoulder or fully ceramic crowns can be cemented onto them.

They are designed to add ceramic for titanium in the titanium abutments. Fusion interval: 1020° - 1155°.

They are made in commercially pure titanium.

This component for the Universal and Universal Plus platforms can also be used in implants from the Wide platform (Interna), provided they are placed juxta- or subcrestally, to improve the biological space.

d. Cylinders for straight and angled MULTI-IM®

Gold cylinders: UCLA-type gold cylinders with plastic casting cylinder for wax-ups. They are cast-on with precious alloys for ceramic, for later ceramisation, or with gold for creating stumps. They are used for the final restoration on MULTI-IM® abutments. Their bevelled design enables a better union between the porcelain and the gold cylinder. We recommend using an alloy with a high gold content for cast-on (CET 25-500°C: 14.6x10-6K-1; 25-600°C: 14.9x10-6K-1). They incorporate a casting chimney, to facilitate the wax-up work in the laboratory. The parts are rotational. If needed, the casting cylinder can be separated from the gold body by traction. They are manufactured in gold alloy. (1400 °C).

e. Transepithelials and accessories DALBO®

Consisting of the male of the Dal-Bo implant retention

for overdentures for implants. They are screwed into the implant at 20 Ncm and the female of the Dal-Bo retention system. This has two parts: The exterior shell in titanium for insertion in acrylic or for welding to a metal structure, and the inside with lamellae that can be adjusted using a screwdriver to graduate the retention of the prosthesis between 200 and 1200g L.

Abutment with two parts, designed to stabilise and retain the partial or complete prosthesis.

Screw in directly to the implant with the handle MDPT and the tip PTDO at a torque of 20 Ncm. It needs a minimum intermaxillar distance (available in different heights). The impression is taken directly from the implant.

They are indicated for partial or complete removable prostheses, with a maximum divergence of 10°. It enables the restoration to emerge even at only 2 mm over the implant.

Made in commercially pure Titanium.

This component for the Universal and Universal Plus platforms can also be used in implants from the Wide platform (Interna), provided they are placed juxta- or subcrestally, to improve the biological space.

The correct seating of the prosthesis on the mucous membrane must be checked at least once a year, if necessary, it must be relined to prevent swinging movements (overloads). We recommend checking the prosthesis at regular intervals of approximately three months and replacing the retention inserts if necessary

f. Screws

Gold screws: They are manufactured in a precious alloy coated with a layer of pure gold measuring 4 microns and with a square slot for the screwdriver. Used to secure the definite restoration with abutment posts, ceramisable posts, gold posts, UCLA abutments, etc. Their design and material make it possible to achieve improved pre-load of the screw. Only available for direct-to-implant prosthetics.

Titanium screws: They are manufactured in grade 5 titanium with a hexagonal slot for the screwdriver. Partial or complete restoration of the teeth in partially or fully edentulate patients, respectively. The titanium screws are used as to secure the prosthesis to the dental implants whether provisionally or definitively, or in the laboratory. Available for direct-to-implant prostheses, angled prostheses and transepithelial restorations.

The following table shows the recommended torques for each type of screw:

TYPE OF SCREW	RECOMMENDED TORQUE (N.cm)
Definitive direct-to-implant screws	35
Definitive screws for angulations	20
Trans-epithelial screws - TTMIR	20
Trans-epithelial screws - TTMIUPA	30
Trans-epithelial screws on DALBO dentures	20

 Apply the torque shown on the label to each screw. If a torque greater than that indicated is applied, it may cause the screw to snap. However, if a lower torque is applied, the prosthetic component may become loose

g. Titanium and Bio-abutment Provisional Cylinders

Titanium cylinders with exterior relief to prevent the acrylic rotating. They are screwed into the implant or transepithelial abutment.

They are multiple/single provisional prostheses in acrylic for screw-mounting on implants or transepithelials.

They have different emergence profiles depending on the implant platform for which they are indicated.

They are screwed into the implant with the gold square screw or the titanium hexagonal screw at 35 Ncm and at 20 Ncm on the MULTI-IM® transepithelial or UNIT® with the titanium screw.

The Bio-abutments® are indicated if aesthetics are required in the anterior sector. They are used when we wish to reduce the transgingival profile, for a horizontal transfer of the biological space in cases of limited interdental space.

h. Titanium straight and angled abutment posts and Bio-abutments®

Straight and angled titanium abutments compatible with the different connections and platforms of BTI implants.

They are the support for multiple/single prostheses for screw-mounting on implants or transepithelials.

They are screwed directly into the implant and used with cemented prostheses.

The impression is taken directly from the implant (indirect technique).

They are screwed into the implant with the square gold or hexagonal titanium screw at 35 Ncm.

In cases juxta-gingival implants, they can be micromilled to the limit (leaving a minimum margin of metal on the vestibular side).

The Bio-abutments® are indicated if aesthetics are required in the anterior sector. They are used when we wish to reduce the transgingival profile, for a horizontal transfer of the biological space in cases of limited interdental space.

This component for the Universal and Universal Plus platforms can also be used it in implants from the Wide platform (Interna), provided they are placed juxta- or subcrestally, to improve the biological space.

i. UCLA gold abutments

Machined gold cylinders used as an alternative in the construction of the conventional UCLA abutment. The non-rotational one is used in the construction of individual teeth (screw-mounted) or stumps (cemented). The rotational one is used in cases of multiple screw-mounted restorations.

The impression is taken directly from the implant (indirect technique).

They are screwed into the implant with the square gold or hexagonal titanium screw at 35 Ncm.

They are indicated for multiple and single restorations. It also solves problems of angle and customisation of the emergence, as angled abutments can be custom made. They are cast-on with precious alloys for ceramic, for subsequent ceramisation, or with gold for manufacturing stumps. (CET 25-500°C: 14.6x10-6K-1; 25-600°C: 14.9x10-6K-1). They are made with a high gold content and a high melting point (1400o).

They are composed of two parts, a gold cylinder and a plastic casting chimney.

This component for the Universal and Universal Plus platforms can also be used it in implants from the Wide platform (Interna), provided they are placed juxta- or subcrestally, to improve the biological space.

9. LABELLING OF THE CONTAINER / SYMBOLS USED

Symbol	Symbol Description	Symbol source
	Catalogue number	ISO 15223-1 (5.1.6)

	Batch code	ISO 15223-1 (5.1.5)
	Manufacturer	ISO 15223-1 (5.1.1)
	Date of manufacture	ISO 15223-1 (5.1.3)
	Medical device	ISO 15223-1 (5.7.7)
	Unique device identifier	ISO 15223-1 (5.7.10)
	CE marking is the manufacturer's declaration that the product meets the requirements of the applicable EC legislation	MDD 93/42/EEC (2007/47/EC) /
	Use-by date	MDR (EU) 2017/745
	Consult instructions for use Follow the link to the eIFU: ifu.bti-biotechnologyinstitute.com/	ISO 15223-1 (5.1.4)
	Caution	ISO 15223-1 (5.4.4)
	Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.	21 CFR 801.109(b)(1)
	Do not re-use	ISO 15223-1 (5.4.2)
	Do not re-sterilize	ISO 15223-1 (5.2.6)
	Sterilized using irradiation	ISO 15223-1 (5.2.4)
	Single sterile barrier system	ISO 15223-1 (5.2.11)
	Do not use if package is damaged and consult instructions for use	ISO 15223-1 (5.2.8)
	MR conditional	FDA Guidance- Testing and Labeling Medical Devices for Safety in the Magnetic Resonance Environment-2021.

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