

BTI Surgical Drills

en	4
es	9
de	14
fr	19
it	24
pt	29
cs	34
ro	39
hu	44
hr	49
sv	54
et	59
fi	64
el	69

BTI Surgical Drills

1. PRODUCT DESCRIPTION

The BTI drill system is used for carrying out surgery, through advanced surgical techniques, related to the placement or extraction of an implant. All BTI drills are made with surgical grade stainless steel. They have a series of depth markings depending on the type of drill and some also include coloured rings for easier identification.

2. LIST OF COMPONENTS AND INDICATIONS FOR USE

FAMILY
INITIAL DRILLS
Indication: Drills with a very active apex and with excellent penetration capacity in the cortical bone. It achieves precise drilling and an exact location of the starting point for the drilling, particularly in narrow crest bones, with good directional control. Designed in two lengths, one short to have better access in posterior areas, another long for cases where the space between teeth conditions the drilling. They also allow lateral drilling should you wish to modify the location of the neo-alveolus by half a millimetre or to considerably alter the angle of the insertion axis.
DIAMETER DRILLS
Indication: With increasing diameters, they are used for widening the socket prior to inserting the implant. They are used progressively, without irrigation, depending on the diameter of the implant to be placed and on the bone quality. They have a very retentive helical design that allows the recollection of drilled bone that is trapped in its threads.
FRONT CUTTING DRILLS
Indication: Designed to achieve maximum settlement for extra-short implants in both maxillae and to work the cortical bone in transalveolar sinus elevations and in proximity to the dental nerve. The apical stem of the drill is inserted in the internal thread of the implant for better guiding to maintain the same axis and optimise function.
PROFILING DRILLS

Indication: BTI profiling drills are profiling tips, as they must be used with the socket wrench (Ref. LLEC / LLEC1) in order to profile the bone around the implant. Manual use is recommended to have greater control and sensitivity.

TREPHINE DRILLS

Indication: Trephine explantation drills are drills designed for the removal of implants or bone blocks for particulate grafting. These drills can work up to a depth of 15 mm, as they have a cylindrical body without a step.

Indication: Trephine drills for implant extraction have different diameters depending on the implant to be extracted. They enable the release of the area of maximum osseointegration (2-3 mm), preserving the surrounding bone and thus reducing the amount of bone eliminated in the extraction. The residual alveolar ridge is respected as much as possible, in addition to reducing the temperatures reached during its use. They have an active limit of 8 mm that acts as a guide for the correct drilling depth. Their use is indicated for cases in which extraction with the reduced protocol with the implant extractors involves jumping the 200 Ncm spanner on three consecutive occasions. The particularly fine profile of these drills minimises the bone loss of the resulting socket after implant extraction.

Indication: The trephine drills FTEX27-6 and FTEX32-6 have a smaller diameter in the first millimetres and are suitable for biopsies, for the removal of Ø2.5mm and Ø3mm expander implants and for the removal of Tiny Ø2.5mm and Ø3mm implants that have broken in the cylindrical area.

TREPHINE DRILL FOR BIOPSIES

Indication: FT2-6 is indicated for biopsies and in the first millimetres it has a more reduced diameter.

COUNTERSINK DRILLS

Indication: Inverted cone-shaped drills that allow adapting the bone cortical to the implant neck, avoiding compression and unwanted resorptions at that level, particularly in very dense cortical bones.

DISPOSABLE DRILL KITS

Indication: The kit is designed for use during the usual drilling protocol in implantology, for both internal and external connections, and for all platforms in the implants with a smaller diameter. The kit is sold in the form of single-use sterile and blister packed drills for use with one patient only.

EXTENSION DRILL

Indication: This adapts to the drills, lengthening them, facilitating drilling in areas where access is hindered by long adjacent teeth.

3. INSTRUCTIONS FOR USE

All of the BTI drills except for bone profiler drills, are designed to be connected to a surgical motor; you must ensure the drill is properly connected to the unit before use, by inserting it as far as possible into the handle until it is clipped in.

The extension drill is used in the same way by connecting the piece to the surgical motor up to clipping and the drill to the extension drill up to clipping.

Recommendations for use:

	Make sure the handle is in perfect technical and hygienic conditions.
	You should check the rotation of the drill before applying it to the work surface.
	Prevent heat damage caused by the rotation at all times; apply a low revolution count and irrigate thoroughly; insufficient irrigation at high revolutions may cause irreversible harm to the adjacent tissue.
	Inappropriate use results in a risk and can reduce the efficacy of the BTI drills; observe the recommendations for use and the speeds defined for each type of drill.
	Avoid an excessive working pressure, that can cause damage to the ridges of the instruments, increase the temperature and in extreme cases, break the instrument.

The following table shows the drilling speed and irrigation conditions recommended for each type of drill.

Drills	No. of Revolutions	Irrigation
Initial drill	800-1000 rpm	Yes
Diameter drills	50-75 rpm	No
Countersink drill	150-200 rpm	No
Front cutting drill	50-75 rpm	No
Trephine drill	500 rpm	Yes
Bone profiler	Manual	No

4. CONTRAINDICATIONS

Do not use in patients with a pre-existing and known allergy to stainless steel. The general contraindications, warnings and precautions in dental implantology must be observed before using this product.

In the case of the extension drill, there are no contraindications related to the device.

5. WARNINGS

	Do not apply more than 80 Ncm with the extension drill.
	The use of these products does not entail additional precautions in children, pregnant women and lactating women to those already contemplated in general implantology.
	The use of drills and other components of the surgical box is restricted to odontologists, stomatologists and maxillofacial surgeons.
	A lack of hygiene or cooperation by the patient, and systemic diseases (diabetes, smoking, etc.) are potential causes that can hinder the subsequent recovery from surgical intervention.
	Reusable products should be cleaned and sterilised before use to avoid risks of infection and cross-contamination.
	In case of ingestion of the piece by the patient, refer the patient to the hospital emergency department for appropriate treatment.
	The direction of rotation of these products is clockwise (normal surgical motor mode).

6. PRECAUTIONS**• EXPIRY AND NUMBER OF USES**

KFDx references (DISPOSABLE DRILL KIT) are products that are marketed sterile and have a shelf life of 5 years from the date of manufacture as long as the storage and preservation conditions of the blister are maintained.

 KFDx: The products contained in these references are for a single use (for one single patient); they must not, under any circumstances, be reused or re-sterilised. Their reuse can lead to the risk of infection and cross contamination.

For reusable products (THE REST OF REFERENCES), BTI recommends a maximum of 13 uses. This is simply a guide, as the true useful life may be different depending on the type of application and/or characteristics of material you are working with (hardness of the bone or mould).

In any event it is recommended that instruments with damaged, deformed and bent blades are disposed of, as they cause

vibrations and can produce imperfections in the preparation margins, as well as irregular surfaces.

It is the user's responsibility to examine the drills and check that they are in an appropriate condition for their intended use. If the instrument shows no signs of wear it may be used.

 Reusable products should be cleaned and sterilised before use in accordance with point 8 to avoid risks of infection and cross-contamination.

7. ADVERSE EFFECTS

The use of BTI Drills as such does not lead to the emergence of adverse effects. Nevertheless, the use of drills is directly linked to dental implantology techniques, so the adverse effects that this can cause.

8. CLEANING AND STERILISATION PROCEDURE

Recommendations for cleaning and sterilization procedure:

	Use protecting gloves
	Before they are used for the first time they must be sterilised, except for the KFDx and KFEx.
	Clean the instrumental after each surgical process
	IMPORTANT: do not clean different materials at the same time (chrome cover, aluminium, stainless steel).
	The autoclave cycle should be allowed to finish, without interrupting the drying stage to avoid the risk of oxidation.
	Make sure that no oxidated instrument or material is sterilised, as it increases the risk of the rest of the instruments becoming oxidised when sterilised together

STAGES OF THE PROCEDURE

PRE-TREATMENT

After each surgical process, every surgical instrument should be placed into a container that will have a cleaning solution (0.5% neodisher MediZym enzymatic, Dr. Weigert; indications for use: 5 mL in 1 L of water) at room temperature, in order to avoid that blood or bodily fluids could dry before decontaminating. Keep the instrument in this way the shortest possible time (within a maximum of 15 minutes) before cleaning and disinfecting.

Immediately after this pre-treatment, they must be manually or automatically cleaned and disinfected, then sterilisation in an autoclave.

MANUAL CLEANING AND DISINFECTION

1. Remove the majority of the impurities with a soft-bristle brush using plenty of running water at ambient temperature for at least 3 minutes, paying attention on areas where access is difficult such as seals, angles, blind/narrow holes, etc. as well as rough surface or containing laser engraving.
2. Place instruments in an ultrasonic bath (Bandelin, HF power 300 Weff, frequency 35 kHz) with purified water which contains a solution of neutral detergent with enzymatic capacity and no foam development, (0.5% neodisher MediZym enzymatic, Dr. Weigert; indications for use: 5 mL in 1 L of water).

Be sure that detergent is suitable for use with instruments and ultrasonic bath

3. Keep instruments on the ultrasonic bath for 20 minutes at 35°C. Be sure that immersion time does not start until the last instrument is placed on the bath and it is completely immersed. The instruments must not have contact among them to prevent possible physical damage. Do not exceed temperatures over 45°C (risk of coagulation proteins). The solution should be replaced as soon as turbidity or suspended particles being detected due to cleaning activity could be inhibited.

 **IMPORTANT:** Cleaning instruments by ultrasonic bath could not be suitable for instruments that contain carbon of tungsten in their composition. It is also not recommended for cleaning of rubber or plastic items as these tend to absorb ultrasounds.

4. Remove instruments from the ultrasonic bath and brush them again with a soft bristle nylon brush rinsing with tap water at room temperature for at least 1 min/piece, paying attention on areas where access is difficult such as seals, angles, blind/narrow holes, etc. as well as rough surface or containing laser engraving
5. Rinse 5 times with at least 10 mL each of deionized water at room temperature using a single-use syringe.
6. Perform a visual inspection of instruments paying special attention on areas where access is difficult such as seals, angles, blind holes, etc. If trace of impurities were observed, cleaning operation from point 1 to 6 must be repeated. Be sure that every instrument is sufficiently covered and there is no contact among them.
7. Once clean, place instruments in disinfection container which contains disinfectant solution at room temperature (Cidex OPA Solution; indications for use: soak instruments for 15 minutes at room temperature).
8. Ensure all the instruments are submerged and are not touching each other.

- Remove instruments from the disinfectant container and rinse them soaking into 8 L of distilled water for 1 minute; rinse all lumens 5 times with at least 20 mL each of distilled water at room temperature using a single-use syringe. Repeat this process two more times to get three independent rinses.

 **IMPORTANT:** If Sterilization process could not be possible or recommended after disinfection, be sure to use sterilized water for rinsing in order to avoid recontamination.

- Dry the instruments at room temperature with filtered compressed air until it is completely dry.

AUTOMATIC CLEANING AND DISINFECTION

 Use always a washer-disinfector certified in accordance with ISO 15883 (for example: Miele G 7836 CD; mobile injector unit E 450 for MIC instruments (Miele); accessory parts E 142 mesh tray with lid (Miele) or E 451 mesh tray with lid for small parts (Miele)).

- Place the instruments in the washer-disinfector according to the manufacturer's instructions. Load instruments such that holes can drain.
- Add a neutral detergent with enzymatic power - according to the manufacturer's instructions -- and select a suitable programme for the cleaning and final rinse according to the instructions of the washer/disinfector manufacturer.
- After cleaning and rinsing, proceed to thermal disinfection: 5 minutes at 90°C.
- Dry the instruments at room temperature with filtered compressed air until it is completely dry.

INSPECTION

Check instruments and ensure that there is no damage, wear or rust, and all of them work properly.

STERILISATION IN AUTOCLAVE

- Introduce instruments, immediately after the inspection, into sterilization pouches with sterilization indicator and seal the bag as specified by the supplier. Ensure that the pack is large enough to contain the instrument without stressing the packaging.
- If sterilization container is used (for example autoclavable BTI surgical kit box) introduce it into a sterilization bag of sufficient size without stressing the packaging.

 **IMPORTANT:** Only use wrapping materials and systems (sterile barrier systems) that comply with the DIN EN ISO 11607-1 standard. The correct application of suitable

wrapping is critical for the success of the sterilization.

- Introduce sterilization pouches into a pre-vacuum autoclave (for example: CELITRON MEDICAL TECHNOLOGIES, STING 11) with paper up or sideways. Be sure not to exceed the maximum load of the autoclave specified by the supplier and not to stack up bags.
- Sterilize at 134°C (pressure under autoclave supplier specification, for example: 2 bar) for 6 minutes. The changing colour of the indicator in the bag shows if the sterilization process has been carried out correctly.
- Be sure the sterilization requirements (pressure, temperature and time) are the ones described in this instruction for use.
- After completing the sterilization cycle, let dry and cool bags inside the autoclave for at least 20 minutes, before removing and storing them.

 **IMPORTANT:** Only use steam sterilizers that comply with the EN 13060 standard. Be sure that sterilization program is established according to the temperature, pressure and time indicated. In case of doubt, contact the steam sterilizer supplier.

INSPECTION

Check instruments and ensure that there is no damage, wear or rust, and all of them work properly.

STORAGE

Place the sealed bags in a closed and dry place. Make sure there is enough space between the instruments and the floor, walls and ceiling for circulate enough air.

The maximum storage time of the product inside sterilisation pouches shall be as indicated by the manufacturer of the sterilisation pouches.

DISPOSAL

At the end of their shelf life, the drills should be disposed of in the biological waste container.

9. LABELLING OF THE CONTAINER / SYMBOLS USED

	Batch code
	Catalogue number
	Do not re-use

	Use-by date
	Sterile using irradiation
	CE marking
Rx only	Caution: Federal law restricts this device to sale by or on the order of a health professional.
	Manufacturer
	Do not use if package is damaged and consult instructions for use
	Do not resterilize
	Caution
	Consult instructions for use or consult electronic instructions for use

NOTE: The devices should be accepted only if the factory packaging and labeling arrive intact.

Contact your distributor if the package is opened or altered.

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