

# BTI Surgical Drills

|           |    |
|-----------|----|
| <b>en</b> | 4  |
| <b>es</b> | 7  |
| <b>de</b> | 10 |
| <b>fr</b> | 13 |
| <b>it</b> | 16 |
| <b>pt</b> | 19 |
| <b>cs</b> | 22 |
| <b>lt</b> | 25 |
| <b>hu</b> | 28 |
| <b>nb</b> | 31 |
| <b>sv</b> | 34 |
| <b>da</b> | 37 |
| <b>fi</b> | 40 |
| <b>el</b> | 43 |
| <b>bg</b> | 46 |

# BTI Surgical Drills

## 1. PRODUCT DESCRIPTION

The BTI drill system is used for carrying out surgery, through advanced surgical techniques, related to the placement or extraction of an implant. All BTI drills are made with surgical grade stainless steel. They have a series of depth markings depending on the type of drill and some also include coloured rings for easier identification.

## 2. LIST OF COMPONENTS AND INDICATIONS FOR USE

| FAMILY                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>INITIAL DRILLS</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Indication:</b> Drills with a very active apex and with excellent penetration capacity in the cortical bone. It achieves precise drilling and an exact location of the starting point for the drilling, particularly in narrow crest bones, with good directional control. Designed in two lengths, one short to have better access in posterior areas, another long for cases where the space between teeth conditions the drilling. They also allow lateral drilling should you wish to modify the location of the neo-alveolus by half a millimetre or to considerably alter the angle of the insertion axis. |
| <b>DIAMETER DRILLS</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Indication:</b> With increasing diameters, they are used for widening the socket prior to inserting the implant. They are used progressively, without irrigation, depending on the diameter of the implant to be placed and on the bone quality. They have a very retentive helical design that allows the recollection of drilled bone that is trapped in its threads.                                                                                                                                                                                                                                          |
| <b>FRONT CUTTING DRILLS</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Indication for FRONT CUTTING DRILLS FOR SHORT IMPLANTS:</b> Indicated for flattening the bone to achieve maximum settlement for extra-short implants and to work the cortical bone in transalveolar sinus elevations and in proximity to the dental nerve. They are also used to collect bone during drilling. These drills should be used when the appropriate socket for their diameter has been prepared, thus making an apical instrumentation at the base of the sinus and/or                                                                                                                               |

allowing further advancement of the implant, avoiding apical compression.

**Indication for FLATTENING FRONT CUTTING DRILLS:** Indicated for flattening the bone in irregular ridges prior to implant placement. They are also used to collect bone during drilling. Before using these drills, in order to flatten the ridge, the 1.8 mm diameter drill must be inserted 3 mm deeper than the height to be reduced. This allows the drill to be guided through the existing protrusion. Recommended to be used without irrigation.

### PROFILING DRILLS

**Indication:** BTI profiling drills are profiling tips, as they must be used with the socket wrench (Ref. LLEC / LLEC1) in order to profile the bone around the implant. Manual use is recommended to have greater control and sensitivity.

### TREPHINE DRILLS

**Indication for EXPLANTATION TREPHINE DRILLS (FTEXxx):** Drills indicated for the extraction of implants or bone blocks for particulate grafting.

**Indication for TREPHINE DRILLS FOR IMPLANT EXTRACTION (FTEXxxx-8):** Drills indicated for the extraction of implants when the extraction of the implant with the extractor fails 3 consecutive times at 200Ncm. They have an active limit of 8 mm that acts as a guide for the correct drilling depth.

**Indication for TREPHINE DRILLS (FTEXxxx-6):** Drills indicated for the extraction of implants or bone blocks for particulate grafting. The first 6 mm of the drill have a smaller diameter.

**Indication for TREPHINE DRILL (FTxx-6):** Drills indicated for the extraction of bone blocks. The first 6 mm of the drill have a smaller diameter.

**Indication for PUNCH DRILL (FPxx):** Punch drill to cut the gingiva and place the implant without raising the flap.

### COUNTERSINK DRILLS

**Indication:** Inverted cone-shaped drills that allow adapting the bone cortical to the implant neck, avoiding compression and unwanted resorptions at that level, particularly in very dense cortical bones.

**DISPOSABLE DRILL KITS**

**Indication:** The kit is designed for use during the usual drilling protocol in implantology, for both internal and external connections, and for all platforms in the implants with a smaller diameter. The kit is sold in the form of single-use sterile and blister packed drills for use with one patient only.

**EXTENSION DRILL**

**Indication:** This adapts to the drills, lengthening them, facilitating drilling in areas where access is hindered by long adjacent teeth.

**3. INSTRUCTIONS FOR USE**

All of the BTI drills except for bone profiler drills, are designed to be connected to a surgical motor; you must ensure the drill is properly connected to the unit before use, by inserting it as far as possible into the handle until it is clipped in.

The extension drill is used in the same way by connecting the piece to the surgical motor up to clipping and the drill to the extension drill up to clipping.

**Recommendations for use:**

|                                                                                     |                                                                                                                                                                                                                |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Make sure the handle is in perfect technical and hygienic conditions.                                                                                                                                          |
|                                                                                     | You should check the rotation of the drill before applying it to the work surface.                                                                                                                             |
|                                                                                     | Prevent heat damage caused by the rotation at all times; apply a low revolution count and irrigate thoroughly; insufficient irrigation at high revolutions may cause irreversible harm to the adjacent tissue. |
|                                                                                     | Inappropriate use results in a risk and can reduce the efficacy of the BTI drills; observe the recommendations for use and the speeds defined for each type of drill.                                          |
|                                                                                     | Avoid an excessive working pressure, that can cause damage to the ridges of the instruments, increase the temperature and in extreme cases, break the instrument.                                              |
|                                                                                     |                                                                                                                                                                                                                |

The following table shows the drilling speed and irrigation conditions recommended for each type of drill.

| Drills              | No. of Revolutions | Irrigation |
|---------------------|--------------------|------------|
| Initial drill       | 800-1000 rpm       | Yes        |
| Diameter drills     | 50-75 rpm          | No         |
| Countersink drill   | 150-200 rpm        | No         |
| Front cutting drill | 50-75 rpm          | No         |
| Trephine drill      | 500 rpm            | Yes        |
| Profiling drill     | Manual             | No         |
| Punch drill         | 50-75 rpm          | No         |

**4. CONTRAINDICATIONS**

Do not use in patients with a pre-existing and known allergy to stainless steel. The general contraindications, warnings and precautions in dental implantology must be observed before using this product.

In the case of the extension drill, there are no contraindications related to the device.

**5. WARNINGS**

|                                                                                     |                                                                                                                                                                                           |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Do not apply more than 80 Ncm with the extension drill.                                                                                                                                   |
|                                                                                     | The use of these products does not entail additional precautions in children, pregnant women and lactating women to those already contemplated in general implantology.                   |
|                                                                                     | The use of drills and other components of the surgical box is restricted to odontologists, stomatologists and maxillofacial surgeons.                                                     |
|                                                                                     | A lack of hygiene or cooperation by the patient, and systemic diseases (diabetes, smoking, etc.) are potential causes that can hinder the subsequent recovery from surgical intervention. |
|                                                                                     | Reusable products should be cleaned and sterilised before use to avoid risks of infection and cross-contamination.                                                                        |
|                                                                                     | In case of ingestion of the piece by the patient, refer the patient to the hospital emergency department for appropriate treatment.                                                       |
|                                                                                     | The direction of rotation of these products is clockwise (normal surgical motor mode).                                                                                                    |

**6. PRECAUTIONS****• EXPIRY AND NUMBER OF USES**

KFDx references (DISPOSABLE DRILL KIT) are products that are marketed sterile and have a shelf life of 5 years from the date of manufacture as long as the storage and preservation conditions of the blister are maintained.

 KFDx: The products contained in these references are for a single use (for one single patient); they must not, under any circumstances, be reused or re-sterilised. Their reuse can lead to the risk of infection and cross contamination.

For reusable products (THE REST OF REFERENCES), BTI recommends a maximum of 13 uses. This is simply a guide, as the true useful life may be different depending on the type of application and/or characteristics of material you are working with (hardness of the bone or mould).

In any event it is recommended that instruments with damaged, deformed and bent blades are disposed of, as they cause vibrations and can produce imperfections in the preparation

margins, as well as irregular surfaces.

It is the user's responsibility to examine the drills and check that they are in an appropriate condition for their intended use. If the instrument shows no signs of wear it may be used.

Reusable products should be cleaned and sterilised before use in accordance with point 8 to avoid risks of infection and cross-contamination.

NOTE: The devices should be accepted only if the factory packaging and labelling arrive intact. Contact your distributor if the package is opened or altered.

## 7. ADVERSE EFFECTS

The use of BTI Drills as such does not lead to the emergence of adverse effects. Nevertheless, the use of drills is directly linked to dental implantology techniques, so the adverse effects that this can cause.

## 8. CLEANING AND STERILISATION PROCEDURE

BTI surgical drills and Dental drill extender are reusable products. They are supplied non-sterile and must be reprocessed before use.

BTI only recommends the reprocessing protocol described in the CAT246 cleaning, disinfection, and sterilisation guide.

## 9. LABELLING OF THE CONTAINER / SYMBOLS USED

For a description of the symbols that appear on product labels and in these instructions, refer to the MA087 guide.

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