
BTI DENTAL IMPLANTS

bg	Български	ЗЪБНИ ИМПЛАНТИ BTI _____	(Страница 3)
cs	Čeština	ZUBNÍ IMPLANTÁTY BTI _____	(Stránka 7)
da	Dansk	BTI-TANDIMPLANTATER _____	(Side 11)
de	Deutsch	BTI-ZAHNIMPLANTATE _____	(Seite 15)
el	Ελληνικά	ΟΔΟΝΤΙΚΑ ΕΜΦΥΤΕΥΜΑΤΑ BTI _____	(Σελίδα 19)
en	English	BTI DENTAL IMPLANTS _____	(Page 23)
es	Español	IMPLANTES DENTALES BTI _____	(Página 27)
fi	Suomen	BTI-HAMMASIMPLANTIT _____	(Sivu 31)
fr	Français	IMPLANTS DENTAIRE BTI _____	(Page 35)
hu	Magyar	BTI FOGIMPLANTÁTUMOK _____	(Oldal 39)
it	Italiano	IMPIANTI DENTALI BTI _____	(Pagina 43)
lt	Lietuvių	BTI DANTŲ IMPLANTAI _____	(Puslapis 47)
nb	Norsk Bokmål	BTI-TANNIMPLANTATER _____	(Side 51)
pt	Português	IMPLANTES DENTÁRIOS BTI _____	(Página 55)
sv	Svenska	BTI-TANDIMPLANTAT _____	(Sida 59)

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BTI Dental Implants

1. PRODUCT DESCRIPTION

BTI dental implants are endosseous implants made of pure grade 4 titanium, available in different platforms, diameters and lengths.

Their surface is obtained through selective acid treatment. The UNICCA surface, chemically modified with calcium and with triple roughness, protects the surface and makes the implants hydrophilic until they are placed. Under normal temperature and relative humidity conditions, it has a moist appearance.

BTI dental implants are available with two types of connection: INTERNA and EXTERNA. The EXTERNA connection consists of a hexagon whose size varies depending on the platform, while the INTERNA connection can have a tetra-, hexa- or octolobular design, depending on the corresponding platform.

The implants are supplied sterile in a screw-cap container, which includes the implant holder or transporter and the cover screw located at the bottom of the inner support.

MATERIALS

Component	Material	Chemical composition	Composition % (mass/mass)
BTI dental implant	Grade 4 titanium	Total residues (Fe, O, C, N, H)	≤1.05
		Titanium (Ti)	Balanced
	Coating	CaCl ₂	0.07 ≤ x ≤ 0.7 (*)
Implant holder Closing screw	Grade 4 titanium	Titanium (Ti)	≤1.05
		Total residues (Fe, O, C, N, H)	Balanced

*Depending on the geometry of the implant.

ACCESSORIES

EXTENSORS (CPIXX and CAIXXX): accessories for the handpiece that facilitate the transport and insertion of implants and other instruments. They are made of surgical grade stainless steel and designed for connection to a surgical motor. Detailed information on their use can be found in the corresponding Instructions for Use (MA290).

COMPATIBLE DEVICES

The following devices can be used in combination with BTI dental implants, and their Instructions for Use should be consulted for general warnings and detailed information.

Compatible device	Associated IFU
Extenders and Keys	CAT246
BTI surgical drills	MA140
MULTI-IM and UNIT ABUTMENTS	MA159
Abutments with Sirona technology	MA161
SQUARE BTI Abutments and Interfaces	MA211
Prosthetic components	MA137
Abutments, screws and healing caps	MA077
KEXIM implant removal kit	MA061

2. INTENDED USE

The BTI dental implant system, consisting of implants, prosthetic components and accessory tools, is used in oral implantology to provide a support structure for the partial or total replacement of teeth in edentulous patients.

3. INDICATIONS FOR USE

Dental implants are indicated for single, partial or complete prosthetic restorations of the upper or lower jaw. Once attached to the bone, the implants act as anchors for various fixed or removable prosthetic solutions that can be used to improve or restore the patient's chewing function. BTI dental implants can be used under the following conditions:

2.5 mm Diameter Implants on platforms: EXTERNA TINY, INTERNA 3.0 and INTERNA 2.5

Indicated for use with splints in multiple restorations screwed to the specific MULTI-IM abutment. Not to be used for single restorations or multiple prostheses screwed directly onto the implant.

3.0 mm Diameter Implants on platforms: EXTERNA TINY and INTERNA 3.0

Indicated for use with splints in multiple restorations screwed to the MULTI-IM abutment and prostheses screwed directly onto the implant. They can be used for single restorations in lower or upper lateral incisors with no occlusion.

3.3 and 3.5 mm Diameter Implants on platforms: INTERNA 3.0

Indicated for splinted use in multiple restorations screwed to the MULTI-IM abutment and prostheses screwed directly onto the implant. They can be used for single restorations in upper lateral incisors with or without occlusion and lower incisors with no occlusion.

4.5 mm length implants on all platforms.

Indicated for use always splinted in multiple restorations screwed to the MULTI-IM abutment. Do not use for single restorations, multiple prostheses screwed directly onto the implant, or for angled or removable prostheses. Recommended for partial or complete prosthetic restorations with severe atrophy of the

upper or lower jaw. In the case of the lower jaw, they require a residual bone of at least 5 mm above the dental nerve.

4. INTENDED USER AND INTENDED PATIENT GROUP

Intended user: Dental healthcare professionals in hospital or clinic settings with specific training in oral implantology and in BTI dental products.

Intended patient group: Partially or totally edentulous male and female adults who do not have one or more conditions or pathologies listed in the contraindications.

5. CONTRAINDICATIONS

Absolute: Situations in which the product should not be used under any circumstances because the risk is too high:

- Severe immunosuppression, whether due to disease or drug treatment.
- Coagulation disorders.
- Drug abuse.
- Ongoing cancer or undergoing antineoplastic treatment.
- Pregnant or breastfeeding women, as implantation usually requires X-rays.
- Bone crest is not high or wide enough to ensure initial fixation and stability of the implant.
- Known allergy or hypersensitivity to commercially pure titanium.
- Children who have not completed bone development.

Relative: Situations in which the use of the product may be possible under certain conditions but with caution and medical supervision. The benefit must always be weighed against the risk:

- Uncontrolled diabetes; in controlled cases, individual assessment is required.
- Autoimmune diseases.
- Patients with high alcohol or tobacco consumption.
- Treatment with bisphosphonates, antidepressants, proton pump inhibitors (PPIs) or other drugs that affect bone metabolism.
- Poor oral hygiene.
- Lack of patient cooperation or failure to comply with medical instructions.

6. GENERAL WARNINGS

- BTI dental implants are intended for single clinical use. They are supplied sterile and must not be reused or re-sterilised, as this carries a risk of infection or cross-contamination.
- Do not use after the expiry date indicated.
- Use only original BTI components; parts from other manufacturers may damage the device or cause malfunctions.
- The failure of prostheses, implants or loss of surrounding bone can result from inadequate surgical technique or incorrect patient selection, negatively affecting the osseointegration process and treatment stability. Factors that can compromise recovery include: insufficient bone quality or quantity, infections, periodontitis, bruxism, poor oral hygiene, poor patient cooperation, drug abuse, alcohol dependence, smoking, medications (antidepressants, antiresorptive therapies, anticoagulants, proton pump inhibitors), mental disorders, coagulation disorders, chemotherapy or radiotherapy treatments, immunosuppressed patients (due to immune deficiencies or immunosuppressive treatments), and systemic diseases (such as diabetes or hypertension).
- Dental nerve injuries can cause anaesthesia, paraesthesia, or dysesthesia.
- In case of accidental ingestion of the piece by the patient, go immediately to a hospital emergency department.

7. PRECAUTIONS

- Perform a complete clinical examination of the patient. It is recommended to supplement the evaluation with panoramic and periapical radiographs and cone beam computed tomography (CBCT) to assess anatomical characteristics, periodontal status, and bone suitability. All candidates should be carefully evaluated prior to surgery.
- Handle all biological samples and sharp objects in strict accordance with the centre's biosafety policies.

8. ADVERSE EFFECTS

The implantation procedure carries certain risks that may be temporary or, in exceptional cases, the loss of sensitivity may be permanent. Possible adverse effects include inflammation of the implant area, local pain, bruising, bleeding, nosebleeds, loose sutures, trismus, sinusitis, peri-implantitis, periodontitis, gingivitis, mucositis, anaesthesia, paresthesia, dysesthesia, fistulas, risk of swallowing or aspiration of pieces during the procedure, temporary sensitivity in the area, loss of sensation in the lower lip or chin, lack of osseointegration of the implant, loss of surrounding bone, bone resorption of the maxillary or mandibular ridge, and implant fracture due to fatigue.

Infection may also occur around the implant, although this can usually be controlled with local treatment.

All materials used in BTI implants are biocompatible. However, in rare cases, contact with titanium may trigger hypersensitivity reactions, such as local or generalised eczema, pruritus, pain, inflammation, erythema, gingival hyperplasia or peri-implant granulomatous reaction. Exceptionally, yellow nail syndrome may occur, characterised by thickened, yellowish, slow-growing nails with lymphoedema, possibly related to the release of titanium ions due to corrosion. These reactions are very rare compared to other metals and usually resolve when the metal component is removed.

If you notice any adverse effects, please notify the manufacturer at the following email address: technicalmanager@bti-implant.es. By reporting adverse effects, you can help provide more information about the safety of this medical device.

9. MAGNETIC RESONANCE IMAGING (MRI) COMPATIBILITY

The radiofrequency (RF) safety of the BTI Implant System, which includes implants and prosthetic components, has not been tested. However, a non-clinical review conducted under worst-case conditions has determined that the materials used are safe in terms of the displacement force and torsion induced by MRI for a given spatial gradient.

Therefore, a patient with this device can be safely scanned in an MRI system provided that the following conditions are met and images are taken by marking a reference point at a minimum distance of 30 cm from the implant or ensuring that the implant is outside the range of the RF coil.

The specific conditions for safe scanning are:

- Static magnetic field strength (B_0) up to 3.0 T.
- Maximum spatial field gradient of 30 T/m (3,000 gauss/cm).
- RF excitation with circular polarisation (CP).
- Use of a body transmission coil only if a reference point is marked at least 30 cm from the implant or if the implant is outside the coil.
- Use of limb T/R coils permitted; use of head T/R coil excluded.
- The examination must be performed in normal operating mode within the permitted imaging zone.
- The maximum SAR for the whole body must not exceed 2 W/kg.
- Safety in relation to SAR in the head has not been evaluated, therefore cranial scanning is not recommended.
- There are no specific limitations on the duration of the scan due to heating of

the implants.

10. INSTRUCTIONS FOR USE

WARNINGS

- Visually inspect the implant before use to verify that the sterile barrier system is not damaged (scratches, breaks, or deformities). Do not use the implant if the plastic seal is broken.
- Do not exceed the insertion torques recommended by BTI, as this could cause bone necrosis.
- The use of trans-epithelial abutments other than those indicated here may compromise the stability of the system and the safety of the treatment:
 - INTERNA PLATFORM 3.0 implants (Ø2.5 mm): use only MULTI-IM trans-epithelial abutments, ref. INTMIP32535xx and INTMIP325xx.
 - INTERNA PLATFORM 2.5 implants (Ø2.5 mm): use only MULTI-IM trans-epithelial abutments, ref. INTMIP241xx and INTMIP235xx.
 - 4.5 mm length implants: use only MULTI-IM trans-epithelial abutments, ref. INTMIPE45xx, INTMIPU45xx, INTMIP345xx, INTMIPEEX45xx and INTMI-PU45xx.

PRECAUTIONS

- Be careful when removing the implant from the inner holder: always pull sideways, never from the front of the upper end, to avoid damaging the surface.
- The implant holder indicates the position of the implant. It is recommended to align the vestibular side with the implant connection.
- If the tightening torque of the prosthetic component is greater than the insertion torque of the implant, delayed loading must be performed.
- Tightening torque of cover screws: 10 Ncm for TINY implants and 20 Ncm for all others.

PROCEDURE:

1. Plan the surgery considering the patient's scan and medical history, and select the most suitable implant.
2. The implant is supplied in a rigid vial with a tamper-evident plastic seal and three labels for recording in the patient's medical records. Check the integrity of the seal before use to ensure that the sterile packaging has not been opened or damaged.
3. Complete the patient details and follow the recommended surgical protocol. For the drilling sequence, refer to:
 - CAT221: drilling sequences for BTI EXTERNA implants.
 - CAT218: drilling sequences for BTI INTERNAL implants.
4. Remove the plastic seal and open the vial cap by turning it anticlockwise until the safety ring comes off.

5. Pour the contents of the vial onto a sterile tray. Hold the central part of the inner holder and insert the implant holder connector (CPIX) into the contra-angle hand piece. Slide your fingers towards the bottom of the holder, at the level

of the cover screw, pressing to remove the implant.

6. Once the implant is attached to the implant holder, transfer it face up to the prepared alveolus.
7. Insert the implant with a surgical motor at 25 rpm, with a maximum torque of 25 Ncm.
8. Complete the insertion manually:
 - For all implants EXCEPT INTERNA PLATFORM 2.5:
 - Use the multi-torque surgical key (LLMQ) without exceeding 70 Ncm.
 - Remove the implant holder by holding the body with the open key (LLAV).
 - For INTERNA PLATFORM 2.5 implants:
 - Remove the implant holder manually with the open-ended wrench (LLAV).
 - Complete the insertion with the multi-torque wrench (LLMQ) connected to the implant extender (EAI2XX), without exceeding 60 Ncm.
- Immediate loading: BTI recommends an insertion torque of 40–50 Ncm to ensure primary stability.
9. If the implant is to be buried under the gum, insert the cover screw: Remove the closure screw from the inner support with the hexagonal conical tip, ensuring that it is securely fastened before removal.
 - For headless cover screws (supplied with TINY implants), use a screwdriver or small hexagonal torque tip.
 - For all other cover screws, use a large hexagonal screwdriver.

11. TRANSPORT AND STORAGE

- Caution: Do not expose the product to direct sunlight.
- Warning:
 - Keep the product in its original packaging until ready for use; do not handle it beforehand.
 - Inspect the container before use; do not use if damaged.
 - If the container is open or has been tampered with, contact your distributor.

12. WASTE DISPOSAL

- Dispose of all used material as biological waste.
- Wear gloves during the process.
- Follow current regulations on infectious waste.

13. NOTIFICATION OF SERIOUS INCIDENTS

In accordance with Regulation (EU) 2017/745 (MDR), si if a serious incident occurs during the use of this product or as a result of its use, the user or patient must immediately notify the following address:
technicalmanager@bti-implant.es.

14. ADDITIONAL INFORMATION

The surface of BTI dental implants may show white calcium deposits or appear less moist due to exposure to extreme temperatures during transport or storage. These signs do not affect the proper functioning or safety of the product.

BTI provides regular training courses to optimise product performance.

The Summary of Safety and Clinical Performance (SSCP) for this device is available in the European EUDAMED database (<https://ec.europa.eu/tools/eudamed>). To consult it, you must use the unique device identifier (UDI-DI) indicated on the label. If you are unable to locate it, you can request a copy directly from the manufacturer or, where applicable, through its authorised distributor.

The UDI-DI code corresponds to a series of numeric or alphanumeric charac-

ters that enable the traceability of the medical device and is displayed on the label in ICAD (Identification and Automatic Data Capture) format and in HRI (human-readable interpretation) format.

15. SYMBOLS

For a description of the symbols that appear on product labels and in these instructions, please refer to guide MA087.

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