
# BTI DENTAL IMPLANTS

|           |    |
|-----------|----|
| <b>en</b> | 4  |
| <b>es</b> | 6  |
| <b>de</b> | 8  |
| <b>fr</b> | 10 |
| <b>it</b> | 12 |
| <b>pt</b> | 14 |
| <b>cs</b> | 16 |
| <b>ro</b> | 18 |
| <b>lv</b> | 20 |
| <b>lt</b> | 22 |
| <b>da</b> | 24 |
| <b>hu</b> | 26 |
| <b>sv</b> | 28 |
| <b>et</b> | 30 |
| <b>fi</b> | 32 |
| <b>el</b> | 34 |
| <b>no</b> | 36 |

# BTI DENTAL IMPLANTS

## 1. PRODUCT DESCRIPTION

BTI dental implants are endosseous implants made of pure grade 4 titanium with different platforms, diameters and lengths that have a surface obtained through selective acid-etching.

The BTI UNICCA® surface has a chemical modification with calcium, on triple roughness, that protects the surface and makes the implants hydrophilic until they are placed. The visual appearance of the UNICCA® surface is wet when in normal conditions of temperature and relative humidity.

BTI dental implants are available with two different connection types: INTERNA and EXTERNA and with a wide variety of diameters and lengths to provide solutions for the many different clinical cases that can arise in clinical practice. The EXTERNA connection is a hexagon that changes dimensions depending on the implant platform and the INTERNA connection has a tetra or hexalobular design, depending on the platform.

The BTI dental implants are sterile, they are supplied in a screw-cap container that also includes the implant mount or transporter and the corresponding cover screw in the lower part of the inner holder.

## 2. INTENDED USE

The BTI Dental Implant system for oral implant surgery is to be used for the partial or total replacement of teeth in edentulous patients. Once osseointegrated, the implants act as an anchor for various fixed or removable prosthetic solutions that can be used to improve or restore a patient's mastication function.

## 3. INDICATIONS FOR USE

BTI Dental Implants can be used for single or multiple restorations in upper and lower jaw, taking into account the following recommendations:

The recommended use of implants with a diameter of 2.5 mm is always splinted. Do not use this implant for multiple prostheses screwed directly onto the implant.

The recommended use of INTERNA 3.0 implants with a diameter of 3 mm is always splinted, except in the lower or lateral upper incisors with no occlusion in which they can be used separately as single implants.

The recommended use of BTI implants with a diameter of 3.3 mm and TINY® implants with a diameter of 3 mm is always splinted, except in lower incisors or upper lateral incisors with no occlusion in which they can be used separately as single implants.

Implants of 4.5 mm length are indicated for partial or complete prosthetic restorations with severe atrophy of upper or lower jaw. In the case of the lower jaw with a residual bone of at least 5 mm above the dental nerve. They must always be splinted.

## 4. INTENDED USER AND INTENDED PATIENT GROUP

It must be used by dental healthcare professionals with specific training in dental implantology and in BTI dental products.

Patients partially or totally edentulous who are going to be subjected to oral implantology techniques and therefore to oral or maxillofacial surgery.

Permanent implants are not indicated in children who have not completed their development phase.

## 5. CONTRAINDICATIONS

BTI dental implants should not be used if the bone crest is not high or wide enough to ensure initial fixation and stability of the implant.

BTI dental implants are not intended for use in children who have not concluded their growth and development stages.

BTI dental implants are not suitable for use in patients who are allergic or hypersensitive to the material used (commercially pure titanium).

## 6. WARNINGS

- BTI dental implants should be accepted only if the factory packaging and labelling arrive intact. Avoid setting the implant in the event of the plastic seal being broken. Contact your distributor if the package is opened or altered.

- The abutments indicated for the BTI L4.5 INTERNA implants are exclusively MUL-TI-IM ref. INTMIPE45XX and INTMIPU45XX.
- Do not use BTI dental implants after the expired date indicated in the packaging of the device.
- The re-use and/or re-sterilisation of these products may cause risk of infection and/or cross contamination.
- Lack of osseointegration of the implant and loss of surrounding bone can occur due to inadequate technique or local or systemic patient factors (e.g. insufficient bone available, bone of poor quality, lack of oral hygiene, smoking, blood disorders, diabetes, infection, periodontitis, bruxism, hypertension, antidepressant drugs, antiresorptive therapies, mental disorders, some chemotherapy, radiation therapy or anticoagulation therapy).
- Damage to the dental nerve can result in anaesthesia, paraesthesia or dysaesthesia.
- Do not exceed the insertion torques recommended by BTI as this may cause bone necrosis.
- In case of ingestion of the piece by the patient, refer the patient to the hospital emergency department for appropriate treatment.

## 7. PRECAUTIONS

- An examination is essential, and also panoramic, periapical x-rays and cone beam CT Scan are recommended in order to determine the patient's anatomical features, the periodontal condition and to determine the bone suitability. All potential implant candidates should be suitably assessed and reviewed.
- Pregnancy and breastfeeding: No additional precautionary measures to those observed in general oral implant surgery are considered, nevertheless, these special conditions shall be taken into account by the clinician, prior to implant placement. Care should also be taken if the patient has any removable dental prostheses or parafunctional tongue habits.
- BTI dental implants should only be used with the corresponding BTI original components.
- If the tightened torque recommended for the prosthetic component is greater than the insertion torque of the implant, it must be carried out delay loading.

## 8. INSTRUCTIONS FOR USE

**Step 1:** Plan the surgery based on the patient's scan and medical history.

**Step 2:** The vial contains a plastic seal as a tamper evident and three labels describing the features of the implant intended to be attached to the patient's records and for other registration purposes.

**Step 3:** Once the implant has been selected and the patient data has been filled in, proceed with the recommended surgical protocol. See section 13.1.

**Step 4:** Remove the plastic seal.

**Step 5:** Open the vial cap by turning counter clockwise until the tamper-resistant ring is detached.

**Step 6:** Pour inner content onto a sterile tray.

**Step 7:** Hold the central part of the inner holder and fit the implant mount connector (CPIX) inserted in the counter angle piece in the implant mount. Once it has taken hold, slide your fingers towards the lower section of the inner holder at the height of the cover screw whilst at the same time pressing to extract the implant.

 Take care when extracting the implant from the inner holder; extracting it through the front of the upper end may damage the surface of the implant. It must be done laterally with respect to the inner holder.

**Step 8:** Once the implant has been fixed to the implant mount connector, it should be carried, face up, to the prepared alveolus.

**Step 9:** Insert the implant with the surgical drill tech. The recommended insertion speed is 25 rpm at a maximum torque of 25 Ncm.

**Step 10:** Complete the insertion process by hand, using the LLMQ Surgical Multitorque Wrench. Complete the insertion of the implant without exceeding the maximum recommended insertion torque of 70 Ncm. For immediate loading BTI recommends an insertion torque of 40-50 Ncm since the implant must have sufficient primary stability.

**⚠** The implant mount will show the position of the implant. It is recommended that the vestibular side lines up with one side of the connection.

**Step 11:** Remove the implant mount manually. During this process, the body of the implant mount should be supported using the open-end wrench (LLAV).

**Step 12:** If the implant is going to be buried beneath the gum, use the cover screw that can be found in the other end of the inner holder, using the large hex screwdriver or the large hex screwdriver tips (for all implants other than the TINY® range). Use it to carry it to the patient's mouth. If using a TINY® implant, use a small hex screwdriver or screwdriver tip for the headless cover screws that are provided with the implant.

**⚠** The cover screws that come with the TINY® implants should be tightened at 10 Ncm; all others can be tightened at 20 Ncm.

## 9. ADVERSE EFFECTS

The implantation procedure carries certain risks, including inflammation of the implant area, loose sutures, short-term sensitivity in the area, bruising, bleeding, loss of sensation in the lower lip and chin, nose bleeds, trismus, sinusitis, peri-implantitis, periodontitis, gingivitis, fistula, mucositis, anesthesia, paresthesia, dysesthesia, local pain and maxillary or mandibular ridge bone resorption. Although these effects are most likely to only be temporary, in some isolated cases the loss of sensation has been permanent.

An infection may also occur around the implant, but this can usually be overcome with local treatment.

## 10. MRI SAFETY INFORMATION

**⚠** Warning: The RF safety of the device has not been tested. The patient may only be imaged by landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the RF coil.

Non-clinical worst case configuration review has determined that materials used in the BTI Dental Implant Systems are safe from MR-induced displacement force for a given spatial gradient and also from MR-induced torque perspective. Therefore, a patient with this device can be scanned safely in an MR system under the following conditions:

|                                       |                                                                                                                                                                                           |
|---------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Device Name                           | BTI Dental Implant Systems – Prosthetic Components.                                                                                                                                       |
| Static Magnetic Field Intensity (B0)  | ≤ 3.0T                                                                                                                                                                                    |
| Maximum spatial gradient of the field | 30 T/m (3,000 gauss/cm)                                                                                                                                                                   |
| RF Excitation                         | Circular polarisation (CP)                                                                                                                                                                |
| RF Transmission Coil Type             | For body transmit coil, landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the coil.<br>Extremity T/R coils permitted.<br>Excludes Head T/R coil. |
| Operating Mode                        | Normal Operating Mode in the allowed imaging zone                                                                                                                                         |
| Maximum Whole-Body SAR                | 2 W/kg (Normal Operating Mode)                                                                                                                                                            |
| Maximum Head SAR                      | Not evaluated for head landmark                                                                                                                                                           |
| Scan Duration                         | No specific constraints due to implant heating                                                                                                                                            |

## 11. MAINTENANCE, STORAGE AND DISPOSAL

Do not expose to direct sunlight.

It is recommended to discard safely potentially contaminated or no longer usable medical devices as dental health care waste in accordance with local regulation and government legislation.

## 12. NOTICE REGARDING SERIOUS INCIDENTS

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority. The contact information for the manufacturer is as follows: <https://www.qualifiedperson@bti-implant.es>

## 13. FURTHER INFORMATION

- It is the user's responsibility to examine the product and check that its components are in an appropriate condition for their intended use.
- The BTI dental implants surface may show signs of white calcium deposits or lack of wet appearance due to having been exposed to extreme temperatures during transport and storage. These signs have no influence on the correct working and safety of the product.

### EUDAMED

When the European database on Medical Devices becomes available, the summary of safety and clinical performance of this medical device will be available via Eudamed at <https://ec.europa.eu/tools/eudamed>. In order to do so, the UDI-DI code, available on the label of the implant, must be entered.

The UDI code corresponds to a series of numeric or alphanumeric characters that allow the traceability of the medical device and is displayed on the label in ICAD format and in HRI (human readable interpretation) format.

### 13.1 DRILLING PROTOCOL

For the appropriate drilling sequence for the implant you are using, see CAT221 for recommended drilling sequences for BTI EXTERNA implants and CAT218 for recommended drilling sequences for BTI INTERNA implants.

## 14. SYMBOLS

For a description of the symbols that appear on product labels and in these instructions, refer to the eMA087 guide.

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