
BTI DENTAL IMPLANTS: UnicCa[®]

en	BTI DENTAL IMPLANTS: UNNICA®.....	(Page 4)
es	IMPLANTES DENTALES BTI: UNNICA®.....	(Página 7)
de	BTI-ZAHNIMPLANTATEN: UNNICA®.....	(Seite 10)
fr	IMPLANTS DENTAIRES BTI: UNNICA®.....	(Page 13)
it	IMPIANTI DENTALI BTI: UNICCA®.....	(Pagina 16)
pt	IMPLANTES DENTARIOS BTI: UNICCA®.....	(Pagina 19)
ro	IMPLANTURILE DENTARE BTI: UNICCA®.....	(Pagină 22)
cs	DENTALNI IMPLANTATY BTI: UNICCA®.....	(Stránka 25)

BTI DENTAL IMPLANTS: UNICCA®

1. PRODUCT DESCRIPTION

BTI dental implants are endosseous implants made of pure grade 4 titanium with different platforms, diameters and lengths that have a surface obtained through selective acid-etching. The BTI unicCa® surface has the same topography as the traditional BTI implants but, in addition, the BTI unicCa® surface has a chemical modification with calcium ions that protects the surface and makes them hydrophilic. The visual appearance of the surface of BTI unicCa® implants is wet in normal conditions of temperature and relative humidity.

The BTI unicCa® dental implants are available with two different connection types: Interna and Externa and with a wide variety of diameters and lengths to provide solutions to the different clinical situations that can be found in clinical practice.

The BTI unicCa® dental implants are sterile, in a screw-cap container that also includes the implant mount or transporter and the corresponding healing screw in the lower part of the grip.

2. INDICATIONS FOR USE

The BTI unicCa® Dental Implant system for oral implant surgery is to be used for the partial or total replacement of teeth in edentulate patients. Once osseointegrated, the unicCa® implants act as an anchor for various fixed or removable prosthetic solutions that can be used to improve or restore a patient's mastication function.

The recommended use of Tiny unicCa® implants with a diameter smaller than 3.3 mm is splinted to 1 and, better yet, 2 implants.

Tiny unicCa® implants measuring 7.5-6.5 mm in length, and those measuring 6.5-5.5 mm in length, must always be used splinted to at least another 2 implants, and the surgical procedure shall be carried out in two phases without immediate loading. They must not be used for angled or removable prostheses.

The recommended use of INTERNA 3.0 unicCa® implants

with a diameter of 3.3 mm is always splinted, except in upper lateral incisors, and lower incisors or upper lateral incisors with no occlusion in which they can be used separately as single implants. Do not use this implant for multiple prostheses screwed directly onto the implant

The recommended use of INTERNA 3.0 unicCa® implants with a diameter of 3 mm is always splinted, except in the lower or lateral upper incisors with no occlusion in which they can be used separately as single implants. Do not use this implant for multiple prostheses screwed directly onto the implant.

The recommended use of INTERNA 3.0 unicCa® implants with a diameter of 2.5 mm is always splinted. Do not use this implant for multiple prostheses screwed directly onto the implant.

Implants of 4.5mm length are indicated for partial or complete prosthetic restorations with severe atrophy of upper or lower jaw. In the case of the lower jaw with a residual bone of at least 5mm above the dental nerve. They must always be splinted. The transepithelial abutments indicated for these implants are exclusively Multi-Im ref. INTMIPE45XX and INTMIPU45xx.

3. INSTRUCTIONS FOR USE

The new unicCa® threaded implant carrier system makes it possible to transfer the sterile implant from the packaging to the implantation site and to position it with accuracy.

Step 1: Select an implant based on the patient's orthodontic/prosthetic programme and x-rays.

Step 2: The vial contains labels describing the features of the implant intended to be attached to the patient's records and for other registration purposes.

Step 3: Open the screw-cap container and tip its content on the sterile sheet.

Step 4: Hold the central part of the pliers and fit the implant-holder connector (CPIX) inserted in the CAHEX

counter angle piece in the implant-holder. Once it has taken hold, slide your fingers towards the lower section of the pliers at the height of the locking screw whilst at the same time pressing to extract the implant.

 Take care when extracting the implant from the grip; extracting it through the front of the upper end may damage the surface of the implant. It must be done laterally with respect to the grip.

Once the implant has been fixed to the connector, it should be carried, face up, to the prepared alveolus. The recommended insertion speed is 25 rpm at a maximum torque of 25 Ncm.

Step 5: The implant should have sufficient primary stability; the recommended fixation torque is 40-50 Ncm.

 We recommend finishing off the insertion process by hand, using the LLMQ surgical multitorque wrench (30-70 Ncm). The implant mount will show the position of the implant. It is recommended that the vestibular side lines up with one side of the connection. Next, remove the mount using a large hex screwdriver (DHG17 - short; DHG24 - long; DHG36 - extra long) or a mechanical screwdriver tip: PHG17 short - PHG24 long). These tips can be inserted into a fit-in wrench (LLEC), a right-angle screwdriver (DAR) or in the LLMTP prosthetic multitorque wrench. During this process, the body of the implant mount should be supported using the openend wrench (LLAV).

Step 6: If the implant is going to be buried beneath the gum, use the healing screw that can be found in the other end of the grip, using the screwdriver tip.

Carry it to the patient's mouth using the large hex screwdriver or the large hex screwdriver tips (for all implants other than the Tiny range). If using a Tiny implant, use a small hex screwdriver or screwdriver tip for the headless healing screws that are provided with the implant.

The healing screws that come with the Tiny implants should be tightened at 10 Ncm; all others can be tightened at 20 Ncm.

If the implant is going to remain exposed, a suitable healing abutment (height and diameter) should be selected. : The healing abutment should be tightened at 20 Ncm.

4 CONTRAINDICATIONS

BTI unicCa® implants should not be used if the bone crest is not high or wide enough to ensure initial fixation and stability of the implant. Implants are not intended for use in children who have not concluded their growth and development stages.

These products are not suitable for use in patients who

are allergic or hypersensitive to the material used (commercially pure titanium)

5. WARNINGS

The implants surface may show signs of white calcium deposits due to having been exposed to extreme temperatures during transport and storage. These signs have no influence on the correct working and safety of the product.

An unsuitable patient or a poor technique can prevent the implants' bone integration and lead to a loss of the surrounding bone. Re-use and/or re-sterilisation of these products may cause risk of infection and/or cross contamination. Lack of implant bone integration may occur if there is insufficient bone available, the bone is of poor quality, there is a lack of oral hygiene in smokers or if the patient has certain medical conditions such as blood disorders or uncontrolled diabetes, infection or bruxism.

 Avoid setting the implant in the event of the seal being broken (the seal moves both laterally and longitudinally). Please note that the seal can also twist round under normal conditions when the pack is closed. The loss of the sterility barrier can cause irritation and toxicity in the patient.

6. PRECAUTIONS

All potential implant candidates should be suitably assessed and reviewed. An examination is essential, and also panoramic, periapical x-rays and CAT scan are recommended in order to determine the patient's anatomical features, the periodontal condition and to determine the bone suitability.

Pregnancy and breastfeeding: No additional precautionary measures to those observed in general oral implant surgery are considered, nevertheless, these special conditions shall be taken into account by the clinician, prior to implant placement. Care should also be taken if the patient has any removable dental prostheses or parafunctional tongue habits.

7. ADVERSE EFFECTS

Detachment of the implant (bone-integration failure) and loss of the prosthesis can occur after surgery.

8. SURGICAL COMPLICATIONS

The implantation procedure carries certain risks, including inflammation of the implant area, loose sutures, short-term sensitivity in the area, bruising, bleeding, loss of sensation in the lower lip and chin or nose bleeds. Al-

though these effects are most likely to only be temporary, in some isolated cases the loss of sensation has been permanent.

An infection may also occur around the implant, but this can usually be overcome with local treatment.

9. MRI SAFETY INFORMATION

MR Conditional

 **Warning:** The RF safety of the device has not been tested. The patient may only be imaged by land-marking at least 30 cm from the implant, or ensuring the implant is located outside of the RF coil.

Non-clinical worst case configuration review has determined that materials used in the BTI Interna Dental Implant UnicCa® Systems are safe from MR-induced displacement force for a given spatial gradient and also from MR-induced torque perspective. Therefore, a patient with this device can be scanned safely in an MR system under the following conditions:

Device Name	BTI Dental Implant System UnicCa®
Static Magnetic Field Strength (B0)	≤ 3.0T
Maximum Spatial Field Gradient	30 T/m (3,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	For body transmit coil, land-marking at least 30 cm from the implant, or ensuring the implant is located outside of the coil. Extremity T/R coils permitted. Excludes Head T/R coil.
Operating Mode	Normal Operating Mode in the allowed imaging zone
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	Not evaluated for head landmark
Scan Duration	No specific constraints due to implant heating

10. LABELING OF THE CONTAINER / SYMBOLS USED

Symbol	Symbol Title and description
	BATCH CODE.

	CATALOGUE NUMBER
	DO NOT REUSE.
	USE-BY DATE.
	STERILIZED USING IRRADIATION.
	SINGLE STERILE BARRIER SYSTEM.
	CE MARKING
	Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner
	MANUFACTURER.
	DO NOT USE IF THE CONTAINER IS DAMAGED. AND CONSULT INSTRUCTIONS FOR USE.
	DO NOT RE-STERILISE.
	CAUTION.
	CONSULT INSTRUCTIONS OF USE OR CONSULT ELECTRONIC INSTRUCTIONS FOR USE.
	DATE DE MANUFACTURE.
	Medical device.
	Unique device identifier
	Mr conditional.
	NON-PYROGENIC
	KEEP AWAY FROM SUNLIGHT

NOTE: The devices should be accepted only if the factory packaging and labeling arrive intact. Contact your distributor if the package is opened or altered.

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