
MULTI-IM and UNIT ABUTMENTS

bg	Български	ТРАНСЕПИТЕЛНИ АБАТМЪНТИ UNIT И MULTI-IM _____	(Страница 3)
cs	Čeština	TRANSEPITELÁLNÍ PILÍŘE UNIT A MULTI-IM _____	(Stránka 6)
da	Dansk	TRANSEPITELIALE ABUTMENTS UNIT OG MULTI-IM _____	(Side 9)
de	Deutsch	MULTI-IM und UNIT DISTANZHÜLSE (ABUTMENT) _____	(Seite 12)
el	Ελληνικά	ΔΙΕΠΙΘΕΛΙΑΚΟ ΣΤΗΡΙΓΜΑ UNIT ΚΑΙ MULTI-IM _____	(Σελίδα 15)
en	English	MULTI-IM and UNIT ABUTMENTS _____	(Page 18)
es	Español	TRANSEPITELIAL MULTI-IM y UNIT _____	(Página 21)
fi	Suomen	EPITEELIN LÄPÄISEVÄ ABUTMENTTI UNIT JA MULTI-IM _____	(Sivu 24)
fr	Français	PILIER MULTI-IM et UNIT _____	(Page 27)
hu	Magyar	TRANSZEPITHELIÁLIS PILLÉR UNIT ÉS MULTI-IM _____	(Oldal 30)
it	Italiano	PILASTRO TRANSEPITELIALE UNIT E MULTI-IM _____	(Pagina 33)
lt	Lietuvių	UNIT IR MULTI-IM TRANSEPITELINIAI ATRAMA _____	(Puslapis 36)
nb	Norsk Bokmål	MULTI-IM og UNIT ABUTMENTS _____	(Side 39)
pt	Português	PILAR TRANSEPITELIAL UNIT E MULTI-IM _____	(Página 42)
sv	Svenska	STÖD ABUTMENT UNIT OCH MULTI-IM _____	(Sida 45)

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MULTI-IM and UNIT abutments

1. PRODUCT DESCRIPTION

Range of transepithelial abutments with different gingival heights (and prosthetic emergence in the case of BTI UNIT and MULTI-IM transepithelial abutments) to suit the needs of each individual. They protect soft tissues and allow for conventional impression taking and intraoral digital scanning.

- UNIT and MULTI-IM TRANSEPITHELIAL ABUTMENTS: used for screw-retained restorations, both single (UNIT) and multiple (MULTI-IM), and compatible with different BTI implant connections and platforms.
- MULTI-IM COMPATIBLE TRANSEPITHELIAL ABUTMENTS: these are screwed into the implants of a specific implant system on the market, converting them into a MULTI-IM, for restoration as such.
- UNIT and MULTI-IM ANGLED TRANSEPITHELIAL ABUTMENTS: with an angulation of 17° or 30°, they correct misalignments in the case of tilted implants, making it easier to construct prosthetic structures.

MATERIALS

Component	Material	Chemical composition	Composition % (mass/mass)
UNI / MULTI-IM BTI abutment MULTI-IM BTI compatible pillars Angled BTI pillars	Grade 4 titanium	Total residues (Fe, O, C, N, H)	≤1.05
		Titanium (Ti)	Balanced
	Coating	Titanium nitride	N/A
Screw	Grade 5 titanium alloy (Ti6Al4V)	Aluminium (Al)	5,5 – 6,5
		Vanadium (V)	3,5 – 4,5
		Total impurities (Fe, O, C, N, H)	≤0,522
		Titanium (Ti)	Balanced
	Coating	WC/C	N/A

Transepithelial abutments are supplied with PEEK (UNIT/MULTI-IM and MULTI-IM COMPATIBLE) or stainless steel (ANGLED UNIT and MULTI-IM) knobs to facilitate aseptic placement during the clinical procedure.

COMPATIBLE DEVICES

The following devices may be used in combination with BTI Transepithelial Abutments, and their Instructions for Use should be consulted for general warnings and detailed information.

Compatible device	Associated IFU
Extenders and Wrenches	CAT246
BTI surgical drills	MA140
BTI Dental Implants	MA145
Articulated bars for temporary restorations	MA131
SQUARE BTI Abutments and Interfaces	MA161
Prosthetic components	MA137
Abutments, screws and healing caps	MA077

* Articulated bars for temporary restorations are not compatible with UNIT transepithelial abutments or angled transepithelial abutments. BTI dental implants are not compatible with compatible MULTI-IM transepithelial abutments.

2. INTENDED USE

BTI Dental Implant Systems consisting of dental implants, prosthetic components and dental accessory tools used in oral implantology to provide a support structure for the partial or total replacement of teeth in edentulous patients.

3. INDICATIONS FOR USE

Dental Abutments fixed to endosseous dental implants used as support for restorations:

- UNIT and MULTI-IM TRANSEPITHELIAL ABUTMENTS: screw-retained, single (UNIT) and multiple (MULTI-IM) in patients with BTI implants. These are the abutments of choice for immediate screw-retained loading in both single (UNIT) and multiple (MULTI-IM) restorations.

Expanded MULTI-IM transepithelial abutments offer a wider prosthetic emergence than the one provided by the inserted implant platform.

- UNIT and MULTI-IM ANGULATED TRANSEPITHELIAL ABUTMENTS: single (UNIT) and multiple (MULTI-IM) to correct misalignments in the case of angled BTI implants, where angulation correction is required.
- MULTI-IM COMPATIBLE TRANSEPITHELIAL ABUTMENTS: multiple screw-retained (MULTI-IM) in patients with implants from different implant systems.

4. INTENDED USER AND INTENDED PATIENT GROUP

Intended user: Dental healthcare professionals in hospital or clinical settings with specific training in oral implantology and in BTI dental products.

Intended patient group: Partially or totally edentulous male and female adults who do not have one or more conditions or pathologies listed in the contraindications.

5. CONTRAINDICATIONS

Absolute: Situations in which the product should not be used under any circumstances because the risk is too high:

- Children who have not completed bone development.

Relative: Situations in which the use of the product may be possible under certain conditions but with caution and medical supervision. The benefit must always be weighed against the risk:

- Pregnant or breastfeeding women, as implantation usually requires X-rays.
- Patients with known allergy or hypersensitivity to commercially pure titanium.

6. GENERAL WARNINGS

- Transepithelial abutments are intended for single clinical use. They are supplied sterile and must not be reused or re-sterilised as this carries a risk of infection, cross-contamination or alteration of their mechanical properties.
- Do not use after the expiry date indicated.
- Use only original BTI components; parts from other manufacturers may damage the device or cause malfunctions.
- Improper technique or incorrect patient selection may result in prosthesis failure, implant failure, or loss of surrounding bone.
- Factors that can compromise recovery include: insufficient bone quality, infections, bruxism, poor hygiene, poor patient cooperation or systemic diseases (such as diabetes).
- In case of accidental ingestion of the piece by the patient, go immediately to a hospital emergency department.

7. PRECAUTIONS

- Handle all biological samples and sharp objects in strict accordance with the centre's biosafety policies.

8. ADVERSE EFFECTS

All materials used in the manufacture of BTI transepithelial abutments are biocompatible. However, in rare cases, contact with titanium may trigger hypersensitivity reactions, such as local or generalised eczema, pruritus, pain, inflammation, erythema, gingival hyperplasia or peri-implant granulomatous reaction. In exceptional cases, yellow nail syndrome may occur, characterised by thickened, yellowish, slow-growing nails with lymphoedema, possibly related to the release of titanium ions due to corrosion. These reactions are very rare compared to other metals and usually resolve when the metal component is removed.

If you notice any adverse effects, please notify the manufacturer at the following email address: technicalmanager@bti-implant.es. By reporting adverse effects, you can help provide more information about the safety of this medical device.

9. MAGNETIC RESONANCE IMAGING (MRI) COMPATIBILITY

The radiofrequency (RF) safety of the BTI Implant System, which includes implants and prosthetic components, has not been tested. However, a non-clinical review conducted under worst-case conditions has determined that the materials used are safe in terms of the displacement force and torsion induced by MRI for a given spatial gradient.

Therefore, a patient with this device can be safely scanned in an MRI system provided that the following conditions are met and images are taken by marking a reference point at a minimum distance of 30 cm from the implant or ensuring that the implant is outside the range of the RF coil.

The specific conditions for safe scanning are:

- Static magnetic field strength (B_0) up to 3.0 T.
- Maximum spatial field gradient of 30 T/m (3,000 gauss/cm).
- RF excitation with circular polarisation (CP).
- Use of a body transmission coil only if a reference point is marked at least 30 cm from the implant or if the implant is outside the coil.
- Use of limb T/R coils permitted; use of head T/R coil excluded.
- The examination must be performed in normal operating mode within the permitted imaging zone.
- The maximum SAR for the whole body must not exceed 2 W/kg.
- Safety in relation to SAR in the head has not been evaluated, therefore cranial scanning is not recommended.
- There are no specific limitations on the duration of the scan due to heating of the implants.

10. INSTRUCTIONS FOR USE

WARNINGS

- Visually inspect the product before use to verify that the sterile barrier system is not damaged (scratches, breaks, or deformities). Do not use the product if the blister pack is damaged.
- A radiological examination is necessary to verify the correct placement of the UNIT abutment on the corresponding implant.
- The specified torque must be applied to the screw. Insufficient or excessive torque may result in the fracture of the abutment or screw, failure of the prosthesis or implant, or loss of the surrounding bone.
- The use of transepithelial abutments other than those indicated here may compromise the stability of the system and the safety of the treatment:
 - INTERNA PLATFORM 3.0 implants (Ø2.5 mm): use only MULTI-IM trans-epithelial abutments, ref. INTMIP32535xx and INTMIP325xx.
 - INTERNA PLATFORM 2.5 implants (Ø2.5 mm): use only MULTI-IM trans-epithelial abutments, ref. INTMIP241xx(*) and INTMIP235xx(*).
 - 4.5 mm length implants: use only MULTI-IM trans-epithelial abutments, ref. INTMIPE45xx(*), INTMIPU45xx(*), INTMIP345xx(*), INTMIPEEX45xx(*) and INTMIPUEX45xx(*)

PROCEDURE:

1. Select the appropriate transepithelial abutment taking into account the implant system, gingival height, size of the required prosthetic restoration, and implant connection.
2. Open the blister pack, hold the transepithelial abutment by the handle, carefully place it on the implant and gently tighten it by hand to secure it in place. Remove the handle.

- To ensure correct alignment and proper seating of the UNIT abutment, use the positioner to place the abutment on the implant connection.

3. Apply the final tightening torque of 35 Ncm using the LLMTTP torque wrench and its tip.

11. TRANSPORT AND STORAGE

- Caution: Do not expose the product to direct sunlight.
- Warning:
 - Store in a dry environment. Exposure to moisture may affect the integrity of the container.
 - Keep the product in its original packaging until ready for use; do not handle it beforehand.
 - Inspect the container before use; do not use if damaged.
 - If the container is open or has been tampered with, contact your distributor.

12. WASTE DISPOSAL

- Dispose of all used material as biological waste.
- Wear gloves during the process.
- Follow current regulations on infectious waste.

13. NOTIFICATION OF SERIOUS INCIDENTS

In accordance with Regulation (EU) 2017/745 (MDR), si if a serious incident occurs during the use of this product or as a result of its use, the user or patient must immediately notify the following address:

technicalmanager@bti-implant.es.

14. ADDITIONAL INFORMATION

BTI provides regular training courses to optimise product performance.

The Summary of Safety and Clinical Performance (SSCP) for this device is available in the European EUDAMED database (<https://ec.europa.eu/tools/eudamed>). To consult it, you must use the unique device identifier (UDI-DI) indicated on the label. If you are unable to locate it, you can request a copy directly from the manufacturer or, where applicable, through its authorised distributor.

The UDI-DI code corresponds to a series of numeric or alphanumeric characters that enable the traceability of the medical device and is displayed on the label in ICAD (Identification and Automatic Data Capture) format and in HRI (human-readable interpretation) format.

15. SYMBOLS

For a description of the symbols that appear on product labels and in these instructions, please refer to guide MA087.

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