
MULTI-IM[®] and UNIT[®] ABUTMENTS

en	MULTI-IM and UNIT ABUTMENTS	(Page 4)
es	TRANSEPITELIAL MULTI-IM y UNIT	(Página 7)
de	MULTI-IM und UNIT DISTANZHÜLSE (ABUTMENT)	(Seite 10)
fr	PILIER MULTI-IM et UNIT	(Page 13)
it	TRANSEPITELIALE MULTI-IM e UNIT	(Pagina 16)
pt	TRANSEPITELIAL MULTI-IM e UNIT	(Pagina 19)
ro	DISPOZITIV TRANSEPITELIAL MULTI-IM și UNIT	(Pagină 22)
cs	MULTI-IM a UNIT ABUTMENT	(stránka 25)
sv	MULTI-IM och UNIT FÖRANKRING	(sida 28)
no	MULTI-IM og UNIT BROSTØTTE	(Side 31)
lv	MULTI-IM un UNIT ABATMENTS	(34. lpp)
lt	MULTI-IM ir UNIT GALVUTĖ	(Puslapis 37)
hu	MULTI-IM és UNIT FELÉPÍTMÉNY	(Oldal 40)
da	MULTI-IM og UNIT ABUTMENT	(Side 43)

MULTI-IM and UNIT ABUTMENTS

1. PRODUCT DESCRIPTION

Abutments are prosthetic attachments in commercially pure titanium. They may be straight or angled, to solve all kinds of angulation problems and to improve the gingival interphase.

The abutments are sold with the coating Ti Black (screw); the Ti Golden coating (sleeve) is only applied in those that have the BTI connection.

2. INDICATIONS FOR USE

Abutments are screwed into dental implants and are used in single (UNIT®) and multiple (MULTI-IM®) restorations, or screw-mounted bridges, to make ceramic screw-mounted prostheses, overdentures with bars or hybrid prostheses.

3. INSTRUCTIONS FOR USE

Abutments are manufactured in titanium and they can be compatible with a wide variety of dental implants, by both BTI and other brands.

(*)The abutments INTMIPE45xx and INTMIPU45xx are for the exclusive use of BTI 4.5mm length implants. They must not be used with BTI implants of other lengths.

Abutments should be screw-retained to the implant at 35 Ncm torque.

Abutments are supplied sterile and are single use products.

(*)These abutments are not cleared in Canada

4. CONTRAINDICATIONS

These products are not suitable for use in patients who are allergic or hypersensitive to the material used (commercially pure titanium)

Equally, other contraindications are those derived from the surgical technique itself (anaesthesia, antibiotic therapy...) and the dental implantology technique: A lack of adequate quality and quantity of bone, infections, a lack of hygiene or cooperation by the patient, and systemic

diseases (diabetes, etc.) are potential causes that can hinder the subsequent recovery from the surgical intervention.

5. WARNING

The reuse and/or re-sterilisation of abutments may cause risk of infection and/or cross contamination.

6. PRECAUTIONS

Abutments must be used only by professionals with sufficient training. The prosthetic restoration techniques required to use these products correctly are very specialised and complex. Inadequate technique can cause the prosthesis to fail.

Abutments are sold sterile; they are for one single use and must be used before the expiry date printed on the label of the product. Do not use the product if the packaging is damaged or has already been opened

7. ADVERSE EFFECTS

BTI abutments are made in biocompatible commercially pure titanium, nevertheless some patients can present allergic reactions or hypersensitivity to this material.

8. STORAGE, CONSERVATION AND HANDLING

The product must not be removed from its original packaging or handled until the moment it is used. Doing so would cause the sterility of the product to be compromised.

NOTE: The devices should be accepted only if the factory packaging and labelling arrive intact.

Contact your distributor if the package is opened or altered.

9. MRI SAFETY INFORMATION

MR Conditional

Warning: The RF safety of the device has not been tested. The patient may only be imaged by land-marking at least 30 cm from the implant, or ensuring the implant is located outside of the RF coil.

Non-clinical worst case configuration review has determined that materials used in the BTI Dental Implant UnicCa® Systems are safe from MR-induced displacement force for a given spatial gradient and also from MR-induced torque perspective. Therefore, a patient with this device can be scanned safely in an MR system under the following conditions:

Device Name	BTI Dental Implant System UnicCa®
Static Magnetic Field Strength (B0)	≤ 3.0T
Maximum Spatial Field Gradient	30 T/m (3,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	For body transmit coil, land-marking at least 30 cm from the implant, or ensuring the implant is located outside of the coil. Extremity T/R coils permitted. Excludes Head T/R coil.
Operating Mode	Normal Operating Mode in the allowed imaging zone
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	Not evaluated for head landmark
Scan Duration	No specific constraints due to implant heating

	Sterilized using irradiation.
	Single sterile barrier system.
	CE marking
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner
	Manufacturer.
	Do not use if the container is damaged. And consult instructions for use.
	Do not re-sterilise.
	Caution.
	Consult instructions of use or consult electronic instructions for use.
	Date de manufacture.
	Medical device.
	MR conditional.

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