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
SQUARE BTI ABUTMENTS AND INTERFACES

1. PRODUCT DESCRIPTION

Titanium posts and interfaces with different gingival heights for BTI dental implants or MULTI-IM / UNIT transepithelial abutments, in rotating and non-rotating versions. They incorporate an anti-rotational system in the upper part, which allows and ensures the adjustment of the prosthetic structure to be cemented on this base.

MATERIALS

Component	Material	Chemical composition	Composition % (mass/mass)
Abutment / SQUARE BTI Interface	Titanium grade 4	Residues (Fe, O, C, N, H) in total	≤1.05
		Titanium (Ti)	Balanced
	Coating	Titanium nitride	N/A
Screw	Titanium Alloy (Ti6Al4V) grade 5	Aluminium (Al)	5.5 - 6.5
		Vanadium (V)	3.5 - 4.5
		Residues (Fe, O, C, N, H) in total	≤0.522
		Titanium (Ti)	Balance
	Coating	WC/C	N/A

2. INTENDED USE

The BTI dental implant system consists of the dental implants themselves, prosthetic components and accessory tools used in oral implant dentistry to provide a support structure for partial or total tooth replacement in edentulous patients.

3. INDICATIONS FOR USE

- With castable cylinder: Used for conventional casting techniques in conjunction with the plastic castable cylinder for single or multiple prostheses.
- Without castable cylinder: They are used as a machined titanium base on which structures are cemented using manual milling or CAD-CAM systems.

In both cases, the titanium abutment and interface and the prosthetic structure, depending on the technique, are cemented in the laboratory.

- Rotational version: Indicated for multiple prostheses screwed to the implant.
- Non-rotating version: Suitable for single prostheses (screw-retained and cemented) and multiple prostheses cemented to the implant.

4. INTENDED USER AND INTENDED PATIENT GROUP

Intended user: Professional healthcare personnel in hospital or clinical settings with specific training in oral implantology and/or dental prosthetics and training in BTI dental products.

Intended Patient Group: Partially or fully edentulous adult males and females who do not have one or more of the conditions or pathologies listed in the contraindications.

5. CONTRAINDICATIONS

Situations in which the product should not be used under any circumstances because the risk is too high.

- Children who have not completed their stage of growth and development.

Situations where use of the product may be possible under certain conditions but with caution and medical supervision. It should always be assessed whether the benefit outweighs the risk.

- Pregnant or lactating women.
- Patients with allergy to titanium.

6. GENERAL WARNINGS

- SQUARE BTI abutments/interfaces are intended for clinical use only. The product should not be reused or re-sterilised after use in a clinical procedure as there is a risk of infection, cross-contamination or alteration of its mechanical properties.
- If the product has been unpacked but not used during the surgical procedure, it may be reprocessed up to a

maximum of 6 cycles, following the cleaning, disinfection and sterilisation protocol described in the CAT246 guide. For reprocessing, the original protective packaging must first be removed.

- Use only original BTI components. The use of non-original components may cause damage to the part or functional failure.
- Do not use in laboratory techniques that involve oven or exposure to high temperatures, either on permanent or temporary components.
- Component or prosthetic structure breakage may be due to poor occlusal definition, excessive abutment/interface modifications or passive misalignment.
- Improper technique or patient selection can lead to prosthesis failure, implant failure or loss of surrounding bone.
- Factors such as insufficient bone quality, infections, bruxism, poor hygiene, poor patient cooperation or systemic diseases (e.g. diabetes) can compromise healing.
- In case of accidental ingestion of a component by the patient, immediately go to a hospital emergency department.
- Visually check the part before fitting, making sure that there is no material residue in the connection.

7. PRECAUTIONS

- Some references allow the use of an alternative screw to modify the chimney angle. See INSTRUCTIONS FOR USE for more information.
- Handle all biological specimens and sharps in accordance with the facility's biosafety policies.

8. ADVERSE EFFECTS

All materials used in the manufacture of SQUARE BTI abutments/interfaces are biocompatible. However, in rare cases, contact with titanium may trigger hypersensitivity reactions such as local or generalised eczema, pruritus, pain, inflammation, erythema, gingival hyperplasia or peri-implant granulomatous reaction. Exceptionally, yellow nail syndrome (thickened, sluggish, yellowish nails with lymphoedema) may occur, possibly related to the release of titanium ions by corrosion.

These reactions are very rare compared to other metals, and usually resolve with removal of the metal component. If you notice any adverse effects, please report them to the manufacturer at the following e-mail address: **technical-**

manager@bti-implant.es. By reporting adverse effects you can help to provide more information on the safety of this medical device.

9. COMPATIBILITY WITH MAGNETIC RESONANCE IMAGING (MRI)

The radiofrequency (RF) safety of the BTI Implant System, which includes implants and prosthetic components, has not been tested. However, non-clinical review of the worst-case configuration has determined that the materials used are safe from the point of view of MRI-induced displacement force and torsion for a given spatial gradient.

Therefore, a patient carrying this device can be safely scanned in an MRI system provided that the following conditions are met and images are taken by marking a reference point at a minimum distance of 30 cm from the implant or by ensuring that the implant is out of range of the RF coil.

Specific conditions for safe scanning are:

- Static magnetic field strength (B_0) of up to 3.0 T.
- Maximum spatial field gradient of 30 T/m (3,000 gauss/cm).
- RF excitation with circular polarisation (CP).
- Use of body coil transmission only if a reference point is marked at least 30 cm from the implant or if the implant is outside the coil.
- Use of extremity T/R coils allowed; use of head T/R coil excluded.
- The examination must be performed in normal operating mode within the permitted imaging area.
- The maximum whole body SAR should not exceed 2 W/kg.
- Safety has not been assessed for SAR at the head and therefore cranial scanning is not recommended.
- There are no specific limitations on the duration of scanning due to heating of the implants.

10. INSTRUCTIONS FOR USE



PRELIMINARY WARNINGS AND PRECAUTIONS

- Process SQUARE BTI abutments/interfaces prior to use following only the protocol described in the CAT246 cleaning, disinfection and sterilization guide.
- Apply the indicated torque to the screw.
 - If more than the recommended torque is applied, the screw may break.

- If a lower torque is applied, the abutment may loosen.
- The angulation of the chimney of the prosthetic framework is for correcting esthetic exits, not for correcting implant angulations.
- SQUARE BTI abutments/interfaces can be cut to suit prosthetic needs. The cutting limit is 4 millimetres from the prosthetic platform to the marked groove.

- Use of improper cement and/or improper cementation of the prosthesis may result in fracture of the prosthesis.

1. Select the abutment or interface according to the technique and type of prosthesis required. In the following table you will find the different options with the recommended angulations and torques. Some abutment and interface references allow the use of another screw to obtain a greater chimney angulation. The table indicates which references allow this, the screw to be used, the maximum angulation and the recommended torque.

DIRECT TO IMPLANT						
Implant Connection	SQUARE Abutment	Maximum angulation* (Ncm)	Recommended torque* (Ncm)	Angled screw and maximum angulation	Recommended torque with angled screw (Ncm)	Prosthetic Ø (mm)
INTERNA 3.0	INPSP30	15°	0°: 35 ≤15°: 20	-	-	4
	INPSP31	15°	0°: 35 ≤15°: 20	-	-	4
	INPSP32	14°	0°: 35 ≤14°: 20	-	-	4
	INPSP33	7°	0°: 35 ≤7°: 20	-	-	4
	INPSP3R0	15°	0°: 35 ≤15°: 20	-	-	3.8
NARROW INTERNA	INPSPE0	15°	0°: 35 ≤15°: 20	-	-	4
	INPSPE1	15°	0°: 35 ≤15°: 20	-	-	4
	INPSPE2	7°	0°: 35 ≤7°: 20	-	-	4
	INPSPE3	5°	0°: 35 ≤5°: 20	-	-	4
	INSPPER0	15°	0°: 35 ≤15°: 20	-	-	3.8
	INSPPER1	15°	0°: 35 ≤15°: 20	-	-	3.8
	INSPPER2	7°	0°: 35 ≤7°: 20	-	-	3.8
	INSPPER3	5°	0°: 35 ≤5°: 20	-	-	3.8
	INPSPEX0	15°	0°: 35 ≤15°: 20	-	-	4.4
	INPSPEX1	15°	0°: 35 ≤15°: 20	-	-	4.4
	INPSPEX2	10°	0°: 35 ≤10°: 20	-	-	4.4
	INPSPEX3	6°	0°: 35 ≤6°: 20	-	-	4.4
UNIVERSAL INTERNA	INSPU0	15°	0°: 35 ≤15°: 20	INTTUHA: 30°	20	4.8
	INSPU1	15°	0°: 35 ≤15°: 20	INTTUHA: 24°	20	4.8
	INSPU2	14°	0°: 35 ≤14°: 20	-	-	4.8
	INSPU3	10°	0°: 35 ≤10°: 20	-	-	4.8
	INSPUR0	15°	0°: 35 ≤15°: 20	INTTUHA: 30°	20	4.4

EXTERNA TINY	PSPT1	15°	0°: 35 ≤15°: 20	TTUHA: 30°	20	4.8
	PSPTR1	15°	0°: 35 ≤15°: 20	TTUHA: 30°	20	4.8
EXTERNA UNIVERSAL	PSPU1	15°	0°: 35 ≤15°: 20	TTUHA: 30°	20	4.8
	PSPU2	15°	0°: 35 ≤15°: 20	TTUHA: 30°	20	4.8
	PSPUR1	15°	0°: 35 ≤15°: 20	TTUHA: 30°	20	4.8

* For screw supplied.

OVER TRANSEPITHELIAL						
UNIT	SQUARE Interface	Maximum angulation* (Ncm)	Recommended torque* (Ncm)	Angled screw and maximum angulation	Recommended torque with angled screw (Ncm)	Prosthetic Ø (mm)
Emergency Ø 3.5	ISMIUPE	15°	0°: 20 ≤15°: 20	-	-	3.5
	ISMIUPEX	15°	0°: 20 ≤15°: 20	-	-	4.5
Emergency Ø 4.1	ISMIUPU	15°	0°: 20 ≤15°: 20	-	-	4.1
	ISMIUPUX	15°	0°: 20 ≤15°: 20	TTMI30: 30°	20	5.1
Emergency Ø 5.5	ISMIUPA	15°	0°: 30 ≤15°: 20	-	-	5.5
	ISMIUPAX	15°	0°: 30 ≤15°: 20	-	-	6.5
Angled / Emergency Ø 5	ISMIUPUA	15°	0°: 20 ≤15°: 20	TTMI30: 30°	20	5.5
MULTI-IM	SQUARE interface	Maximum angulation* (Ncm)	Recommended torque* (Ncm)	Angled screw and maximum angulation	Recommended torque with angled screw (Ncm)	Prosthetic Ø (mm)
Emergency Ø 3.5	ISMIPE	15°	0°: 20 ≤15°: 20	-	-	4
	ISMIPEX	15°	0°: 20 ≤15°: 20	-	-	4.5
Emergency Ø 4.1	ISMIPU	15°	0°: 20 ≤15°: 20	-	-	4.4
	ISMIPUX	15°	0°: 20 ≤15°: 20	TTMI30: 30°	20	5.1
Emergency Ø 5.5	ISMIPA	15°	0°: 20 ≤15°: 20	TTMI30: 30°	20	5.5
	ISMIPAX	15°	0°: 20 ≤15°: 20	TTMI30: 30°	20	6.5
Angled / Emergency Ø 5	ISMIPUA	15°	0°: 20 ≤15°: 20	TTMI30: 30°	20	5.1

* For screw supplied

2. Please select the technique for the prosthesis placement.
- With castable cylinder: Used for conventional casting techniques in conjunction with the plastic castable cylinder for single or multiple prostheses.
 - Without castable cylinder: These are used as a machined titanium base on which frameworks are cemented using manual milling or CAD-CAM systems.

In both cases, the titanium component and the prosthetic structure, depending on the technique, are cemented in the laboratory.

11. TRANSPORT AND STORAGE

- Do not expose to direct sunlight.
- Store in a dry environment. Exposure to moisture may affect the integrity of the packaging.
- The product should not be removed from its original packaging or handled until ready for use.
- Do not use after expiry date.
- Inspect the product before use to verify that the container is not damaged. Do not use if container is damaged.
- Contact your distributor if the container is opened or modified.

12. WASTE DISPOSAL

- Dispose of all used material as biological waste.
- Wear gloves during processing.
- Follow current regulations for infectious waste.

13. NOTIFICATION OF SERIOUS INCIDENTS

In accordance with Regulation (EU) 2017/745 (MDR), if a serious incident occurs during the use of this product or as a result of its use, the user or patient should immediately report it to the following address: **technicalmanager@bti-implant.es**

14. ADDITIONAL INFORMATION

BTI conducts regular training courses to optimise product performance.

The Summary of Safety and Clinical Performance (SSCP) for this device is available in the European EUDAMED database (<https://ec.europa.eu/tools/eudamed>). To consult it, you must use the unique device identifier (UDI-DI) indicated on the labelling. If you are unable to access EUDAMED, you can request a copy directly from the manufacturer or, where applicable, through its authorised distributor.

The UDI-DI code corresponds to a series of numeric or alphanumeric characters that allow traceability of the medical

device and is displayed on the label in ICAD (Automatic Identification and Data Capture) format and in HRI (Human Readable Interpretation) format.

15. SYMBOLS

For a description of the symbols that appear on product labels and in these instructions, see MA087.

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