
SQUARE ABUTMENTS AND INTERFACES

en	3
es	9
de	15
fr	21
it	27
pt	33
cs	39
ro	45
lv	51
lt	57
da	63
hu	69
sv	75
et	81
fi	87
el	93
no	99

SQUARE ABUTMENTS AND INTERFACES

1. PRODUCT DESCRIPTION

SQUARE BTI abutments / interfaces are attachments made of pure titanium grade 4 with Ti Golden treatment, which allow the elaboration of unitary and multiple definitive dental prostheses.

They are available in different gingival heights and are compatible with the different connections and platforms of the BTI implant system, in rotating and non-rotating versions. The abutments and interfaces intended for unitary use incorporate in the upper part an anti-rotational system, which allows and ensures the adjustment of the prosthetic structure that will be cemented on that basis.

SQUARE abutments / interfaces are supplied together with the corresponding titanium screw.

2. INTENDED USE

SQUARE BTI abutments / interfaces are designed to be placed on dental implants or MULTI-IM/UNIT abutments for partial or total tooth replacement.

3. INDICATIONS FOR USE

SQUARE BTI abutments / interfaces are used in the partial or total restoration of teeth in partially or fully edentulous patients, respectively.

SQUARE BTI abutments / interfaces are placed on dental implants or abutments, either as a support or fixation element, of the different fixed prosthetic solutions that allow to improve or restore the masticatory function of the patient.

4. INTENDED USER AND INTENDED PATIENT GROUP

SQUARE BTI abutments / interfaces should be used by health-care professionals with specific training in oral implantology and/or dentures and training in BTI dental products.

The group of patients envisaged are partially or totally toothless patients who are going to undergo oral implantology techniques and therefore oral or maxillofacial surgery.

SQUARE BTI abutments / interfaces should not be used in children who have not completed their growth and development stage.

5. CONTRAINDICATIONS

There are no particular contraindications except those derived from the surgical technique itself (anesthesia, antibiotic therapy...) and oral implantology. Lack of adequate bone quality and quantity, infections, lack of hygiene or cooperation on the part of the patient, and generalised diseases (diabetes, etc.) are potential causes that can aggravate the subsequent recovery of the surgical intervention. See instruction for use of BTI implants.

6. CAUTION

-  SQUARE BTI abutments / interfaces should be sterilised prior to use in patients regardless of prior prosthetic laboratory phase. Laboratory sterilisation will not be necessary. See section 8.2 STERILISATION.
-  Reuse and/or re-sterilisation of SQUARE abutments may result in risk of infection and/or cross-contamination.
-  Reuse and/or re-sterilisation may cause a decrease in mechanical characteristics and/or dimensions leading to failure of prosthetic restoration.
-  There are certain factors that may favour the breaking of the SQUARE BTI abutment / interface or prosthetic structure, such as poor occlusion definition, poor cementing, excessive abutment / interface modification or passive maladjustment.
-  Apply the torque indicated to the screw. If a torque greater than that indicated is applied, it may cause the screw to snap. Conversely, if a lower torque is applied, the abutment may become loose.
-  SQUARE BTI abutments / interfaces are sold with the corresponding screw, although some abutment and interface references allow the use of another screw to obtain a higher chimney angle. See section 8.1 MODE OF USE.
-  SQUARE BTI abutments / interfaces should not be subjected to high temperatures, therefore, they should not be introduced into any type of oven.
-  In case of ingestion of the piece by the patient, refer them to the emergency department of the hospital for appropriate treatment.

7. PRECAUTIONS

Prosthetic and laboratory components should only be used by suitably trained professionals. The prosthetic restoration

techniques that are required to use these products correctly are very specialised and complex. The use of an inadequate technique can cause the failure of the prosthesis.

The lack of bone integration of the implants can occur for different reasons, for example, there is not enough bone available, the bone is of poor quality, there is a lack of oral hygiene, smoking, blood disorders, diabetes, infection, periodontitis, bruxism, hypertension, antidepressant drugs, antiresorptive therapies, mental disorders, some chemotherapy, radiotherapy or anticoagulant therapy. See instruction for use of BTI implants.

8. INSTRUCTIONS FOR USE

8.1 MODE OF USE

SCREW:

In the following table we can see the maximum angle of the chimney with the screw supplied and the recommended torque, both for a structure with a straight chimney, and for a structure with the maximum allowed angle.

Some references of abutments and interfaces allow the use of another screw to obtain a greater chimney angle. The table indicates which references allow it, the screw to be used, the maximum angulation we can obtain and the recommended torque of the screw.

Finally, the prosthetic diameter of the SQUARE BTI abutment / interface reference is indicated in the table.

CONNECTION			REFERENCE	DESCRIPTION	MAXIMUM SCREW ANGU- LATION SUPPLIED	RECOMMEN- DED SCREW TORQUE SUPPLIED (Ncm)	ANGULAR SCREW AND MAXIMUM ANGULA- TION	RECOM- MENDED TORQUE ANGULAR SCREW (Ncm)	Ø PROTE- SIC (mm)
SQUARE ABUTMENT / INTERFACE	DIRECT IMPLANT	INTERNA /Universal Platform	INPSPU0	SQUARE ABUTMENT; INTERNA UNIVERSAL PLATFORM	15°	0°: 35 ≤15°: 20	INTTUHA 30°	20	4.8
			INPSPU1	SQUARE ABUTMENT; INTERNA UNIVERSAL PLATFORM; 1MM	15°	0°: 35 ≤15°: 20	INTTUHA 24°	20	4.8
			INPSPU2	SQUARE ABUTMENT; INTERNA UNIVERSAL PLATFORM; 2MM	14°	0°: 35 ≤14°: 20	-	-	4.8
			INPSPU3	SQUARE ABUTMENT; INTERNA UNIVERSAL PLATFORM; 3MM	10°	0°: 35 ≤10°: 20	-	-	4.8
			INPSPU0	SQUARE ABUTMENT; INTERNA UNIVERSAL PLATFORM	15°	0°: 35 ≤15°: 20	INTTUHA 30°	20	4.4

CONNEC- TION		REFEREN- CE	DESCRIPTION	MAXIMUM SCREW ANGU- LATION SUPPLIED	RECOMMEN- DED SCREW TORQUE SUPPLIED (Ncm)	ANGULAR SCREW AND MAXIMUM ANGULA- TION	RECOM- MENDED TORQUE ANGULAR SCREW (Ncm)	Ø PROTE- SIC (mm)	
SQUARE ABUTMENT / INTERFACE	DIRECT IMPLANT	INTERNA /Narrow Platform	INPSPE0	SQUARE ABUTMENT; INTERNA NARROW PLATFORM	15°	0°: 35 ≤15°: 20	-	-	4
			INPSPE1	SQUARE ABUTMENT; INTERNA NARROW PLATFORM; 1MM	15°	0°: 35 ≤15°: 20	-	-	4
			INPSPE2	SQUARE ABUTMENT; INTERNA NARROW PLATFORM; 2MM	7°	0°: 35 ≤7°: 20	-	-	4
			INPSPE3	SQUARE ABUTMENT; INTERNA NARROW PLATFORM; 3MM	5°	0°: 35 ≤5°: 20	-	-	4
			INSPPER0	ROTATING SQUARE ABUTMENT; INTERNA NARROW PLATFORM	15°	0°: 35 ≤15°: 20	-	-	3.8
			INSPPER1	ROTATING SQUARE ABUTMENT; INTERNA NARROW PLATFORM; 1MM	15°	0°: 35 ≤15°: 20	-	-	3.8
			INSPPER2	ROTATING SQUARE ABUTMENT; INTERNA NARROW PLATFORM; 2MM	7°	0°: 35 ≤7°: 20	-	-	3.8
			INSPPER3	ROTATING SQUARE ABUTMENT; INTERNA NARROW PLATFORM; 3MM	5°	0°: 35 ≤5°: 20	-	-	3.8
			INPSPEX0	EXPANDED SQUARE ABUTMENT; INTERNA NARROW PLATFORM	15°	0°: 35 ≤15°: 20	-	-	4.4
			INPSPEX1	EXPANDED SQUARE ABUTMENT; INTERNA NARROW PLATFORM; 1MM	15°	0°: 35 ≤15°: 20	-	-	4.4
			INPSPEX2	EXPANDED SQUARE ABUTMENT; INTERNA NARROW PLATFORM; 2MM	10°	0°: 35 ≤10°: 20	-	-	4.4
			INPSPEX3	EXPANDED SQUARE ABUTMENT; INTERNA NARROW PLATFORM; 3MM	6°	0°: 35 ≤6°: 20	-	-	4.4
			INTERNA / PLATFORM 3.0	INPSP30	SQUARE ABUTMENT; INTERNA 3.0 PLATFORM	15°	0°: 35 ≤15°: 20	-	-
		INPSP31		SQUARE ABUTMENT; INTERNA 3.0 PLATFORM; 1MM	15°	0°: 35 ≤15°: 20	-	-	4
		INPSP32		SQUARE ABUTMENT; INTERNA 3.0 PLATFORM; 2MM	14°	0°: 35 ≤14°: 20	-	-	4
		INPSP33		SQUARE ABUTMENT; INTERNA 3.0 PLATFORM; 3MM	7°	0°: 35 ≤7°: 20	-	-	4
		INPSP3R0		ROTATING SQUARE ABUTMENT; INTERNA 3.0 PLATFORM	15°	0°: 35 ≤15°: 20	-	-	3.8
		EXTERNA /UNIVERSAL PLATFORM	PSPU1	SQUARE ABUTMENT; EXTERNA UNIVERSAL PLATFORM; 1MM	15°	0°: 35 ≤15°: 20	TTUHA 30°	20	4.8
			PSPU2	SQUARE ABUTMENT; EXTERNA UNIVERSAL PLATFORM; 2MM	15°	0°: 35 ≤15°: 20	TTUHA 30°	20	4.8
			PSPUR1	ROTATING SQUARE ABUTMENT; EXTERNA UNIVERSAL PLATFORM; 1MM	15°	0°: 35 ≤15°: 20	TTUHA 30°	20	4.8

CONNECTION			REFERENCE	DESCRIPTION	MAXIMUM SCREW ANGU- LATION SUPPLIED	RECOM- MENDED SCREW TORQUE SUPPLIED (Ncm)	ANGULAR SCREW AND MAXIMUM ANGULA- TION	RECOM- MENDED TORQUE ANGU- LAR SCREW (Ncm)	Ø PRO- TESIC (mm)	
SQUARE ABUTMENT / INTERFACE	DIRECT IMPLANT	EXTERNA / TINY Narrow Platform	PSPT1	SQUARE ABUTMENT; EXTERNA TINY PLATFORM; 1MM	15°	0°: 35 ≤15°: 20	TTUHA 30°	20	4.8	
			PSPTR1	ROTATING SQUARE ABUTMENT; EXTERNA TINY PLATFORM; 1MM	15°	0°: 35 ≤15°: 20	TTUHA 30°	20	4.8	
	ABUTMENT	MULTI-IM	PLAT- FORM Ø3.5	ISMIPE	SQUARE INTERFACE; MULTI-IM; Ø3.5 PLATFORM	15°	0°: 20 ≤15°: 20	-	-	4
				ISMIPEX	EXPANDED SQUARE INTERFACE; MULTI-IM; Ø3.5 PLATFORM	15°	0°: 20 ≤15°: 20	-	-	4.5
			PLAT- FORM Ø4.1	ISMIPU	SQUARE INTERFACE; MULTI-IM; UNIVERSAL PLATFORM	15°	0°: 20 ≤15°: 20	-	-	4.4
				ISMI- PUX	EXPANDED SQUARE INTERFACE; MULTI-IM; UNIVERSAL PLATFORM	15°	0°: 20 ≤15°: 20	TTMI30 30°	20	5.1
			PLAT- FORM Ø5.5	ISMIPA	SQUARE INTERFACE; MULTI-IM; WIDE PLATFORM	15°	0°: 20 ≤15°: 20	TTMI30 30°	20	5.5
			ANGULAR / PLAT- FORM Ø5	ISMI- PUA	SQUARE INTERFACE; MULTI-IM ANGLED; UNIVERSAL PLATFORM	15°	0°: 20 ≤15°: 20	TTMI30 30°	20	5.1
		UNIT	PLAT- FORM Ø3.5	ISMIU- PE	SQUARE INTERFACE; UNIT; Ø3.5 PLATFORM	15°	0°: 20 ≤15°: 20	-	-	3.5
				ISMIU- PEX	EXPANDED SQUARE INTERFACE; UNIT; Ø3.5 PLATFORM	15°	0°: 20 ≤15°: 20	-	-	4.5
			PLAT- FORM Ø4.1	ISMIU- PU	SQUARE INTERFACE; UNIT; UNIVERSAL PLATFORM	15°	0°: 20 ≤15°: 20	-	-	4.1
				ISMIU- PUX	EXPANDED SQUARE INTERFACE; UNIT; UNIVERSAL PLATFORM	15°	0°: 20 ≤15°: 20	TTMI30 30°	20	5.1
			PLAT- FORM Ø5.5	ISMIU- PA	SQUARE INTERFACE; UNIT; WIDE PLATFORM	15°	0°: 30 ≤15°: 20	-	-	5.5
			ANGULAR / PLAT- FORM Ø5	ISMIU- PUA	SQUARE INTERFACE; UNIT ANGLED; UNIVERSAL PLATFORM	15°	0°: 20 ≤15°: 20	TTMI30 30°	20	5.5

! The chimney angulation of the prosthetic structure is to correct aesthetic outlets, not to correct implant angulations.

CYLINDERS:

It can be used with two types of cylinders:



- **With casting cylinder:** They are used for conventional techniques of manufacturing prostheses by casting together with the plastic casting cylinder for unitary prostheses.
- **Without casting cylinder:** They are used as a machined base in titanium on which zirconium bushes can be modelled with manual milling systems or CAD-CAM systems for the elaboration of prostheses or abutments with aesthetic materials such as Zirconium.

In both cases, the titanium component and the one made of zirconium will be cemented in the clinic.

ABUTMENT / INTERFACE

SQUARE BTI abutments / interfaces can be cut to fit prosthetic needs.

! The cut-off limit of the SQUARE BTI abutments / interfaces is 4 millimetres from the prosthetic platform, to the marked slot.

8.2 STERILISATION

BTI only recommends the sterilisation protocol described in the CAT-246 cleaning, disinfection and sterilisation guide.

9. SIDE EFFECTS

All materials used in the manufacture of BTI prosthetic and laboratory components are biocompatible, however, some patients may have allergic reactions or irritation from materials used in SQUARE BTI abutments / interfaces (pure titanium grade 4).

10. RESIDUAL COMPLICATIONS AND RISKS

Complications other than those that may arise during the prosthetic technique are not expected.

11. MRI SAFETY INFORMATION

! The RF safety of the BTI Implant System: UnicCa® implants + prosthetic components has not been tested. Images of the patient can only be taken by marking a landmark at a minimum distance of 30 cm from the implant, or by making sure that the implant is outside the RF coil.

Non-clinical review of the worst-case configuration has determined that the materials used in the BTI Interna Dental Implant UnicCa® Systems are safe from the standpoint of MRI-induced displacement force for a given spatial gradient and also from the perspective of MRI-induced torsion. Therefore, a patient with this device can be safely scanned on an MRI system under the following conditions:

Static Magnetic Field Intensity (B0)	≤ 3.0T
Maximum spatial gradient of the field	30 T/m (3,000 gauss/cm)

RF Excitation	Circular polarisation (CP)
RF Transmission Coil Type	In the case of the body transmission coil, a reference point must be marked at least 30 cm from the implant, or ensure that the implant is outside the coil. Allowed end T/R coils. Excluding the head T/R coil.
Operating mode	Normal operating mode in image zone allowed
Maximum SAR throughout the body	2 W/kg (Normal operating mode)
Maximum SAR in the head	Head Benchmark Not Assessed
Scan Duration	There are no specific limitations due to the heating of the implants

It is the responsibility of the user to inspect the product and verify that its components are in proper condition for the intended use.

15. LABELLING OF THE PACKAGING / SYMBOLS USED

For a description of the symbols that appear on product labels and in these instructions, refer to the MA-087 guide.

12. MAINTENANCE, STORAGE AND REMOVAL

BTI prosthetic and laboratory components do not require special storage conditions.

WASTE DISPOSAL:

1. Wear gloves to reduce the risk of infections (latex or nitrile).
2. Medical devices containing contaminated blood must be disposed of in containers suitable for potentially infectious waste.
3. General directives on hygiene and legal regulations governing appropriate disposal of infectious material must be complied with.

13. SERIOUS INCIDENT WARNING

If a serious incident occurs while or after using this product, please notify the manufacturer and your national authority. The manufacturer's contact details are: <https://www.qualifiedperson@bti-implant.es>.

14. OTHER INFORMATION

When the European database on medical devices is available, the summary of safety and clinical performance of this medical device will be available via Eudamed at <https://ec.europa.eu/tools/eudamed>. To do this, it will be necessary to enter the UDI-DI code, available on the implant label.

The UDI code corresponds to a series of numerical or alphanumeric characters that allow the traceability of the medical device and is shown on the label in ICAD format and in a humanly readable manner.

This product should only be used with the corresponding original BTI components.

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