STOP & GUIDE 3.0 KIT

EN
ES
FR
DE
IT
PT
CS
DA
FI
HU
LT
NB
SV
BG
EA

STOP & GUIDE 3.0 KIT

INSTRUCTIONS FOR USE

1. PRODUCT DESCRIPTION

Stop & Guide 3.0 is a surgical system with components to guide the placement of dental implants. It is a versatile system that can be used as a "drill stopper kit" in conventional (analog) surgery, or as a "guided surgery kit" to perform piloted, semi-guided, and fully guided surgery protocols.

2. INTENDED USE

Surgical instrument system intended for the placement of dental implants.


3. INDICATIONS FOR USE

The system is designed to allow fully guided placement of implants with a diameter of less than 4 mm whose drilling protocol has a drill with a diameter of less than 3.4mm as the largest diameter drill and whose implant holder does not exceed 4mm in diameter. Using semi-guided surgery protocols, it is possible to place larger diameter implants using this system, as long as the largest diameter drill is less than 3.4 mm.

4. STOP & GUIDE 3.0 SYSTEM COMPONENTS

The *Stop & Guide 3.0* kit includes the following components:

- A set of 12 drivers of different lengths.
- A connector for the drivers to the surgical contra-angle.
- A punch drill for flapless surgery.
- A driver tester, which is compatible with all 12 drivers.
- Specific implant mounts for fully guided implant placement protocol with the 3.0 internal connection implants range.
- A *Stop & guide 3.0* box.

PUNCH DRILL Ø 3.5, INTERNA 3.0	STOP & GUIDE DRIVER TESTER
FP335	CMCG
	

STOP & GUIDE DRIVER CONNECTOR	DRIVER STOP & GUI- DE (N°XX) INTERNA 3.0	STOP & GUIDE IMPLANT MOUNT; SLEEVE 3.0; INTERNA IMPLANT; 3.0 PLATFORM
CDCG	D3CGxx (followed by driver number)	P3CG3 P3CG325
		
STOP & GUIDE 3.0 BOX		
CCG30		
		

5. GUIDED SURGERY PROCEDURE WITH THE STOP & GUIDE 3.0 SYSTEM

To perform any of the 3 types of static guided surgery (piloted, semi-guided, or fully guided) it is necessary to have:

- A radiographic study of the patient (CBCT) and a scan of the surface of the patient's mouth.
- A surgical guide for guided surgery prepared prior to surgery, designed using surgical planning software.
- Sleeve to prepare the surgical guide for Guided Surgery. Ref: S3CG. Manufactured by BTI. The sales unit contains 5 units.
- A *Stop & guide 3.0* guided surgery kit.
- Punch drill for flapless surgery 3.5mm diameter ref. FP335. Manufactured by BTI.
- Flattening front cutting drill, ø4.1. Ref. FCFA41. Manufactured by BTI.

- A drill kit containing the required drill diameters corresponding to the dimensions of the planned implants. The recommended drilling sequence can be found in the CAT218 catalog for reference.
- A WI-75 E/KM surgical contra-angle compatible with the connector to which the drivers are coupled. Manufactured by W&H. Distributed by BTI.

5.1. SURGICAL GUIDE FEATURES

Surgical guides for guided surgery are designed and created using various software platforms that allow the merging of surface scans of the patient's mouth (STL file) with radiographic data from a CBCT scan (DICOM file). After case planning, the guide design should accommodate metal master sleeves with the following specifications: 5 mm in height and 4.12 mm in internal diameter (Ref. S3CG). The sleeves must be positioned with a 9 mm offset, measured from the implant platform to the upper edge of the sleeve.

5.2. SURGICAL GUIDED SLEEVE

CP4 Titanium cylindrical sleeves with the following dimensions:

Ref. S3CG:

INTERNA 3.0 SURGICAL GUIDE SLEEVE BTI

Outside diameter (D): 4.9mm

Outer diameter between flat surfaces (minimum external diameter): 4.6 mm

Inside diameter (d): 4.12 mm

Length (L): 5.0mm

5.3. CONTRA-ANGLE W&H WI-75 E/KM

Contra-angle compatible with any micromotor following ISO 3964 standards for handpiece connector dimensions. The connection is standard, non-illuminated, and requires no specific configuration. For more details on the WI-75 E/KM surgical contra-angle, refer to the W&H Dentalwerk Bürmoos GmbH 'Contra-angle handpieces WI-75 E/KM' user instructions.

5.4. PUNCH DRILL

The punch drill has a 3.5 mm cutting edge (Ref. FP335) designed to perform a gingival incision in a flapless surgical approach. Its outer diameter allows it to be inserted and guided through the S3CG guide sleeves.

Ref: FP335:

PUNCH DRILL Ø3,5 INTERNA 3.0

Indications for use: It should be connected to the surgical motor, does not require irrigation, and the recommended drilling speed is 50-75 rpm.

It should not be used when the driver connector is mounted on the contra-angle, as it will not be able to reach the proper cutting depth.

These products rotate in a clockwise manner (in normal surgical motor mode).

5.5. FLATTENING FRONT CUTTING DRILL

The flattening front cutting drill is a drill with a cutting edge of 4.1mm-diameter (Ref. FCFA41) for flattening the bone in irregular ridges prior to implant placement.

Indications for use: It should be connected to the surgical motor, does not require irrigation, and the recommended drilling speed is 50-75 rpm.

Before using this drill, in order to flatten the ridge, the 1.8 mm diameter drill must be inserted at least 3 mm subcrestally. This allows the initial part of the drill to be inserted in the bone to properly guide the crest flattening.

Ref: FCFA41:

FLATTENING FRONT CUTTING DRILL, Ø4.1 INTERNA 3.0

5.6. DRIVER TESTER

The driver tester consists of a handle and a head into which a driver can be screwed. This attachment facilitates the insertion of the driver through the sleeve, ensuring that the initial, smaller-diameter portion of the driver passes completely through the sleeve. This initial comprobation determines whether the implant was planned for a juxta-crestal or subcrestal position before the surgical guide is manufactured.

In cases of mesio-distal or bucco-lingual height discrepancies of the crest, or when the implant is planned for a subcrestal position, the narrow part of the driver will not contact the upper portion of the master sleeve. This discrepancy should be measured using the lateral laser markings on the drivers, and the appropriate driver must be selected to achieve the desired drilling depth.

Indications for use: This attachment to the driver aids in its insertion through the master sleeve to assess the discrepancy in the inclination of the bone crest margins and the depth at which the implant is planned to be positioned, especially in subcrestal placements.

5.7. DRIVER CONNECTOR

An element that fits on the head of a compatible contra-angle, WI-75 E/KM, with a part on which it is possible to attach the drivers.

5.8. DRIVERS STOP & GUIDE INTERNA 3.0

The drivers are cylindrical components that, once attached to the driver connector at the top, concentrically surround the drills used to prepare the implant site. Their outer diameter allows smooth insertion and guidance through the guide sleeves, ensuring the drill rotates freely inside without friction against the sleeve.

The surgical box has 12 *numbered* drivers (from 2 to 13) to facilitate their identification and with depth marks (laser marking) that help in the placement of implants when they are planned in a subcrestal position. All drivers have an identical 9 mm length initial part that is inserted through the sleeve of the surgical guide and that allows the drills to be guided. The outer diameter of this initial part of the driver is adjusted to the inner diameter of the sleeve in the surgical guide. Beyond these 9 mm, the portion of the driver closest to the contra-angle has a larger diameter, and the length of this section varies incrementally (1 mm) between two consecutive drivers.

The internal diameter of the drivers is 3.45 mm, which allows drilling with any drill with a smaller diameter.

Indications for use: When a driver and drill are placed in the contra-angle, the portion of the drill extending beyond the driver determines the effective drilling depth. The initial end of the driver stops further depth when it contacts the bone crest, ensuring only the exposed part of the drill advances into the bone. Using a shorter driver increases the drilling depth, and vice versa. The laser markings on the drivers indicate this effective drilling depth.

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5.9. GUIDED SURGERY IMPLANT MOUNT

Guided surgery implant mounts are designed with an outer diameter that matches the inner diameter of the guide sleeve, allowing precise guidance of the implant during insertion. With a height of 9 mm, matching the sleeve offset, they also provide control over the final implant positioning depth.

Indications for use: The guided surgery implant mount is attached to the implant, when the aim is to perform a fully-guided surgery protocol.

This implant mount is compatible with the prosthetic connection of BTI 3.0 internal implants (3.0 mm prosthetic connection).

6. INTENDED USER AND INTENDED PATIENT GROUP

The *Stop & Guide 3.0* surgical system should be used by clinicians with specific training in dental implantology, guided surgery for implants, and in BTI dental products.

The group of patients envisaged are partially or totally toothless patients who are going to undergo oral implantology techniques and therefore oral or maxillofacial surgery.

7. CONTRAINDICATIONS

The *Stop & Guide 3.0* surgical system is not suitable for use in patients who are allergic or hypersensitive to titanium.

The *Stop & Guide 3.0* surgical system is not indicated when there is a limitation of oral opening in the patient and/or interdental space of less than 4.6mm.

The *Stop & Guide 3.0* surgical system should not be used in

children who have not completed their growth and development stage.

For more information on the WI-75 E/KM surgical contra-angle, see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

8. WARNINGS

 The components of the Stop & Guide 3.0 surgical system as well as the S3CG sleeve, the FCFA41 drill and FP335 punch drill must be sterilised prior to use in patients. See REPROCESSING section.

 Lack of adequate bone quality and quantity, infections, lack of hygiene or cooperation on the part of the patient, teeth grinding and generalised diseases (diabetes, etc.) are potential causes that can aggravate the subsequent recovery from the surgical intervention.

 In case of ingestion of the piece by the patient, refer them to the emergency department of the hospital for appropriate treatment.

 All products must only be used with the corresponding original BTI components. Otherwise, the parts could be damaged.

 The system is optimised for BTI drill as in other systems compatibility with the marks on the drill cannot be ensured.

For more information on the WI-75 E/KM surgical contra-angle, see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

9. CAUTIONS

NUMBER OF USES

BTI recommends a maximum of 13 uses for all drills including punch drill. This value serves as a reference, as the true service life may be different depending on the type of application and/or characteristics of the material being worked with (hardness of the bone or mould). For the rest of the kit components manufactured by BTI, a maximum of 50 uses is recommended.

In any case, instruments should be discarded as soon as any physical deterioration that may affect their proper functionality is noticed, e.g., instruments with damaged, deformed, and nicked blades, as they cause vibrations and imperfections in the preparation margins, as well as uneven surfaces.

The sleeve (S3CG) is single-use product and therefore should not be reused.

For more information on the WI-75 E/KM surgical contra-angle, see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

10. ADVERSE EFFECTS

The surgical procedure carries certain risks, such as inflammation of the implant site, loose sutures, short-term tenderness of the area, bruising, bleeding, loss of sensation in the lower lip and chin, nosebleeds, trismus, sinusitis, peri-implantitis, periodontitis, gingivitis, fistula, mucositis, anesthesia, paresthesia, dysesthesia, local pain, and bone resorption of the maxillary or mandibular crest. Although these effects are most likely only temporary, in some isolated cases the loss of sensation has been permanent. An infection can also occur around the implant, but can usually be resolved with local treatment.

For more information on the WI-75 E/KM surgical contra-angle see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

11. INSTRUCTIONS FOR USE

1. Ensure that the surgical guide you have prepared after the case planning sits correctly in the patient's mouth.
2. For flapless surgery, insert the appropriate punch drill into the contra-angle without attaching the driver connector. Begin drilling clockwise at a low speed (50-75 rpm) until reaching the bone crest. If a flap is opened during the procedure, this step can be skipped.
3. Insert a driver (preferably long) into the driver tester and guide it through the sleeve. When the bone crest is flat and the implant is planned for juxta-crestal positioning, the initial part of the driver will pass completely through the sleeve until it reaches the top. If the crest is uneven or the implant is planned for subcrestal placement, the driver will not stop at the sleeve's top. The laser markings, spaced every 0.5 mm, will help estimate this distance.
 - The distance measured with the laser marking in this situation is the measurement to be added to the effective drilling depth, by selecting a shorter Driver (with a greater number), as explained in greater detail below.
 - For flattening the bone in irregular ridges, use the flattening front cutting drill (Ref. FCFA41), if needed. Connect it to the surgical motor, without irrigation, and start drilling in a clockwise manner at a low speed, 50-75 rpm.
- 3.1 Before using this drill, in order to flatten the ridge, the 1.8 mm diameter drill must be inserted 3 mm deeper than the height to be reduced. This allows the drill to be better guided.
4. Attach the driver connector to the contra-angle (Ref: WI-75 E/KM), ensuring proper alignment. Verify that it is securely mounted and exhibits no lateral or vertical movement. The friction felt during insertion into the drive head

is normal and ensures accurate guidance.

5. Insert the longest driver (No. 2) along with the largest diameter long drill that fits through it (Ø3.25 BTI drill). Rotate the drill at low speed and verify the correct concentric alignment of all components.
6. Remove the Ø3.25 drill and driver no. 2.
7. Place the starting drill (BTI long drills) and attach the driver corresponding to the desired drilling depth.
8. Proceed with the drilling protocol appropriate for the implant being placed, taking into account the patient's bone quality. Refer to CAT218 for details on BTI INTERNA implants. To achieve the desired drilling depth, you will need to alternate between drills and drivers accordingly.

	D3CG2	D3CG3	D3CG4	D3CG5	D3CG6	D3CG7
LONG DRILLS	2.5mm	3.5mm	4.5mm	5.5mm	6.5mm	7.5mm
SHORT DRILLS						

	D3CG8	D3CG9	D3CG10	D3CG11	D3CG12	D3CG13
LONG DRILLS	8.5mm	9.5mm	10.5mm	11.5mm	12.5mm	13.5mm
SHORT DRILLS		2mm	3mm	4mm	5mm	6mm

For **yuxta-crestal implant positioning**, select the driver that corresponds to the length of the implant to be placed. For example, if placing a 6.5 mm implant, driver No. 6 is selected. The driver system is optimized for use with **BTI Long drills**, so this numbering correspondence may not apply when using other drill types, depending on their length. After placing the driver, ensure the drill's depth markings align perfectly with the lower end of the driver. Drills are switched from the initial drill to the final diameter drill, following standard drilling protocols adapted to the patient's bone quality. Additionally, the five shortest drivers can be used with **BTI Short drills**, allowing for a drilling depth range of 2 to 6 mm. Drivers numbered 9 to 13 include an additional marking indicating the drilling depth with a short drill. For example, driver 11/4 allows a drilling depth of 11.5 mm with a long drill and 4 mm with a short drill.

 *Since there is some manufacturing tolerance in the length of the drills, it is important to confirm that the lower end of the driver matches the desired depth mark on the drill every time the drill is changed.

For **sub-crestal positioning**, as for yuxta-crestal positioning, the appropriate driver is first selected based on the length of the implant. For example, when a 6.5 mm implant is planned to be placed, driver No. 6 is selected, and it is ensured that the portion of the drill protruding from

the driver (which determines the drilling depth) aligns with the depth marks on the drill.

With a default offset of 9 mm (the distance between the implant platform and the top of the guide sleeve), when the implant is positioned below the ridge margin, the initial end of the driver will contact the crest before being fully inserted into the guide during drilling. The drivers have laser markings every 0.5 mm to measure this discrepancy (verified with the Driver Tester prior to drilling). This measurement is then added to the drilling length (the portion of the drill outside the driver) to accurately drill and position the subcrestal implant.

For example, if a 6.5 mm implant is planned to be positioned 2 mm sub-crestally, the laser mark visible when driver No. 6 contacts the crest will indicate 2 mm. To reach the planned depth, this 2 mm should be added to the implant length when calculating the drilling depth. In this case, using a driver that is 2 mm shorter (No. 8) would achieve the desired depth. In this case, using a driver that is 2 mm shorter (No. 8) would achieve the desired depth.

9. Once the drilling protocol is completed, remove the driver.
10. Proceed with the replacement of the conventional implant mount (pre-assembled in the implant) by the guided surgery implant mount, which can be done in different ways:
 - 10.1. If you wish to perform a half-guided procedure, also remove the surgical guide from the patient's mouth, and proceed with the placement of the implant as if it were a conventional surgery.
 - 10.2. If you want to follow a fully guided protocol, you have several options depending on the length of the implant and bone density (quality):
 - 10.2.1. When inserting implants in favorable bone quality, the first 2-3 turns of the implant can be initiated through the sleeve (without removing the guide) as in conventional surgery, using the pre-assembled implant mount and implant mount connector. After inserting the initial threads, remove the conventional implant holder and attach the guided surgery implant mount, ensuring it fits correctly. Continue placing the implant until the upper end (with a larger diameter) of the guided surgery implant mount contacts the sleeve. At this point, the implant will be positioned at the planned direction and depth (9 mm offset). Finally, remove the guided surgery implant mount and the surgical guide
 - 10.2.2. When placing extra-short implants or in cases of unfavorable bone quality, an alternative approach is to begin by inserting the implant conventionally through the guide sleeve, as described above, until the first 2-3 mm are complete and the implant achieves sufficient stability for the mount change. At this point, carefully remove the guide, detach the conventional implant mount (fixing it with the aid of an open-end wrench), and reinsert the guide. Then, connect the guided surgery implant mount

through the guide (ensuring it is securely connected) and complete the implant insertion. Finally, remove the guided surgery implant holder and the surgical guide.

10.2.3. Another method to exchange the implant holder, in cases of low bone quality or when the previous procedures are difficult, is to begin by following the same initial steps (inserting through the sleeve). After the first 2-3 mm are completed, remove the guide. With the surgical guide out of the mouth, remove the conventional implant mount and replace it with the guided surgery implant mount. Then, disinsert the implant, along with the guided surgery implant mount. Reinsert the surgical guide and fully insert the implant through the sleeve using the guided surgery implant mount. Finally, remove the guided surgery implant holder and the surgical guide.

11. Following these instructions eliminates the need for additional components during the exchange of implant holders and prevents the implant surface from coming into contact with anything other than the patient's bone.

11. 1 REPROCESSING

The components of the *Stop & guide 3.0* surgical system, the S3CG sleeve, flattening front cutting drill and the punch drills are NON-STERILE when sold and must be sterilised prior to use. For products manufactured by BTI, the reprocessing protocol described in the CAT246 cleaning, disinfection and sterilisation guide is recommended.

 The driver connector must be disassembled to be sterilised (head and washer separately).

For more information on the WI-75 E/KM surgical contra-angle, see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

12. MAINTENANCE, STORAGE AND REMOVAL

The components of the *Stop & guide 3.0* surgical system, the S3CG sleeve, flattening front cutting drill and the punch drills do not require special storage conditions.

For more information on the WI-75 E/KM surgical contra-angle see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

WASTE DISPOSAL:

It is recommended to discard safely potentially contaminated or no longer usable medical devices as dental health care waste in accordance with local regulation and government legislation.

13. SERIOUS INCIDENT WARNING

If a serious incident occurs while or after using this product, please notify the manufacturer and your national authority.

The manufacturer's contact details are:
qualifiedperson@bti-implant.es

14. OTHER INFORMATION

You must only accept devices with intact packaging and factory labelling. Contact your dealer if the container is opened or modified.

It is the responsibility of the user to examine the product and check that its components are in a suitable condition for the intended use. If the instrument does not show wear, it can be used.

For more information on the WI-75 E/KM surgical contra-angle, see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

The UDI code corresponds to a series of numerical or alphanumeric characters that allow the traceability of the medical device and is shown on the label in ICAD format and in a humanly readable manner.

For more information on the WI-75 E/KM surgical contra-angle see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

15. SYMBOLS

For a description of the symbols that appear on product labels and in these instructions, please refer to the MA087 guide.

For more information on the WI-75 E/KM surgical contra-angle see the W&H Dentalwerk Bürmoos GmbH "Contra-angle handpieces WI-75 E/KM" user instructions.

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