
Surgical Motor Accessories

bg	Български	ХИРУРГИЧЕСКИ МОТОРНИ АКЦЕСОРИ _____	(Страница 3)
cs	Čeština	CHIRURGICKÉ MOTOROVÉ PŘÍSLUŠENSTVÍ _____	(Stránka 5)
da	Dansk	KIRURGISK MOTORTILBEHØR _____	(Side 7)
de	Deutsch	ZUBEHÖR ZUM CHIRURGISCHEN MOTOR _____	(Seite 9)
el	Ελληνικά	ΑΞΕΣΟΡΙΑ ΧΕΙΡΟΥΡΓΙΚΟΥ ΚΙΝΗΤΗΡΑ _____	(Σελίδα 11)
en	English	SURGICAL MOTOR ACCESSORIES _____	(Page 13)
es	Español	ACCESORIOS DE MOTOR QUIRÚRGICO _____	(Página 15)
fi	Suomen	KIRURGISET MOOTTORITARVIKKEET _____	(Sivu 17)
fr	Français	ACCESSOIRES DU MOTEUR CHIRURGICAL _____	(Page 19)
hu	Magyar	SEBÉSZETI MOTOROS TARTOZÉKOK _____	(Oldal 21)
it	Italiano	ACCESSORI PER MOTORI CHIRURGICI _____	(Pagina 23)
lt	Lietuvių	CHIRURGINIŲ VARIKLIŲ PRIEDAI _____	(Puslapis 25)
nb	Norsk Bokmål	TILBEHØR TIL KIRURGISKE MOTORER _____	(Side 27)
pt	Português	ACESSÓRIOS PARA MOTORES CIRÚRGICOS _____	(Página 29)
sv	Svenska	TILLBEHÖR TILL KIRURGISKA MOTORER _____	(Sida 31)

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Surgical Motor Accessories

1. PRODUCT DESCRIPTION

SURGICAL MOTOR ACCESSORIES (IMPLANT MOUNT CONNECTOR (CPIXX) and IMPLANT APPLICATOR CONNECTOR (CAIxxx) are extenders for the handpiece, which facilitates implant and other instrument transport and insertion. Made of surgical grade stainless steel. They are designed to be connected to a surgical motor.

2. INTENDED USE

BTI Dental Implant Systems consisting of dental implants, prosthetic components and dental accessory tools used in oral implantology to provide a support structure for the partial or total replacement of teeth in edentulous patients.

3. INDICATIONS FOR USE

SURGICAL MOTOR ACCESSORIES are used connected to a surgical motor to facilitate handling of implant system components and instruments.

4. INTENDED USER AND INTENDED PATIENT GROUP

The devices must be used by dental healthcare professionals with specific training in dental implantology and in BTI dental products.

The intended patient group are partially or totally edentulous male and female adults.

5. CONTRAINDICATIONS

There are no contraindications.

6. WARNINGS

- A lack of hygiene or cooperation by the patient, and systemic diseases (diabetes, etc.) are potential causes that can hinder the subsequent recovery from surgical intervention.
- Reusable products should be reprocessing before use to avoid risks of infection and cross-contamination.
- Before applying torque to the CAIxxx verify that it is completely and correctly connected to the implant connection. Otherwise, implant connection could be damaged.
- Avoid an excessive working pressure, that can cause damage to the ridges of the instruments, increase the temperature and in extreme cases, break the instrument
- In case of ingestion of the piece by the patient, refer the patient to the hospital emergency department for appropriate treatment.

7. CAUTIONS

Pregnancy and breastfeeding: No additional precautionary measures to those observed in general oral implant surgery are considered, nevertheless, these special conditions shall be taken into account by the clinician, prior to implant placement. Care should also be taken if the

patient has any removable dental prostheses or para-functional tongue habits.

breakage that could affect the safety or functionality of the product.

8. INSTRUCTIONS FOR USE

8.1. How to use

8.1.1. IMPLANT MOUNT CONNECTOR:

IMPLANT MOUNT CONNECTOR are accessories that are attached at one end to the surgical motor and at the other end to the implant mount or to expander and compactors.

1. Connect the implant mount connector to the contre angle.
2. Program the surgical motor revolutions
3. Connect the corresponding device to the implant mount connector
4. Apply the required torque

8.1.2. IMPLANT APPLICATOR CONNECTOR:

IMPLANT APPLICATOR CONNECTOR are accessories that are attached at one end to the surgical motor and at the other end to the implant inner connection.

1. Connect the implant applicator connector to the contre angle.
2. Program the surgical motor revolutions
3. Connect the implant to the implant applicator connector
4. Apply the required torque

8.2. Reprocessing

SURGICAL MOTOR ACCESSORIES are multiple-use products, are supplied NON STERILE and must be reprocessed before use in patients.

BTI only recommends the reprocessing protocol described in the CAT246 cleaning, disinfection, and sterilisation guide.

For SURGICAL MOTOR ACCESSORIES, the number of uses is not limited by its use or by cleaning, disinfection and sterilisation cycles. The product should be checked after each use and cycle to ensure that it is still fit for use and shows no signs of oxidation, degradation or

9. ADVERSE EFFECTS

The use of these devices is not linked to the onset of adverse effects, although their use is directly linked to the dental implantology techniques used; thus, the adverse effects that these may cause should be taken into account.

10. STORAGE, CONSERVATION AND DISPOSAL

SURGICAL MOTOR ACCESSORIES do not require special storage conditions as long as they are kept in their original packaging.

It is recommended to discard safely potentially contaminated or no longer usable medical devices as dental health care waste in accordance with local regulation and government legislation.

11. NOTICE REGARDING SERIOUS INCIDENTS

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority. The contact information for the manufacturer is as follows: qualifiedperson@bti-implant.es

12. FURTHER INFORMATION

It is the responsibility of the user to make a visually inspection of the product and verify that its components are in proper condition for their intended use.

BTI conducts regular training courses for optimal product performance.

The UDI code corresponds to a series of numeric or alphanumeric characters that allow the traceability of the medical device and is displayed on the label in ICAD format and in HRI (human readable interpretation) format.

13. LABELING OF THE CONTAINER / SYMBOLS USED

For a description of the symbols that appear on product labels and in these instructions, refer to the MA087 guide.

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